
easyUS™ Ultrasonic Surgical System Handpiece Instructions for Use

For handpiece eHP01

Ezisurg Medical Co., Ltd.

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easyUS™ Ultrasonic Surgical System Handpiece

Please read all information carefully before use.

Rev.XH
2026-06-23

Important!

Please read all information carefully.

1. This manual is designed to provide instructions for use of the handpiece eHP01 of Ultrasonic Surgical System. It is not a reference to surgical techniques.
2. The handpiece should be used with the generator and single use shears (SHC14A, SHC23A, SHC36A, SHC45A) of Ultrasonic Surgical System. Please refer to the compatibility of the handpiece eHP01.

Device Description

The Ultrasonic Surgical System converts high frequency current into ultrasonic oscillation to mechanically fragment soft-tissue cells upon contact with a vibrating tip for cutting and coagulating tissue during surgery. It mainly consists of a generator (with software embedded), a footswitch, a handpiece and a single use shears. The shears and the handpiece have a torque wrench and/or a test tip as spare parts. The torque wrench is to secure the test tip or shears to the handpiece. The test tip is to connect the handpiece and is used to help to calibrate and diagnose the handpiece performance. Among them, the generator outputs electrical energy (controlled by the software); the footswitch is used to active MAX or MIN output; the handpiece converts the electrical energy to mechanical motion; the shears connects to the handpiece and pass the mechanical motion to the distal end of the blade in the shears.

The Ultrasonic Surgical System does not include an integrated irrigation/suction system.

The Ultrasonic Surgical System does not support the function of patient physiological parameter monitoring.

Handpiece mainly consists of a transducer, a handpiece hull, a cable and a plug. The handpiece should work with Ultrasonic Surgical System generator. The handpiece converts electrical energy to mechanical motions. There is a fixed cable on handpiece which has a plug on the other side to connect to the receptacle in front of a generator.

Handpiece is packed with a torque wrench and a test tip. All of the handpiece, the torque wrench and the test tip are not sterilized when manufactured and packaged non-sterile. The torque wrench and the test tip packaged with the handpiece are only used outside the operating room.

Handpiece is re-usable. Medical institutions should sterilize the handpiece prior to use. The handpiece performance might decrease during general usage. The handpiece should be replaced when generator notifies errors about the handpiece,

See Figure 1 for the diagram of the handpiece.


Figure 1 Diagram of Handpiece

See Table 1 for the product codes of the handpiece and the compatibility of generators and instruments.

Table 1 Handpiece Specification

PRODUCT CODE	RESONANT FREQUENCY	CABLE LENGTH	COMPATIBLE GENERATOR	COMPATIBLE INSTRUMENTS
eHP01	55.5±0.5KHz	2.85m	ES01	SHC14A SHC23A SHC36A SHC45A

Theory of Operation

The Ultrasonic Surgical System is for hemostatic cutting and/or coagulation. The blade in shears vibrates longitudinally at 55.5 kilohertz. This ultrasonic vibration at the blade enhances its cutting ability. The same vibration seals small vessels with coagulated blood and tissue proteins. Hemostasis occurs when tissue couples with the instrument. This coupling causes collagen molecules within the tissue to vibrate and become denatured, forming a coagulum.

Device performance

The single use shears can output ultrasonic energy stably to incise soft tissue and seal small vessels when activate the hand switch or footswitch after connecting to the handpiece and generator.

Indications for Use

The easyUS™ Handpiece, when used in conjunction with the easyUS™ Instruments, is indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The instruments can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels.

Contraindications

1. The instruments are not indicated for incising bone.
2. The instruments are not intended for contraceptive tubal occlusion.
3. The Ultrasonic Surgical Systems instruments are contraindicated when, in the judgment of the physician, ultrasonic surgery would be contrary to the best interest of the patient.

Point of Use Treatment (in the operating room)

1. Clean device as soon as possible after use.
2. Remove excessive solids using disposable wipes, enzymatic foam, or per hospital procedure. Soiled devices should be separated from unsoiled devices.
3. Soiled devices should be covered with a towel dampened with tap, sterile, or critical water to prevent blood and/or debris from drying. Encrusted soil can increase the time and effort required to complete processing.
4. If the recommendations in the steps above are not possible, the device may be placed in a container and immersed in tap, sterile, or critical water and covered for transport to the processing area.

Caution: To minimize the chance of corrosion, do not use saline for cleaning.

5. Follow the **Cleaning and Sterilization** instructions below.

Transportation:

Safe storage and transportation to the reprocessing area to avoid any damage and contamination to the environment.

Target Patient Group

The intended patient population consists of adult and pediatric patients.

Intended Users

Intended users are medical professionals qualified for the ultrasonic surgical technique and surgical procedure to be performed.

Clinical Benefits

Support for successful cutting and coagulation.

Undesirable Side Effects/ Residual Risks

Potential bleeding, tissue injury, property or environmental damage, unintended harm, extended surgery or altered surgical approach may result from incorrect use or damaged

devices.

Operation

① Diagram (Figure 1)


Figure 2 Construction of Handpiece

② Installation

Prior to use, remove and discard the protective cover over the Shears Mount Surface.

Prior to each surgery, secure the test tip with the handpiece and run system test on generators to examine the handpiece performance. Please refer Troubleshooting step.

To install handpiece with shears, please refer to *easyUS™ Single Use Ultrasonic Shears Instructions for Use* for SHC14A, SHC23A, SHC36A, SHC45A.

To install the handpiece and the generator, please make sure the handpiece connector is inserted into the receptacle along the location line. Fully insert the handpiece connector.

③ Use

Power up the generator, hold the Standby touch key for 3 seconds to switch the generator to Ready mode. Ensure the shears jaw is open and not in contact with any objects. Activate system by hand switch or footswitch. The system will perform Pre-Run Test and play a test tone until the test is passed and system enters Running mode. If test failed the system enters Error mode and prompts failed error information. Please refer *easyUS™ Ultrasonic Surgical System Generator User Manual* for troubleshooting.

④ Finish

After surgery is done, remove shears from patients. Detach handpiece from shears with torque wrench. The single use shears are disposable. Disposal should be performed according to hospital requirements and regulations. The handpiece is able to reuse.

Please refer **Instructions for Reprocessing**

1. Handpiece is packaged non-sterile. Medical institutions should sterilize the handpiece prior to use.
2. Please refer to **Cleaning and Sterilization** section for details.

Cleaning and Sterilization section for details.

⑤ Troubleshooting

When generator prompts handpiece errors during testing or normal use, it is able to use test tip to diagnose the error root cause.

Detach shears from handpiece, use torque wrench to secure test tip with handpiece until to hear at least two notification tone “Ka-”. Power cycle the generator, switch to Ready mode and perform system test. Please refer System Test section on *easyUS™ Ultrasonic Surgical System Generator User Manual*.

If system test failed with test tip, please replace the handpiece as it has been broken. If system test passed with test tip, please replace the shears.

Instructions for Reprocessing

3. Handpiece is packaged non-sterile. Medical institutions should sterilize the handpiece prior to use.
4. Please refer to **Cleaning and Sterilization** section for details.

Cleaning and Sterilization

The user must ensure that cleaning and sterilization are conducted in accordance with the appropriate guidelines, standards, or national Health Authorities’ requirements. The handpiece must be cleaned and sterilized prior to each use. The handpiece has been designed to allow for thorough cleaning and safe sterilization.

Precleaning

All surgical instruments are subject to a degree of wear and tear as a result of normal use. A regular and precise visual check of the instrument should be made before each use.

The handpiece is immersible and may be soaked in a pH neutral detergent or pH neutral enzymatic detergent for a period of up to one hour prior to cleaning.

Note: The use of ultrasonic cleaners is not recommended for the handpiece.

Cleaning Regarding cleaning/disinfection, rinsing and drying, it is to distinguish between manual and automated reprocessing methods. Preference is to be given to automated reprocessing methods, especially due to the better standardizing potential and industrial safety.

Caution:

Follow the instructions and warnings issued by the suppliers of any cleaning agents and equipment used. Avoid exposing devices to hypochlorite solutions, as these will promote corrosion. Do not use cleaning agents with chlorine or chloride as they are corrosive.

Automated Cleaning Method(Preferred Method)

Following the pre-cleaning step, the device must be cleaned/disinfected in a validated automated washer disinfector in compliance to ISO 15883 or equivalent standard.

Automated cleaning program (Neutral enzyme cleaner)

STEP	MEDIUM/ CHEMICAL	HOLDING TIME	TEMPERATURE	DOSAGE
------	---------------------	-----------------	-------------	--------

PreWashing	Water	4 min	Cold	/
Emptying	/			
Washing	Belimed Protect® Neutral Enzyme Cleaner	5 min	37 °C	1.5 ml/l
Emptying	/			
Rinsing	Water	3 min	40 °C	/
Emptying	/			
Secondary rinsing	Water	5 min	40 °C	/
Drying	/	8 min	50 °C	/

Automated cleaning program (Alkaline cleaner)

STEP	MEDIUM/ CHEMICAL	HOLDING TIME	TEMPERATURE	DOSAGE
PreWashing	Water	4 min	Cold	/
Emptying	/			
Washing	neodisher MediClean forte	5 min	55 °C	5 ml/l
Emptying	/			
Rinsing	Water	3 min	40 °C	/
Emptying	/			
Secondary rinsing	Water	5 min	40 °C	/
Drying	/	8 min	50 °C	/

Note:

1. Do not expose the device to a lubrication cycle.
2. Temperatures used during the enzymatic and/or wash phases should meet the requirements of the detergent manufacturer's recommendations.
3. Manually clean the upper and lower contact rings using a cotton swab dipped in isopropyl alcohol until visibly clean.

Manual Cleaning Method

Following conditions of manual cleaning can be performed for the product:

1. 0.5-1% neodisher MediClean forte (Dr. Weigert), 10-30 minutes immersion time at 30-40 °C.
2. 0.2-0.3% Belimed Protect® Neutral Enzyme Cleaner, 10-30 minutes immersion time at 30 °C.

If other cleaners are used, select immersion times, concentrations and temperatures of the solutions as instructed by the manufacturer of the cleaner.

1. Put the device in a suitable cleaning basin.
2. Create a detergent solution as instructed by the manufacturer of the detergent.
3. Completely immerse the device in the detergent solution. Make sure that there are

no air bubbles on the products while immersed.

4. The product is brushed with a soft, nylon bristled brush for 30 seconds until all visible residues are removed. Pay special attention to the product gaps. Repeat the step for 5 times. The cleaning process takes about 10 minutes.
5. Remove the device from the cleaning basin. Rinse the device under running water for 30 seconds at room temperature.
6. Put the cleaned and rinsed device into a suitable clean basin with fresh de-ionized water. Thoroughly rinse and flush the device surface with the water for 5 minutes.
7. Dry the device through lint free towel. Insufflate cavities of products by using sterile compressed air.

Disinfection

Automated thermal disinfection

Thermal disinfection can be included as part of the automatic cleaning cycle.

A disinfection cycle of 5 min disinfection at 90°C has been validated for the device to achieve an A0 value of > 3000. Automated thermal disinfection in washer/disinfector under consideration of national requirements in regards to A0 value (see EN 15883).

Here we suggest a disinfection cycle of 5 min disinfection time at 93 °C.

Automated **disinfection** program (Glutaraldehyde Disinfectant)

The process parameters of automated disinfection with glutaraldehyde disinfectant are specified as follows:

STEP	MEDIUM/ CHEMICAL	HOLDING TIME	TEMPERATURE	DOSAGE
Disinfection	Belimed Protect™ Glutaraldehyde Disinfectant	5 min	55 °C	20 ml/l
Rinsing	Water	1 min	Cold	/
Emptying	/			
2 nd Rinsing	Water	1 min	Cold	/
Drying	/	8 min	/	/

Manual disinfection:

The following are the disinfectants, concentration percentage and contact times. Please follow hospital protocol throughout the disinfection process.

DISINFECTANT	RECOMMENDED CONCENTRATION	MINIMUM CONTACT TIME
MetriCide™ OPA Plus	Ready-to-use solution, undiluted	5 minutes

Note:

1. Any disinfection process, including tools and solutions, may influence the wear and tear on a device or equipment. In some instances, changing to another disinfectant may be required.

2. Within the applied decontamination process, ensure that detergent and disinfectant residuals are sufficiently removed. Purified or deionized water should be used during final rinsing processes, where applicable (multiple rinses may be required). Refer to the manufacturer's recommendations for the removal of disinfectant residuals.

Drying

Drying of outside of instrument through drying cycle of washer/disinfector. If needed, additional manual drying can be performed through lint free towel. Insufflate cavities of instruments by using sterile compressed air.

Functional Testing and Maintenance

Visual inspection for cleanliness of the instruments and reassembling, if required. Functional testing according to instructions of use. If necessary, perform reprocessing process again until instrument is visibly clean.

Before packaging and autoclaving, make sure that these devices have been maintained acc. to manufacturer's instruction.

Packaging

Pack the instruments in an appropriate packaging material for sterilization. The packaging material and system refer to EN ISO 11607.

Sterilization

Following the cleaning and disinfection steps above, the handpiece must be sterilized. The possible method includes steam sterilization and low temperature hydrogen peroxide plasma sterilization. Handling of the handpiece should follow hospital protocol throughout the cleaning and sterilization process.

Sterilization of instruments by applying a fractionated pre-vacuum steam sterilization process (according to EN 285/EN 13060/EN ISO 17665) under consideration of the respective country requirements.

For this instruments, preference is to be given to steam sterilization. H₂O₂ plasma sterilization is validated as alternative sterilization method.

Steam

The following steam sterilization cycle parameters are approved for use.

STERILIZER TYPE	METHOD	STERILIZATION TIME (AT TEMPERATURE)	TEMPERATURE SET POINTS
Prevacuum	Non-woven fabric wrapping or Tyvek/PE sterilization bag	5 minutes	134 °C

Note:

1. The above table includes the minimum temperature and time validated to assure sterility.
2. For steam sterilization, we recommend a drying time of 20 to 40 minutes.

3. The drying time should be determined by the medical institution based on the specific performance of its sterilizer and the configuration of the load.
4. Health Authorities in some regulated regions do not accept unwrapped sterilization methods. Please review the appropriate guidelines, standards and national Health Authorities' guidelines when determining acceptable steam sterilization process parameters for use in each respective country.

Low Temperature Hydrogen Peroxide Plasma

Sterilization is carried out in a low-temperature hydrogen peroxide plasma sterilizer. The sterilizer is the Sterrad® 100S sterilization system from Advanced Sterilization Products, Inc. Please use the certificated device and follow the instruction from device manufacturer. H₂O₂ plasma sterilization (full cycle of the long cycle program):

- H₂O₂ concentration at the beginning after injection: 58%
- Sterilization temperature: 55-60 °C
- Injection time: 6 minutes
- Holding time in diffusion: 10 minutes
- Plasma time: ≥ 5 minutes
- Injection time: 6 minutes
- Holding time in diffusion: 10 minutes
- Plasma time: ≥ 5 minutes
- Vent time: 1~2 minutes

Precautions

1. Minimally invasive procedures should be performed only by persons have adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
2. Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure. Please check carefully on the invasive parts, especially on the surface roughness, sharp edges or outcrop, which might cause safety hazard of the surgery.
3. Handpieces are not sterilized when manufactured. Nonsterile handpiece will contaminate the sterile area and expose the patient to the risk of infection. Medical institutions should sterilize the handpiece prior to use.
4. The torque wrench and test tip are not sterilized ,and are only used in non-sterile environments. They must not be used in sterile environments.
5. A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse electrosurgical instruments in liquid unless the instruments are

designed and labeled to be immersed.

6. Please use the torque wrench provided in the shears package to secure the shears and handpieces. Using torque wrench provided other than Ezisurg Medical Co., Ltd might not be able to secure shears and handpieces in appropriate torque, which might cause products broken or decrease the system performance. Please refer *easyUS™ Single Use Ultrasonic Shears Instructions for Use* for installation between handpiece and shears.
7. Operators should perform system test when shears are first in use. During the system test, please keep the jaw opened, ensure the blade is not in contact with any other instruments, drapes, patients and other objects, Safety measures (in accordance with hospital protocol) taken in the presence of aerosols should be in effect during the system test.
8. When replace shears or re-secure shears and handpiece, please switch generator to Standby mode.
9. Ultrasonic Surgical System is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. It is possible to create sparks by hitting other metal instruments. Sparks may ignite flammable gases such as bowel gas.
10. When generator prompts error messages and/or error tones, it is an indicator that the shears or handpiece is not operating properly. It may be due to the handpiece is beyond its useful life or the shears has not been attached properly, which may result in abnormally high temperatures and user or patient injury. Please refer *easyUS™ Ultrasonic Surgical System Generator User Manual* for the details of error messages.
11. If the handpiece temperature becomes uncomfortable to hold, stop use to prevent burn injury,
12. As with all energy sources (Electrosurgery, Lasers, or Ultrasound) there are concerns about the carcinogenic and infectious potential of the by-products such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
13. To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the *easyUS™* Ultrasonic shears.
14. The aerosols creates by the activations of *easyUS™* Ultrasonic shears in fatty tissue are potentially flammable.
15. To avoid user or patient injury in the event that accidental activation occurs, the shears' blade should not be in contact with the patient, drapes or flammable materials while not in use.
16. During and following activation in tissue, the shears' blade may become hot. Avoid unintended blade contact with tissue, drapes, surgical gowns, or other unintended sites after activation.
17. The handpiece meets the international safety standard EN 60601-1 for user contact and is not intended for patient contact.
18. Handle the handpiece carefully. If the binder bolt is damaged the handpiece resonant

frequency may be shifted and decrease the performance.

19. Do not bang or drop the handpiece.
20. Do not clean the handpiece electrical connector with alcohol.
21. Use handpiece eHP01 only with generator ES01 to avoid potential electric shock hazard.
22. After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
23. In case of system failure, ensure the availability of appropriate back up equipment relevant to the specific procedure.
24. The Ultrasonic Surgical System shall not be serviced or maintained while in use with a patient.
25. Successful hemostasis may require adjunct measures when the device is used on solid organs. Due to the difficulty of visualizing internal structures, proceed slowly and do not attempt to transect large masses of tissue in one activation. Avoid the division of large vascular/biliary bundles when using the instrument under these conditions.
26. Device has not been evaluated in main vessels of the central circulatory system and is not intended for use in the following named vessels: arteriae pulmonales, aorta ascendens, arcus aortae, aorta descendens to the bifurcatio aortae, arteriae coronariae, arteria carotis communis, arteria carotis externa, arteria carotis interna, arteriae cerebrales, truncus brachiocephalicus, venae cordis, venae pulmonales, vena cava superior, and vena cava inferior.

Disposal


Some internal components of the handpiece contain lead. Disposal should be performed according to local requirements and regulations, or contact Ezisurg Medical representative for assistance.

Use Times

Handpiece can be reused 100 times after clean and sterilized. However the actual life depends on many facts like how often it used, the environment, surgical technique and re-sterilization method. To ensure the handpiece is able to work properly, please attach the test tip with the handpiece then run System Test each prior to use. If test failed, please replace the handpiece.

The torque wrench and the test tip packaged with the handpiece can be reused 50 times.

Storage and transportation Conditions:

Storage of sterilized instruments in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.

Ezisurg Medical Co., Ltd.
Site :No. 399, Miaoqiao Rd., Pudong New District, 201315 Shanghai,
P.R. China
Site web : www.ezisurg.com

Shanghai International Holding Corp.GmbH(Europe)
Add: Eiffestrasse 80,20537, Hamburg, Germany
Tel: +49-40-2513175
Fax: +49-40-255726

For a patient/ user/ third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/ or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

For additional information, contact your Ezisurg Medical Co., Ltd. Sales Representative or Ezisurg Medical Co., Ltd. Customer Service.

Symbols

Caution

Affix on the package label

Keep dry

Affix on the transport
package label

This way up

Affix on the device package

Keep away from
sunlight

Affix on the transport
package label

Temperature limit

Affix on the package label
and transport package label

Humidity limitation

Affix on the package label
and transport package label

Use no hooks

Affix on the transport
package label

Stacking limit by
number

Affix on the transport
package label

Model number

Affix on the bar code label,
package label and transport
package label

Serial Number

Affix on the bar code label
and package label

Date of manufacture
Country of
manufacture: CHINA

Affix on the bar code label
and package label

	Manufacturer	Affix on the package label and transport package label
	Authorized representative in the European Community	Affix on the package label and transport package label
	Return waste to a collection system or treatment and recycling facilities. Applicable in EU. Follow decontamination instructions before returning waste.	Affix on the package label
	Device meets the requirements of Medical Device Directive 93/42/EEC, as amended by 2007/47/EC. 0197 is the code of notified body TUV Rheinland.	Affix on the package label
	Non-sterile	Affix on the package label
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Affix on the package label
	Unique device identifier	Affix on the bar code label, package label and transport package label
 www.ezisurg.com	Consult instructions for use or consult electronic instructions for use	Affix on the package label
	Medical device	Affix on the bar code label, package label and transport package label
	Atmospheric pressure limitation	Affix on the transport package label

Refer to instruction
manual/booklet

Affix on the package label

Package quantity

Affix on the transport
package label

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Shanghai International Holding Corp. GmbH (Europe)