EC REP Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffestrasse 80, 20537 Hamburg, Germany
Tel: +49-40-2513175
Fax: +49-40-255726

 Ezisurg Medical Co., Ltd.
399, Miaoqiao Rd., Pudong New District, 201315 Shanghai, P.R. China
Tel: +86-21-50456176
Fax: +86-21-50676156










CE 0197



Ver.B
P.2.0992

easyEndo™ E-Lite

- en** Single Use Powered Endoscopic Staplers and Cartridges
- pl** Jednorazowe zasilane staplery endoskopowe i ładunki
- cs** Jednorázové elektrické endoskopické staplery a zásobníky
- pt** Agrafadores Endoscópicos Alimentados de Uso Único e Cartuchos de Agrafo
- es** Grapadora endoscópica eléctrica de un solo uso y cargas
- it** Cucitrici e Cartucce Endoscopica Elettrica Monouso
- de** Elektrische Einweg-Endoskop-Klammerer und -Magazine
- fr** Agrafeuses endoscopiques motorisées et chargeurs à usage unique
- nl** Aangedreven endoscopische nietmachines en patronen voor eenmalig gebruik van
- sv** drivna endoskopiska klammerapparater och patroner för engångsbruk
- ro** Capsatoare și cartușe endoscopice electrice de unică folosință
- el** Ηλεκτρικά Ενδοσκοπικά Συρραπτικά και Φυσίγγια Μίας Χρήσης
- hu** Egyszer használatos motoros endoszkópos tűzőgépek és patronok

Please read the following instructions carefully before using this product. Failure to follow the instructions may lead to serious surgical consequences, such as anastomotic dehiscence or leakage.

Notice!

- This instruction is intended to guide the operator to use the Single Use Powered Endoscopic Staplers and Cartridges (hereinafter referred to as "The Staplers and Cartridges for Single Use"), and is not a surgical technical reference.
- Powered stapler for single use (hereinafter referred to as "the Stapler"), contained in this product can only be used for a single operation of a single patient. Reuse or reassembly of the instrument may cause potential failure, posing a risk of patient injury, and may also lead to patient infection and contagion. Do not re-use, re-assembly, re-sterilization or reprocessing.
- Cartridges can only be fired once, do not fire repeatedly.
- easyEndo™ is the registered trademark of this product.
- Cartridges from other brands cannot be used with the Stapler of this product.
- The Stapler can only be used with battery pack supplied with this product, and cannot be used with any other battery pack from other products.

Device Description

The Single Use Powered Endoscopic Staplers and Cartridges consist of Powered staplers for single use and Cartridges. Each Powered stapler for single use contains a Battery Pack, which must be installed before use.

The titanium staple is similar to door-shape, and similar to B-shape after formation. Each tissue incision is sutured by three rows of formed staples, which are arranged interlaced to ensure the capillaries are compressed and closed, to prevent erythema after tissue cut. The tissue is cut by Cutting Knife after anastomosis.

In a single surgery, per Stapler can fire at most 13 times by replacing Cartridges, and each Cartridge can be fired only once. The Stapler and Cartridge for Single Use are sterile medical devices for single use, and every part is individually packaged and sterilized.

The Shaft has an outer diameter of 12.3mm, which is fit for passing through a 12mm trocar. Trocar is optional and not contained in this product. Medical institutions should purchase trocar additionally.

Structure of the Stapler and Cartridge for Single Use is shown in Figure 1.

Structure of the Stapler is shown in Figure 2.

Structure of Cartridge is shown in Figure 3.

Figure 1 Structure diagram of The Stapler and Cartridge for Single Use

1. Cartridge 2. Powered stapler for single use

Figure 2 Structure diagram of Powered stapler for single use

- | | | | |
|---------------------------|----------------------|---------------------|------------------------|
| 1. Anvil Jaw | 2. Cartridge Jaw | 3. Shaft | 4. Articulation Knob |
| 5. Rotating Knob | 6. Close/Fire Button | 7. Open/Back Button | 8. Safety Button |
| 9. Handle | 10. Battery Pack | 11. Power Indicator | 12. Forward Indicator |
| 13. Backward Indicator | 14. Safety Indicator | 15. End of Cut line | 16. End of Staple line |
| 17. Knife Blade Indicator | | | |

Figure 3 Structure diagram of Cartridge

1. Staple Retaining Cap 2. Cover Plate 3. Staple Cartridge 4. Staple*

I. Inner II. Middle III. Outer

* Note: Staples material - Titanium Grade 2

Product model/specification

The Staplers and Cartridges have multiple specifications, with the differences shown in Table 1 and Table 2. The Staplers and Cartridges have multiple specifications, which allow surgeons to select appropriate specification and combine according to surgical needs. The Staplers and Cartridges with same length of staple line can be compatible.

Table 1 Spec. of Powered staplers for single use

Model of Stapler	Description	Length (L3)	Shaft Length (L1)	Unit: mm	
				Outer Diameter of Shaft (D1)	
P2UL45	Single Use Powered Endoscopic Stapler, Long, 45mm	440	340	12.3	
P2UM45	Single Use Powered Endoscopic Stapler, Medium, 45mm	350	250	12.3	
P2US45	Single Use Powered Endoscopic Stapler, Short, 45mm	260	160	12.3	
P2UL45G	Single Use Powered Endoscopic Stapler, Long, 45mm, Curved tip	440	340	12.3	
P2UM45G	Single Use Powered Endoscopic Stapler, Medium, 45mm, Curved tip	350	250	12.3	
P2US45G	Single Use Powered Endoscopic Stapler, Short, 45mm, Curved tip	260	160	12.3	
P2UL60	Single Use Powered Endoscopic Stapler, Long, 60mm	440	340	12.3	
P2UM60	Single Use Powered Endoscopic Stapler, Medium, 60mm	350	250	12.3	
P2US60	Single Use Powered Endoscopic Stapler, Short, 60mm	260	160	12.3	
P2UL60G	Single Use Powered Endoscopic Stapler, Long, 60mm, Curved tip	440	340	12.3	
P2UM60G	Single Use Powered Endoscopic Stapler, Medium, 60mm, Curved tip	350	250	12.3	
P2US60G	Single Use Powered Endoscopic Stapler, Short, 60mm, Curved tip	260	160	12.3	

Note: Refer to Figure 2 for L3, L1 and D1.

Table 2 Spec. of Cartridges

Unit: mm

Model of Cartridge	Staple Line Length (L2)	Open staple height (H)	Color of Staple Cartridge	Description
PN45M	45	2.0	Grey	Single Use Cartridge, Only be used with the 45 mm Stapler
PN45W	45	2.5	White	
PN45B	45	3.5	Blue	
PN45C	45	3.8	Gold	
PN45G	45	4.1	Green	
PN45T	45	4.4	Black	
PN45BC	45	3.5	Blue, ArcGrip*	
PN45CC	45	3.8	Gold, ArcGrip	
PN45GC	45	4.1	Green, ArcGrip	
PN45VM	45	Inner: 2.0 middle: 2.5 Outer: 3.5	Brown	Single Use Cartridge, Only be used with the 60 mm Stapler
PN45MT	45	Inner: 3.5 middle: 3.8 Outer: 4.1	Purple	
PN45XT	45	Inner: 3.8 middle: 4.1 Outer: 4.4	Black	
PN60M	60	2.0	Grey	
PN60W	60	2.5	White	
PN60B	60	3.5	Blue	
PN60C	60	3.8	Gold	
PN60G	60	4.1	Green	
PN60T	60	4.4	Black	
PN60BC	60	3.5	Blue, ArcGrip	Single Use Cartridge, Only be used with the 60 mm Stapler
PN60CC	60	3.8	Gold, ArcGrip	
PN60GC	60	4.1	Green, ArcGrip	
PN60VM	60	Inner: 2.0 middle: 2.5 Outer: 3.5	Brown	
PN60MT	60	Inner: 3.5 middle: 3.8 Outer: 4.1	Purple	
PN60XT	60	Inner: 3.8 middle: 4.1 Outer: 4.4	Black	

Note: Refer to Figure 3 for L2 and H.

*ArcGrip: The contact surface of cartridge is arched

WARNING: The 45mm Stapler is only used with the 45 mm Cartridges. The 60mm Stapler for single use is only used with the 60 mm Cartridges.

Type CF Applied part

1. Cartridges (refer to figure1)
2. Parts of the Stapler: Anvil Jaw, Cartridge Jaw, Shaft, Cutting Knife (refer to figure2)

Essential performance

1. No knife movement without activation of the Close/Fire Button and Open/Back Button.
2. When there is a Stapler failure, manual forced reset can return the Cutting Knife and open the jaws of the Stapler.
3. No tissue will be cut without activation of the Safety Button.

Indications /Intended Purpose

The easyEndo™ E- Lite Single Use Powered Endoscopic Staplers and Cartridges, which used together, are intended for transection, resection, and/or creation of anastomoses. The instruments have application in multiple open or minimally invasive general abdominal, gynecologic, thoracic, and pediatric surgical procedures.

The intended purpose of the device is transection and resection of tissue, and/or creation of anastomosis.

Contraindications

1. Grey Cartridges in The Stapler and Cartridge for Single Use cannot be used for tissues and organs whose thickness is less than 0.75 mm after compression or which cannot be appropriately compressed to 1.0 mm;
2. White Cartridges in The Stapler and Cartridge for Single Use cannot be used for tissues and organs whose thickness is less than 1.0 mm after compression or which cannot be appropriately compressed to 1.5 mm;
3. Blue Cartridges in The Stapler and Cartridge for Single Use cannot be used for tissues and organs whose thickness is less than 1.5 mm after compression or which cannot be appropriately compressed to 2.0 mm;
4. Gold Cartridges in The Stapler and Cartridge for Single Use cannot be used for tissues and organs whose thickness is less than 1.7 mm after compression or which cannot be appropriately compressed to 2.2 mm;
5. Green Cartridges in The Stapler and Cartridge for Single Use cannot be used for tissues and organs whose thickness is less than 2.0 mm after compression or which cannot be appropriately compressed to 2.5 mm;
6. Black Cartridges (PN45T/PN60T) in The Stapler and Cartridge for Single Use cannot be used for tissues and organs whose thickness is less than 2.3 mm after compression or which cannot be appropriately compressed to 3.0 mm;
7. Brown Cartridges in The Stapler and Cartridge for Single Use cannot be used for tissues and organs whose thickness is less than 0.88 mm after compression or which cannot be appropriately compressed to 1.5 mm;
8. Purple Cartridges in The Stapler and Cartridge for Single Use cannot be used for tissues and organs whose thickness is less than 1.5 mm after compression or which cannot be appropriately compressed to 2.25 mm;
9. Black Cartridges (PN45XT/PN60XT) in The Stapler and Cartridge for Single Use cannot be used for tissues and organs whose thickness is less than 2.25 mm after compression or which cannot be appropriately compressed to 3.0 mm;
10. It is not suitable to use this product for tissues that are prone to damage under compression;
11. It is not suitable to use this product if it is impossible to observe hemostatic condition after suture;
12. These products are sterile medical device for single use and can only be used for a single procedure. Do not reuse or use it after re-sterilization.
13. Do not use the instruments on the aorta.
14. Do not use the instruments on ischemic or necrotic tissue.

Target Patient Group

- Age: No restrictions;
- Gender: No restrictions;
- Weight & Height: No restrictions;
- Nationality: No restrictions;
- Patient health: Patients who need to perform transection, resection, and/or creation of anastomoses in multiple open or minimally invasive general abdominal, gynecologic, thoracic, and pediatric surgical procedures;
- Nationality: multiple;
- Patient state: Patient is not the OPERATOR.

Intended Users

The Staplers and Cartridges for Single Use should only be used by professional medical staff who have been fully trained in the use of The Staplers and Cartridges for Single Use techniques and surgeries.

Clinical Benefits

- Improve anastomotic success rate;
- Reduced anastomotic bleeding and anastomotic leakage rate;
- Reduced hospital stay and operative time.

Undesirable Side Effects/ Residual Risks

Undesirable side effects and risks associated with powered surgical stapling include the potential for bleeding, tissue injury, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, hypersensitivity, leakage, Biological toxicity, property or environmental damage, electrical shock, and foreign body magnetic resonance incompatibility. Also, unintended harm, extended surgery, or altered surgical approach may result from staple line failures, inability to cut, or damaged devices.

Operation instructions

1. Structure diagram of Powered stapler for single use and Cartridge (refer to Figure 1, Figure 2, Figure 3)
2. Assembling the Stapler
 - 2.1 Using sterile technique, remove the Stapler, battery pack, and Cartridge from their respective packages. To avoid damage, do not flip the Stapler, battery pack, or Cartridge into the sterile field.
 - 2.2 Install the battery pack. The battery pack must be installed prior to use. Align the battery pack with slot on lower back side of the Stapler and insert it (refer to Figure 4). The side of battery pack with buckle is faced upward. Make sure battery pack is fully inserted into the Stapler. When the battery pack is fully inserted, a click can be heard.

Notes:

- The Stapler must be used within 12 hours of inserting the battery pack. Refer to the Battery Pack Disposal section for battery pack disposal instructions.
- Rotate articulation knob, and jaws of the Stapler will rotate accordingly. If the articulation knob cannot control the jaws, the Stapler needs to be replaced.

Figure 4 Battery pack installation and removal

IV. Buckle spring plate

2.3 Assembling Cartridges

- 2.3.1 Verify that the Cartridge size matches the Stapler size to be used (e.g. use an 60 Cartridge with an 60 Stapler).

Cautions:

- Choose Cartridge with appropriate height of staple according to the tissue thickness. If tissue thickness exceeds the height range, there will be a staple malformation.
- Cartridge shall be selected with a staple line length that is the same with that of the Stapler, otherwise it is unable to load or fire. (e.g. loading a Cartridge PN60G into a Stapler P2UL60G, or loading a Cartridge PN45G into a Stapler P2UL45G.)
- Before re-assembling, wash and wipe anvil jaw and cartridge jaw with sterile solution and remove residual staples from the jaw of Stapler (refer to Figure 5). Stapler shall not be used before carrying out visual inspection to confirm that there is no staple on anvil jaw and cartridge jaw.

Figure 5 Wash the jaw of Stapler

2.3.2 Verify jaws of the Stapler are open.

2.3.3 Verify articulation knob is in the neutral position.

Note: The Cartridge cannot be assembled when articulation knob is not in neutral position or the Stapler is not fully reset.

- 2.3.4 Examine the Cartridge for the presence of a staple retaining cap. If the staple retaining cap is not in place, discard the Cartridge.

- 2.3.5 Insert a new Cartridge into the Stapler by sliding it against the bottom of cartridge jaw until the cartridge buckle is aligned with cartridge base slot, refer to Figure 6(a). Press Cartridge firmly into the slot, and there will be a click sound.

Figure 6 Cartridges assembly

Warning: If Cartridges and Stapler are not fully assembled, Cartridges may become loose and fall into body cavity during surgery.

- 2.3.6 Remove staple retaining cap from Cartridges according to the arrow in Figure 6(b).

Warning: After removing the staple retaining cap, observe the surface of the new Cartridge. The Cartridge must be replaced with another one if any colored drivers are visible. (If colored drivers are visible, the Cartridge may not contain staples.)

Note: When Cartridges are removed from sterile package, make sure staple retaining cap is taken out together. The staple retaining cap must not be removed from Cartridges before assembly is completed.

3. Disassembly

- 3.1 Ensure that the jaws of the Stapler are fully open by pressing Open/Back Button completely.
- 3.2 Adjust articulation knob to neutral position so that the jaws of the Stapler are in linear status.
- 3.3 Press the distal end of the Cartridge with thumb, while the other fingers hold the anvil jaw, according to Figure 7 to remove Cartridges from the Stapler.

Note: Do not press your thumb against the cartridge jaw during disassembly.

Figure 7 Cartridges disassembly

4. Using the instrument

- 4.1 Press the Close/Fire Button to close the jaws of the Stapler before inserting the Stapler into trocar. Keep articulation knob in neutral position Insert Stapler into an appropriate trocar or into a larger trocar through a converter. (Refer to Figure 8)

Cautions:

- Make sure anvil jaw fully exposed outside trocar cannula before opening the jaws of the Stapler in body cavity.
- Before inserting into trocar, keep articulation knob in neutral position and make sure the jaws of the Stapler are in linear and closed status.

Figure 8 Inserting into trocar

- 4.2 When the jaws of the Stapler have completely entered body cavity, press Open/Back Button until the jaws of the Stapler open automatically. (Refer to Figure 9)

Note: The shaft can rotate 360° by adjusting rotating knob.

Figure 9 Reset and adjustment of jaws

- 4.3 In body cavity, turn articulation knob to bend the jaws of the Stapler to an appropriate angle before grasping tissue. (Refer to Figure 10)

Warning: Do not turn articulation knob after grasping the tissue.

Note: The Stapler enables articulation up to 60° in either direction, but will only lock at pre-determined angles at 15°, 30°, 45°, and 60°.

Figure 10 Articulation operation

- 4.4 Close the jaws of the Stapler across the tissue to be transected by pressing Close/Fire Button. The jaws of the Stapler may be repositioned on the tissue prior to firing by pressing the Open/Back Button, allowing the jaws of the Stapler to open. (Refer to Figure 11)

Notes:

- Long press or continuous short press of the Open/Back Button can completely reset the jaws of the Stapler.
- The Stapler is equipped with a safety interlock. The Stapler will not fire the staples or cut tissue unless the green safety button is pressed.

Figure 11 Grasping tissue

V. Tissue stop wedge

Cautions:

- When positioning the Stapler on the application site, make sure that there are no obstacles, such as clips that have been inserted. Otherwise, it may result in incomplete cutting and/or staple malformation.
- The Stapler will not cut the tissue outside the end of cut line. If the tissue to be cut exceeds staple line length (45mm or 60mm) of loading cartridges, more than one Cartridge will be required.
- Any tissue beyond end of cut line will not be cut.
- 4.5 In order to fire the staples, press green safety button to deactivate safety. It is ready to be fired when the safety indicator is on. Long press or short press the Close/Fire Button repeatedly until cutting knife positioning block slides to the end of cut line. (Refer to Figures 12-14)
When cutting knife positioning block slides to the end of cut line, forward indicator and safety indicator is off, only power indicator is on.

Notes:

- Make sure the staple is completely fired. The time to complete the firing depends on staple line length (45mm or 60mm) of Cartridges and tissue status. Incomplete firing may result in incomplete cutting and/or staple malformation, which may not achieve hemostasis effectively.
- When a stalling failure occurs during firing, the system will automatically stop. Try to continue firing at this time. If the firing cannot be completed after several attempts, stop firing and long press Open/Back Button to move backward. If it is still unable to move backward, refer to "5. Manual forced reset".

Figure 12 Deactivating firing safety

Figure 13 Firing
Figure 14 Stop firing

4.6 Once the Stapler has been completely fired, release firing button, and cutting knife positioning block will automatically return to the original position before firing. Then open the jaws of the Stapler by pressing Open/Back Button. Gently remove the Stapler from the tissue.
Inspect the staple line for hemorrhage. Minor bleeding can be controlled with electrocautery or manual sutures. (Refer to Figure 15)

Figure 15 Releasing tissue

4.7 Adjust the articulation knob to the neutral position to keep the Cartridge in linear status (Refer to Figure 16); long press the Close/Fire Button to close the jaws of the Stapler; remove the Stapler from the body cavity, and unload the Cartridge from the Stapler.

Warning: For the Cartridge that have been fired once, the safety interlock will take effect to prevent from being fired again. Do not try to force the used Cartridge to fire again, otherwise it may cause damage to the Stapler or tissue.

Caution: During a single surgery, per Stapler only can fire 13 times by replace Cartridges; it is recommended to replace the Stapler after 13 times.

Note: The empty Cartridge cannot be fired again due to a safety interlock in the Stapler. Installing Cartridge that has been fired once or failing to install, Stapler properly can close the jaws of the Stapler, but the safety device will prevent the Stapler from being fired.

Figure 16 Jaws reset to linear status

5. Manual override system

In the following two cases, manual override system may be used to return the cutting knife and open the jaws of the Stapler:

- When the reset force is so great that the Stapler cannot reset;
- When there is a Stapler failure, or a battery pack failure after re-installing the battery pack;

Caution: After the manual override system is used, the Stapler is disabled and cannot be used for any subsequent firings.

There are two steps for manual override:

- Open reset cover which labeled "Manual Override" on the top of the handle. Insert the reset wrench into the first square hole (Refer to 1 in Figure 17). Turn the reset wrench counterclockwise to the end.
- Insert the reset wrench into the second square hole (Refer to 2 in Figure 17), and turn the reset wrench back and forth counterclockwise until the cutting knife completely returns to its home position. At this time, the jaws are opened to release the tissue.
After releasing the tissue, move the jaws of the Stapler away. Then the Stapler can be removed from body cavity. Discard the Stapler.

Figure 17 Manual override
VI. Reset Cover VII. Reset Wrench

Battery Pack Disposal

The battery pack contains a built-in battery consumption function, so it must be installed in the Stapler to trigger the automatic battery drain function. The battery pack does not need to be kept in the Stapler. After the battery pack is removed, it will continue to drain. Once removed, according to local regulations, the battery can be directly put into battery recycling bin or normal waste material logistics.

Before installing it into the Stapler:

If the battery pack needs to be disposed of before being installed in the Stapler (for example, this product has exceeded the expiration date on the package or the battery pack is dropped), first install the battery pack in the Stapler and then remove it to trigger the built in battery drain function.

After using the Stapler:

Before disposal, the battery pack must be removed from the Stapler. If the battery pack needs to be decontaminated before disposal, please follow the hospital's regulations.

To remove the battery pack, press the buckle spring plate and pull out the battery pack backwards (Figure 4). There is no need to disassemble the battery pack.

Standard Conventions Used

The Use of Warning, Caution, and Note Statements

Information relative to the completion of a task in a safe and thorough manner will be supplied in the form of a Warning, a Caution, or a Note statement. These statements are found throughout the documentation.

These statements should be read before continuing to the next step in a procedure.

WARNING: A Warning statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in personal injury or loss of life.

Caution: A Caution statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in damage to or destruction of the equipment.

Note: A Note statement indicates an operating or maintenance problem, practice, or condition that is necessary to accomplish a task efficiently.

Warnings

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Radiation therapy before surgery can cause tissue changes, resulting in tissue thickness exceeding the applicable range of the selected staple. Therefore, when choosing the specification of the staple, patient's various preoperative treatments should be considered.
- Tissue thickness should be carefully evaluated prior to using the Stapler. Refer to the Contraindications for proper Cartridge selection.
- Ensure tissue has not extended (extruded) proximal to the Tissue stop wedge on the jaw. Tissue forced into the Tissue stop wedge on the jaw may be transected without staples. (Refer to Figure 11)
- After removing the staple retaining cap, observe the surface of the new Cartridge. The Cartridge must be replaced with another one if any colored drivers are visible. (If colored drivers are visible, the Cartridge may not contain staples.)
- After firing, check the hemostatic condition of the staple line. Minor hemorrhage can be treated with an electrosurgical knife or manually sutured.
- When the Stapler needs to be used for multiple times in a single operation, the used Cartridges must be unloaded and replaced with a new one. The built-in safety device of the Stapler can prevent the used cartridges from being fired again. Do not attempt to forcefully breach the single-use safety device as this may result in device damage or tissue injury.
- Before re-assembling, wash and wipe anvil and cartridge base with sterile solution and remove residual staples from the jaw of Stapler.
- Make sure that there are no obstacles, such as clamps that have been inserted, between the jaws of the Stapler when the tissue to be cut is positioned; otherwise, incomplete cutting and/or staple malformation may be caused.
- This product has never been verified for cross assembly with cartridges/reloads of similar devices of other brands. Do not cross-assemble the cartridges/reloads of other brands with this product during the operation. Otherwise, it may cause stapler damage, failure to fire, staple malformation, anastomosis failure, etc.
- The Stapler and Cartridges which come into contact with bodily fluids may require special disposal handling to prevent biological contamination. The battery part needs to be handled in accordance with the Battery Pack Disposal section.
- The Staplers and Cartridges are provided sterile and are intended for use in a single procedure only. Discard after use. Do not re-sterilize. Multiple patient use may compromise the device integrity or create a risk of contamination that, in turn, may result in patient injury or illness.
- Do not use when sterile package is damaged.
- Do not use moist heat, ethylene oxide or radiation to re-sterilize or disinfect the battery pack.
- If handled incorrectly, the battery may cause a fire. Do not disassemble, heat to above 100°C, autoclave, squeeze, puncture, short-circuit external contact or charge.
- Any type of battery other than that provided with the Stapler may cause unpredictable performance problems, increased electromagnetic radiation, or reduced electromagnetic immunity of the Stapler.
- Avoid using the Staplers near other equipment or stacked with other equipment. If it is necessary to use it near or stacked with other equipment, please monitor the Stapler and other equipment to ensure normal operation.
- If it is suspected that there is electromagnetic interference between the Stapler and other equipment, reposition the Stapler or remove possible sources of interference from the room (e.g., mobile phones, radios, etc.).
- Do not modify the Stapler.

Cautions

- In a single surgery, Stapler can be fired at most 13 times by replacing Cartridges, and every Cartridge can be fired only once.

2. When endoscopic device and accessories made from different manufacturers are used in a single operation, the compatibility between the products should be verified, and the insulation and grounding should be ensured.
3. In order to prevent accidental cessation of the operation due to functional failure of the Stapler, it is recommended to prepare a full set of products as a backup during the operation.
4. Do not try to load the Cartridges while pressing and holding forward or backward button, or firing device is not reset.
5. Before inserting or removing the Stapler from trocar cannula, jaws of the Stapler should be closed and be in linear status.
6. Make sure anvil fully exposed outside trocar cannula before opening the jaws of the Stapler in body cavity.
7. Do not try to pry this product forcibly when the Stapler is inserted or removed from the trocar cannula, and during the entire operation; otherwise, it may cause Stapler failure or damage.
8. Cobalt classified as a CMR substance is present in the device in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

MRI Safety Information

 Non-clinical testing of the worst-case implants has demonstrated that the implant matrix "Staples" is MR Conditional. A patient with staples can be safely scanned in an MR system meeting the following conditions. Failure to follow these conditions may result in injury to the patient.	
Item Name/Identification	Staples
Static Magnetic Field Strength [T]	1.5 T or 3.0 T
Maximum Spatial Field Gradient [T/m]	55 T/m (5,500 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Whole body transmit coil (short body coil), Head RF receive coil or Extremity RF receive coil
Maximum Whole Body SAR [W/kg]	whole body SAR < 2.0 W/kg at 1.5 T and 3.0 T
Maximum B1 ¹ rms [μT]	B1 ¹ rms < 27.1 μT at 1.5 T Not reported for 3.0 T
Scan Duration	for 15 minutes of continuous RF (a sequence or back to back series/scan without breaks)
MR Image Artifact	The presence of this implant may produce an image artifact of 1.92 mm.
If information about a specific parameter is not included, there are no conditions associated with that parameter.	
Note: All parameters are connected by AND logical conjunction and have to be within the limits always at the same time. For MR systems that use SAR only, SAR values shall be used. For MR systems that use B1 ¹ rms and SAR, B1 ¹ rms can be used.	

Note: Staples are safe as implant components under specific conditions in 1.5T and 3.0T MRI environment.

Electrical safety performance:

1. The general electrical safety performance of the Powered stapler for single use complies with the requirements of IEC 60601-1:2005+A1:2012+A2:2020. See Table 5 below for details.

Table 5 Electrical safety characteristics

This product has following safety characteristics:

Serial number	Safety characteristics	Parameter
1	Classified according to type of protection against electric shock	Internally powered equipment
2	Classified according to degree of protection against electric shock	Type CF applied part
3	Classified according to degree of protection against ingress of water	IPX0

Serial number	Safety characteristics	Parameter
4	Classified according to safety degree when used in the environment with flammable anesthetic mixture with air or flammable anesthetic mixture with oxygen or nitrous oxide	Non-AP/APG equipment
5	Classified according to mode of operation	Continuous operation
6	Input voltage and frequency	12VDC
7	Applied part with protection against the effect of the defibrillation discharge:	No
8	Signal output or input part	No
9	Permanently installed equipment or non-permanently installed equipment	Hand-held equipment

2. Electromagnetic compatibility requirements:

For this instrument, special precautions regarding electromagnetic compatibility (EMC) must be taken, and it must be installed and used in accordance with the electromagnetic compatibility information specified in this Instruction.

Medical electrical equipment requires special EMC precautions and must be installed and used in the EMC environment specified in this document.

Portable and mobile radio frequency communication equipment may affect this instrument.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the MANUFACTURER of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of this equipment, including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic Emissions Guidance:

Guidance and manufacturer's declaration – electromagnetic emission		
The Powered stapler for single use is intended to be used in the electromagnetic environment specified below. The customer or the user should guarantee that it is used in such an environment:		
Emission test	Compliance	Electromagnetic environment - guidance
RF emission CISPR 11	Group 1	The Powered stapler for single use uses radio frequency energy only for its internal function. Therefore, its radio frequency emission is very low and is not likely to cause any interference to nearby electronic equipment.
RF emission CISPR 11	Class A	
Harmonic emission IEC 61000-3-2	Not applicable	The Powered stapler for single use can't be connected to the public low-voltage power supply network
Voltage fluctuation/flicker emission IEC 61000-3-3	Not applicable	

Electromagnetic Immunity Guidance:

Guidance and manufacturer's declaration – electromagnetic immunity		
The Powered stapler for single use is intended to be used in the electromagnetic environment specified below. The customer or the user should guarantee that it is used in such an environment:		
Immunity test	Compliance level	Electromagnetic environment - guidance

Electrostatic discharge IEC 61000-4-2	± 8 kV contact discharge ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air discharge	Floor should be wood, concrete or ceramic tile. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM fields IEC 61000-4-3	3V/m, 80MHz ~ 2.7 GHz 80%AM at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of Staplers
Proximity fields From RF wireless communications equipment	Refer to next page	
Power frequency magnetic field (50 Hz/60 Hz) IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Recommended separation distance between portable and mobile RF communication equipment and the powered endoscopic linear cutting staplers and components for single use

Nowadays, many RF wireless equipment have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. The Powered staplers for single use have been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014+A1:2020. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and The Powered staplers for single use recommended below.

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	430-470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
1845						
1970						

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
2450	2400-2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9
5500						
5785						

Sterilization

Powered staplers for single use and Cartridges was sterilized by irradiation.

Shelf life

Three years since the day of irradiation sterilization.

Intended using environment

The Stapler and Cartridge for Single Use are generally used in professional healthcare environments such as hospitals and clinics.

- Temperature: +10°C ~ +40°C ;
- Relative humidity: 30% ~ 75% ;
- Range of atmospheric pressure: 700hPa ~ 1,060hPa;

Storage conditions

Store at room temperature, and storage environment parameters are as follows.

- Ambient temperature: +5°C ~ +30°C ;
- Relative humidity: 10% ~ 80% ;
- Range of atmospheric pressure: 700hPa ~ 1,060hPa;

Transportation conditions

This product should be prevented from pressure, direct sunlight, rain and snow during transportation.

- Ambient temperature: -20°C ~ +60°C ;
- Relative humidity: 5% ~ 80% ;
- Range of atmospheric pressure: 700hPa ~ 1060hPa;

Contact information

Ezisurg Medical Co., Ltd.
399, Miaojiao Rd., Pudong New District, 201315 Shanghai, P.R. China
Tel: +86-21-50456176
Fax: +86-21-50676156

Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffelstrasse 80, 20537 Hamburg, Germany
Tel: +49-40-2513175
Fax: +49-40-255726

After launch of European Database on Medical Devices/EUDAMED, the summary of safety and clinical performance (SSCP) will be made available to the public at <https://ec.europa.eu/tools/eudamed>.

For a patient/ user/ third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/ or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

For additional information, contact your Ezisurg Medical Co., Ltd. Sales Representative or Ezisurg Medical Co., Ltd. Customer Service Tel: +86-4000906176

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Shanghai International Holding Corp. GmbH (Európa)