enScope™

Video Image Processor Instructions for Use

CE 0197

Dedicated To ESP30 Processor

Ezisurg Medical Co., Ltd.

P.2.1013
Ver.B

Rev.B
2025-07-30

Dedicated To ESP30 Processor

en	Instructions for Use	2
fr	Mode D'emploi.....	23
de	Gebrauchsanweisung	45
it	Istruzioni per l'uso.....	67
es	Instrucciones de uso	89
pt	Instruções de Utilização	111
nl	Gebruiksaanwijzing	133
pl	Instrukcja Obsługi	155
cs	Návod k Použití Obrazového.....	177
sv	Bruksanvisning	199
hu	Használati utasítás	221

en

► Intended use/indications for use	3
► Clinical benefit.....	3
► Intended use population.....	3
► Contraindications	3
► Operation qualifications.....	3
► Information signs.....	3
► Warnings & cautions	3
► Adverse events	5
Basic product information	5
► Name & catalogue number	5
► Structure and composition	5
► Software information	5
► Network security.....	5
► Specification.....	7
✓ Environmental requirements	7
✓ Technical specifications	8
Unpacking inspection & appearance & function	8
► Unpacking inspection.....	8
► Product appearance.....	9
► Monitor image	10
Operation & use	14
► 1. Installation and power connection.....	14
► 2. Connection of the processor to the endoscope and the monitor (if needed) ...	14
✓ Connection of the processor to the monitor	14
✓ Connection of the processor to the endoscope.....	14
► 3. Function inspection	15
► 4. White balance	15
► 5. During surgery.....	15
► 6. After surgery.....	16
Reprocessing	16
► Reuse life	16
► Cleaning process	16
After-sales maintenance, troubleshooting, and disposal	17
► Routine maintenance	17
► Troubleshooting	17
► Disposal	18
Essential performance	18
Electromagnetic compatibility	19
► Manufacturer information	22

enScope™
Video Image Processor, Portable

en

► Intended use/indications for use

The video image processor is suitable for use in medical institutions and designed to be connected with the single use digital flexible endoscopes, and other auxiliary equipment for endoscopic diagnosis, treatment and video observation.

► Clinical benefit

Provide clear images for visualization.

► Intended use population

Adults only.

► Contraindications

No known contraindications.

► Operation qualifications

The product should be used only by professional endoscope physicians or medial staff under their supervision. The operators should receive sufficient clinical skill training. But this document does not cover any explanation of or discussion over clinical skills.

► Information signs

The following information signs are used throughout the document:

WARNING : The improper use or misuse of the product may cause damage, death or other serious adverse reactions.

CAUTION : The improper use or misuse of the product may cause problems, such as fault, failure or damage.

► Warnings & cautions

Follow the warnings and cautions below when using the Video Image Processor Additionally, the sections herein provide separate warnings and cautions.

WARNING

- Do not expose the product to rain or damp environment. Do not install or operate the product with wet hands to prevent fire or electric shock.
- Do not connect the high-voltage cable of the peripheral device to the product port via the connector.
- For safety, use compatible cables or accessories that are delivered together with the product, or use cables that are produced by qualified manufacturers per cable specifications. Otherwise, it may cause fire or electric shock.
- Neither pull, bend or wind cables, nor damage the insulated jackets of cables. Otherwise, it may cause fire or electric shock.
- Neither remove screws and housing, nor replace components inside the product.
- Do not block the exit of the product to prevent any fire hazard.
- Stop using the product in case of cable damage or aging to prevent electric shock and other hazards.
- The customer is not allowed to modify the product. The modification may have negative influence on electric safety and electromagnetic radiation. Contact Ezisurg Medical Co., Ltd. for repair and maintenance.
- Do not install, operate or store the product in the following places. Otherwise, it may result in fire, injury or product damage.
 - (1) A place storing flammable chemicals, alcohol, diluent or benzene;
 - (2) A place where there is fluid nearby or the product is prone to exposure to rain;
 - (3) A place with direct sunlight;
 - (4) A places that is close to the air outlet of the air-conditioner or the ventilating device;
 - (5) A place that is close to heat source, such as heating radiator;
 - (6) A place with unstable power supply;
 - (7) A place with high salt content or under sulfur-containing environment;
 - (8) An extremely cold or hot place;
 - (9) An extremely damp or dusty place;
 - (10) A place with frequent mechanical vibration or instability;
 - (11) A place that is close to a strong electromagnetic source, such as TV and radio equipment.

- The power grid connected to the operation room must be grounded reliably. Otherwise, it may hurt medical staff and patients.
 - Do not use the product when it is placed close to, or is used together with, another device, as it may result in improper operation. If the product has to be used, please observe it and other device to check that it works normally.
 - The use of any accessory, sensor and cable other than those specified or provided by the product manufacturer may cause increased electromagnetic radiation or reduced electromagnetic immunity, and may result in improper operation.
 - The portable radio frequency communication device (including antenna cable, external antenna and other peripheral components) should be located at least 30 cm (12 in) from the video image processor (including cable designated by the manufacturer). Otherwise, it may impair product performance.
 - The manufacturer will not assume liability for any interference or response due to the use of cable other than recommended one or the unauthorized change or modification of the product. The unauthorized change or modification of the product may disable the access of the user to operate the product.
 - It is not advisable to use devices emitting radio waves other than those provided by the Company nearby the product, such as mobile phones, radio transceivers, mobile radio transmitters, radio-controlled toys.
 - The product emits some radio frequency energy. it may generate radio frequency interference with other medical, non-medical and radio communication devices. In order to prevent such interference effectively, the product meets radio frequency emission limits under Class A, Group 1 as specified in IEC 60601-1-2. However, the Company will not guarantee that it is free from interference generation in a specific installation environment.
 - Medical devices used together with the product should meet applicable requirements in IEC 60601.
- CAUTION**
- The product should not be used under the following environmental conditions:
 - (1) A place with the explosion hazard;
 - (2) A place with chemicals including combustibles, explosives and volatile solvents (such as alcohol and benzene).
 - Do not store or use the product in a damp room. The humidity for normal operation should be within 30%~75%. Do not expose the product to water drop, splash and spray. Do not place any container with fluid or gas on the electronic device or instrument.
 - The product should be installed, operated, and maintained by, or under the supervision of, the trained and qualified operator.
 - The product can be modified only by the technical personnel or other authorized personnel of Ezisurg Medical Co., Ltd.
 - The maintenance of the product that is allowed by the customer should be conducted following the Instructions for Use. Any additional maintenance can be conducted only by the technical personnel or the authorized personnel of Ezisurg Medical Co., Ltd.
 - The manufacturer is only responsible for the safety, reliability and performance of the product that meets the following requirements:
 - (1) The electrical installation of the product is conducted per the requirements herein;
 - (2) The product is installed, operated and maintained per the requirements herein.
 - Without authorization, it is not allowed to modify the product. Such modification will deprive the user of the warranty right, and the manufacturer will assume no liability for such modification.
 - The product should be used only together with other auxiliary devices that are provided by Ezisurg Medical Co., Ltd. The auxiliary devices of other brands, if used by the user, should be approved by Ezisurg Medical Co., Ltd. or the manufacturer of these devices in terms of their technical and safety availability.
 - The Instructions for Use should be accessible to the operator and the maintainer at any time.
 - Do not inset the power plug into the socket by force. Check whether the plug matches the socket. If the plug or the socket gets damaged, send them to the authorized service provider for repair.
 - Do not pull the power cable forcibly when powering the product off. Instead, pull the power cable while holding the plug.
 - Inspect whether there exists mechanical damage and abnormal function visually before each operation.

en

- In case of any smoke, spark, strange noise or smell, power the product off immediately by shutting off the switch and pulling out the power cable.
- A medical monitor is recommended.
- Before each use, check the compatibility between the endoscope system and any accessory and/or treatment devices according to any standard for safe use specified in the Instructions for Use.

► Adverse events

The Video Image Processor is a product designed to be used together with the Single-use Electronic Endoscope. For more information of adverse events, refer to the Instruction for Use of the Single-use Electronic Endoscope.

Basic product information

► Name & catalogue number

Product name: Video Image Processor

Catalogue number: ESP30

► Structure and composition

The product is composed of a processor and accessories including power cable ,HDIM cable and Power adapter.

► Software information

The software version is controlled by Ezisurg Medical Co., Ltd. The information of the current version, V1, is available in the menu item of "Product Information".

► Network security

Users employ recommended cybersecurity control measures

Item	Description
Equipment repair	Users are not allowed to disassemble the equipment. If there are any issues with the equipment, they must notify the manufacturer or authorized personnel for maintenance.
User authentication and authorization	To ensure data security, users must enter the correct password for authentication before they can view and export data, the initial password is: 666666, and users are allowed to change the password.
Software upgrade	1. The equipment upgrade must be performed by our technical staff or authorized personnel; users cannot upgrade it on their own. 2. When upgrading the software, you need to enter the correct password for identity verification, and the equipment will also verify the integrity and legitimacy of the upgrade package.
Network environment	1. The device does not have a network connection interface, nor does it have a wireless network card interface, and it does not involve network connections. 2. The device has only USB ports for data transfer externally, and we have disabled other USB functions, leaving only the read and write capabilities of the USB ports.

Network ports

The application uses the following network ports for communication with other applications and/or external devices:

Port	Connection	Type	Description
USB port	Local	USB 2.0	Data can be stored and exported through the USB port, we have disabled other USB functions, leaving only the read and write capabilities of the USB ports.

Removable and portable media

When users need to export data using a USB, they must enter the correct password for identity verification.

The data transmit by USB shall be protected, if there is no operation for 30 minutes, will be logged out of the file system.

When using removable media (USB or others,) be aware of these security issues:

- Inserting removable media can introduce a virus to the medical device.
 - Removing media containing patient data can allow access to the data by unauthorized individuals.
- If the media is to be discarded it must be destroyed or disabled so that the data can no longer be accessed.

CAUTION

Whenever media is inserted into the CT, be sure that the media has not been exposed to potential viruses, worms and trojans that infect desktop PCs.

CAUTION

Removable media that contains images and/or other medical information must be stored in a secure area that is not accessible by unauthorized individuals

SBOM download

To enhance the visibility and transparency of the software supply chain, and to provide better risk management and security, we will provide a machine-readable SBOM (Software Bill of Materials) table. The download address for the SBOM table is: <https://www.jianguoyun.com/p/DfxsGhwQqbf-DBi7ytsFIAA>. The SBOM table includes a variety of elements and cycles of the software application, and we will regularly update the SBOM to ensure it remains current and accurate.

Incident response

Item	Description
User authentication and authorization	1. "Please change your password regularly" is displayed on the password login interface to remind users to update their passwords periodically. 2. If the password is entered incorrectly 5 consecutive times, the device authentication will be disabled for 5 minutes, and the disablement information will be recorded and saved in the log. 3. The data transmit by USB shall be protected, if there is no operation for 30 minutes, will be logged out of the file system.

Log file management

The log files that can be used for forensic analysis are stored on the local file system. The deletion/archiving of those log files are handled as part of maintenance activities. Access control for log files is in place. The log files are not recycled. The log files are structured and allow parsingby an automated process (e.g., for IDS).

The log contains the following information:

1. User login information is recorded.
2. File import and export information is recorded.
3. Software upgrade is recorded.
4. User password change information is recorded.
5. Records of account lock due to incorrect password entry by the user.

CAUTION

The viewing and exporting of logs require the entry of the correct password for identity verification. If the correct password is entered and the system is accessed, the device will exit the system files if no operation is performed for 30 minutes.

Cybersecurity use environment

Cybersecurity considerations were made in accordance with the nature of the device, including the device type and the use environments of the system during its lifetime. The following common best practice security controls ("security hygiene") in the use environment are expected to be followed, independently of the status of the system (powered on / powered off):

● Good physical security to prevent unauthorized physical access to the system.

If the system is powered on, follow these additional security controls:

- Network access controls, such as segmentation, to limit medical device communication.
- Malware protection to prevent unauthorized code execution on adjacent devices in the use environment.
- Ensure that prescribed maintenance is done as required, including installation of security patches.
- Security awareness training of the users.

Using the device outside of the intended operating environment may lead to cybersecurity incidents that may affect safety and effectiveness of the device.

Hospital administration is responsible for limiting access to the system and the patient data to appropriate personnel.

Be aware that the screenshot created by the software contains patient-specific data. Saving the screenshot locally thereby saves patient information to the system.

Cybersecurity vulnerability

If you detect a cybersecurity vulnerability or weakness of the device or if a cybersecurity incident takes place involving the device, contact Ezisurg support. Depending on the actual event, Ezisurg supports taking corresponding measures.

Expected service life

From the installation of the device, Ezisurg will provide software services for 5 years. During this period of time, software updates as well as field support are offered. If the software remains in service after the end of support, the cybersecurity risks may increase over time.

When the device exceeds its service life and is scrapped, users must delete the data on the TF card and format the data.

Throughout the software development lifecycle, Ezisurg have regularly released software upgrades and plan to release an upgrade within a year after the last upgrade is released to correct software bugs for testers and customers.

Device software upgrade

The maintenance of the Ezisurg's equipment must be carried out by Ezisurg.

Users can check the software version by opening the device information in the menu.

The manufacturer's system software upgrade is not realized through the network or other means, but through the manufacturer's professional training service personnel to upgrade and maintain the equipment. The technical service personnel manually upgrade the equipment according to the company's upgrade operation guide.

► **Specification**
✓ **Environmental requirements**

Environmental requirements			
	Temperature	Relative humidity	Atmospheric pressure
Operation	10 °C ~ 40 °C	30% RH ~ 75% RH	700hPa ~ 1060 hPa
Transportation	-20 °C ~ 55 °C	10% RH ~ 80% RH	
Storage			

- Do not expose the product to vapor; do not operate the product with wet hands.
- Do not expose the product to direct sunlight.
- Do not place the product in a place with dramatic temperature or humidity changes. Normal operating temperature: 10 °C ~ 40 °C ; humidity: 30% RH ~ 75% RH.
- Keep the product away from any instrument that may generate heat.
- Do not place the product in a place with high humidity or poor ventilation.
- Prevent the product from extreme vibration or sway.
- Keep the product away from chemicals and do not expose the product to inflammable and explosive gas.

- Operate the product carefully and prevent dust (metal in particular) from entering into the product.
- Do not disassemble the product without authorization. Ezisurg will assume no liability for unauthorized disassembly.
- Do not switch off the fan on the side or the back of the product.
- Do not insert the power plug into the socket before the components of the product are connected in place. Otherwise, it may cause a fault.
- Pull out the power cable by holding the plug (instead of the power cable).

✓ Technical specifications <i>Technical specifications</i>		
No.	Item	Specification
1	Product specification/catalogue number	ESP30
2	The power adapter inputs voltage and frequency	AC 100V~240V 50Hz/60Hz
3	Input voltage and frequency	DC 12V
4	Input power	20 VA
5	Type of protection against electric shock	Type I
6	Degree of protection against electric shock	BF (applied part: insert portion of the endoscope)
7	Degree of protection against harmful ingress liquid	General device (enclosed device with no fluid inflow) IPX0
8	Safety degree when used with flammable anesthetic mixture with air or flammable anesthetic mixture with oxygen or nitrous oxide	A device that should not be used with flammable anesthetic mixture with air or flammable anesthetic mixture with oxygen, or nitrous oxide.
9	Working form	Continuous
10	Video output resolution	1280*800P
11	Video output port	1 HDMI output port
12	Service life	5 years

Unpacking inspection & appearance & function

► **Unpacking inspection**

Inspect product components per the packing list after taking delivery. Do not use the product and contact the manufacturer immediately in case of shortage or damage. Please keep all the packing materials that may be required when claiming from the carrier.

The Video Image Processor consists of the following components (subject to the packing list that the customer uses as a reference to take delivery).

Packing list of Video Image Processor

No.	Name	Quantity
1	Video Image Processor	1
2	Power cable	1
3	Power adapter	1
4	HDMI video cable	1
5	Instructions for Use	1
6	Certificate Card	1

► Product appearance

The schematic diagram of the structure is shown in Figure 1 and Figure 2


Figure 1: Schematic diagram of the frontal structure


Figure 2: Schematic diagram of the structure on the back

The structure is described in Table.

Structure description

No.	Name	Description (under standard state unless otherwise specified)
1	Power switch	Short press the power switch key, the medical endoscopic image processor begins to run the work; Long press the power switch key, the medical endoscopic image processor shuts down and stops working.
2	HDMI interface	Port for outputting HDMI video signals.
3	USB interface	Insert a USB flash drive for video and picture storage and export.
4	Power input connector	Use a power adapter to connect to the grid power.
5	Touch screen	Used to display images and medical endoscopic image processor function settings.
6	Signal interface 1	For connecting the electronic endoscope, when using the electronic endoscope, align the signal plug of the electronic endoscope (snap up) with the signal interface 1 of the medical endoscope image processor.
7	Signal interface 2	For attaching the electronic endoscope. When using the electronic endoscope, align the signal plug of the endoscope with the signal interface 2 of the medical endoscope image processor.
8	Dimming button "-"	Used to reduce the brightness of the lighting LED. Users can select "LED0", "LED1", "LED2", "LDE3" and "LED4"
9	Dimming button "+"	Used to increase the brightness of the lighting LED. Users can select "LED0", "LED1", "LED2", "LDE3" and "LED4"
10	Information label	Label that reflect product name, model, specification, batch, category, etc.
11	Device holder	Used to support the device when placing the platform, adjustable.
12	Heat vents	Used for heat dissipation in medical endoscopic image processor

► Monitor image

For the main display interface, see Figure 3.


Figure 3 . Display the main interface

Interface description

No.	Name	Description (under standard state unless otherwise specified)
1	Device manufacturer information	Display the device manufacturer's LOGO and name information
2	Electronic endoscopic information	Displays the name of the inserted endoscope
3	Electronic endoscopic information	Display image zoom level, duration, and other information.
4	Take photos and Video information	Displays information such as the status of the photo and video, and the duration of the recording
5	Store information	Displays information about the remaining memory of the storage device
6	Date, time information	Displays date and time information
7	Warning message	This message will be prompted when 10 minutes remain in the endoscope use.
8	Image display	Display a real-time view of the electronic endoscope
9	WB	White balance the display
10	Photo	The current image displayed by the electronic endoscope is photographed and saved
11	Freeze	The current electronic endoscope display screen is frozen, and the frozen screen is photographed and saved
12	Video	The current image displayed by the electronic endoscope is recorded and saved
13	Switch	Used to switch the main input signal path of the electronic endoscope. For example, the channel A signal is the main input signal, the channel B signal is the secondary input signal, when the "switch" button is pressed, the main signal channel becomes the channel B signal, and the secondary signal channel becomes the channel A signal.
14	Lighting	Used to adjust the LED brightness of the endoscope (5 levels: 0, 1, 2, 3, 4)
15	ADJ	Image adjustment of the picture displayed by the electronic endoscope; The functions of adjustment include brightness, contrast, structural enhancement, electronic magnification, noise reduction and other functions.
16	Settings	Setting up the image processor; Configurable functions include playback, controller button customization, output settings, language settings, time settings, factory reset, device information, etc

Figure 4 shows the image adjustment display interface


Figure 4: Image adjustment display interface

Description and operation instructions of image adjustment functions

No.	Menu tab	Feature description	Operating instructions
1	Close	Used to turn off the image feature setting	Click "x" to close
2	Brightness	Adjust the display brightness of the image processor	Click the touch screen, pop-up menu interface→ click Image Adjustment→ select the appropriate gear (4 gears: 1, 2, 3, 4) under the brightness function
3	Contrast	Adjust the contrast of the displayed image	Click the touch screen, pop-up menu interface→ click Image Adjustment→ select the appropriate gear (4 gears: 1, 2, 3, 4) under the contrast function
4	Enhancement	Structural enhancements to the displayed image	Click the touch screen, pop-up menu interface→ click Image Adjustment→ select the appropriate gear (4 gears: 1, 2, 3, 4) under the structure enhancement function
5	Zoom	Electronic enlargement of the displayed image, support image enlargement, and zoom out to the original image after zooming in	Click the touch screen, pop up menu interface→ click Image Adjustment→ select the appropriate gear (3 gears: x0, x1, x2) under the electronic magnification function
6	Noise Reduction	Noise reduction adjustment of the displayed image	Click the touch screen, pop up menu interface→ click Image Adjustment→ select the appropriate gear (2 gears: OFF, ON) under the noise reduction function

Figure 5 shows the display screen


Figure 5: Setting up the display interface

Description of setting features and their corresponding actions

No.	Menu tab	Feature description	Operating instructions
1	Close	Close the settings menu	Click "x" to close
2	View	View and export saved images and videos	Click the touch screen, the menu interface will pop up → click Settings→ click Playback→ select the corresponding folder→ select the desired picture or video for viewing and exporting
3	Handle Button	The function of the electronic endoscope handle button (if available) can be customized	Click the touch screen, pop-up menu interface→ click Settings→ click the Customized Key→ select the corresponding key change function (button 1 defaults to the electronic magnification function; Optional photo, video, electronic magnification, freezing function; Button 2 defaults to the photo function; Optional photo, video, electronic magnification, freeze function)
4	Output Setting	Set the border format for image display	Click the touch screen, pop-up menu interface→ click Settings→ click Output Settings→ select the displayed screen format
5	Language Setting	Set the language in which the image is displayed	Click the touch screen, pop-up menu interface→click Settings→click Language Settings→ select the desired language
6	Time Setting	The time and date can be set for the image processor	Click the touch screen, pop-up menu interface → click Settings→click Time Setting→ enter the time that needs to be adjusted, and click Confirm to finish
7	Factory Settings	All options in the menu can be factory reset	Click the touch screen, pop-up menu interface → click Settings→ click Factory Settings→ click Confirm
8	Device Info	You can view the device information and software version of the image processor	Click the touch screen, pop up the menu interface→ click Settings→click Device Information→ view device information
9	Format	Erase all data and restore storage space.	Tap the touchscreen to open menu → Tap Settings → Tap Format → Erase all data

Operation & use

en

► 1. Installation and power connection

- (1) Install the product on a stable plane.
- (2) Connect the power cable attached to the product to the main socket with the protective grounding terminal only after it is confirmed that the power switch is in "OFF" state.

CAUTION

- Do not use the product near the high-power device. Otherwise, the product is prone to image output disturbance.
- Arrange the processor in a ventilated place.
- If the Video Image Processor is placed on the cart, the area available to place the processor on the cart should be equal to or larger than the base area of the processor.
- Inspect whether the power cable matches the power source before turning the mains switch on.
- Connect the product to the preset possible equipotential terminal and keep the product grounded completely after inspecting the power supply.
- Insert the power plug into the socket completely to prevent the plug from being pulled out accidentally during the surgery.
- Use the power cable and power adapter provided by Ezisurg , and do not modify the power cable.
- Keep the power socket which the processor is connected to grounded reliably. When turning the power supply off, do not pull the cable by force. Instead, pull out the power cable by holding the plug.

► 2. Connection of the processor to the endoscope and the monitor (if needed)

✓ Connection of the processor to the monitor

The Video Image Processor comes with a display for image display, and supports the use of HDMI video cable to connect the video output of the Video Image Processor to the video input of the monitor. The monitor must meet the IEC60601 standard.

✓ Connection of the processor to the endoscope

When using an electronic endoscope that matches Signal Interface 1, ensure that the Video Image Processor is powered off, then connect the endoscope to the Video Image Processor and turn on the processor. When attaching the endoscope to the processor, align the bulge on the endoscope signal plug with the notch on the processor signal socket and insert it for connection (Figure 6).


Figure 6. Connect the electronic endoscope signal plug to the Video Image Processor

WARNING

- The convex cannot be inserted into the concave if they are not aligned. Insertion by force may damage the endoscope or the processor.
 - Use the electronic endoscope of the specified model. Otherwise, it may damage the product.
- When using an electronic endoscope that matches signal interface 2, ensure that the Video Image Processor is powered off, then connect the endoscope to the Video Image Processor and turn on the processor. When attaching an endoscope to the processor, snap the endoscope signal board up and insert it to make the connection (Figure 7).


Figure 7. Connect the endoscope signal plug to the Video Image Processor

WARNING

- Insertion will fail if the direction is wrong. Insertion by force may damage the endoscope or the processor.
- Use the electronic endoscope of the specified model. Otherwise, it may damage the product.

► **3. Function inspection**

After proper installation, inspect the following items:

- (1) Inspect whether the product is connected properly;
- (2) Short press the power button and wait for the device to turn on, click the switch button to switch the signal to the desired signal channel, and inspect whether there is any image on the monitor.
- (3) Please check the internal storage space and make sure it is adequate.

CAUTION

- Do not replace the endoscope when the Video Image Processor keeps powered on. Neither remove nor insert the signal plug of the endoscope when the processor is working.
- To ensure smooth surgery, inspect whether the product has mechanical damage or abnormal function visually before each operation, and have a standby device in the operating room as far as possible.
- Switch to the signal input source matching the monitor. Otherwise, there may be no output on the external monitor.
- (4) Inspect whether main buttons function normally.

CAUTION

- Do not expose the tip portion of the endoscope to high light during operation. Otherwise, it may impair image quality.

► **4. White balance**

White balance adjustment: Keep the tip portion of the endoscope facing a white object (white gauze recommended) at a distance of 1-2 cm. Press the "White Balance" button, leaving the monitor to display "WB" in the interface. Keep the head end still until the monitor screen displays "WB OK ". Then move the tip portion away from the white object to complete white balance adjustment.

CAUTION

- Re-adjust white balance if the image color is unsatisfactory after each startup.
- If color quality is unsatisfactory, inspect the monitor first (chroma, contrast, saturation, etc.). For some monitor models, converters for color temperature adjustment are provided on their rear panels to facilitate users to select favorite colors.

WARNING

- Do not look straight at the luminous tip portion of the endoscope with naked eyes as such action may result in eye discomfort or damage.

► **5. During surgery**

Capture the proper image by operating keys and menu items as mentioned in "Unpacking Inspection, Appearance and Function".

WARNING

- Ensure the patient has no access to the Video Image Processor in work.
- Inspect whether the applied parts of the Video Image Processor have rough surfaces, sharp edges and protrusions that may result in safety hazards before each use.

- If it is suspected or confirmed that the product develops a fault, contact Ezisurg Medical Co., Ltd. or its designated service provider for inspection before use.
- Before each use or change of view mode/setting, the operator should check that real-time images (instead of stored images) with correct image directions are viewed via the Video Image Processor.
- Do not operate the endoscope and its accessories when playing the video.
- According to interconnection requirements, Type BF applied parts should be used as the applied parts of the other ME device in the control unit. And the applied parts of some CF devices cannot be used together. Otherwise, the leak current in the applied part of the CF device may exceed limits.
- Inspect the accessories of the product before each use.
- The operator and his/her assistant should wear rubber gloves.

CAUTION

- If photo and video functions are required during operation, ensure the hard disk drive has sufficient remaining memory capacity. If the memory capacity is used up, insufficient memory space will be prompted.

► **6. After surgery**

- (1) Turn off the power switch, and disconnect the power cable from the socket by pulling out the power plug.
- (2) Remove the electronic endoscope from the processor.
- (3) Clean and sterilize the product housing:
 - Wipe off dust on the product surface gently with soft dry cloth first. Then, clean off stubborn dirt with the cloth that is dipped with neutral detergent and wipe up the surface with dry cloth.
 - Clean up the product surface with soft cloth or operating gauze that is dipped with 75% medical alcohol twice, which lasts for 3 min.

WARNING

- Check that the power cable is pulled out and the power switch is turned off before cleaning the processor.

CAUTION

- Do not use abrasive detergent or dissolving agent. Otherwise, it may damage the surface of the product or its label.
- Note that excessive detergent may permeate into the product and consequently, damage the electronic components. Do not immerse the product into liquid.
- Only when the cleaned product and its components are dried completely can they be powered on.
- Clean, inspect and maintain the product following the Instructions for Use every other week. (For periodic maintenance, see the chapter below.)

Reprocessing

► **Reuse life**

The service life of the device is 5 years, the calculation method starts from the completion of the equipment installation.

Perform the visual and functional inspection before each reuse to ascertain whether the device has exceeded its use life.

1. All connecting cables are in good condition and without damage.
2. The signal input interface and output interface are intact and free of dirt.
3. The housing has no unacceptable deterioration such as corrosion, discoloration, pitting, cracked seals.
4. The device can be turned on normally and all buttons can be operated normally.

► **Cleaning process**

1. Wipe off dust on the product surface gently with soft dry cloth first. Then, clean off stubborn dirt with the cloth that is dipped with neutral detergent and wipe up the surface with dry cloth.
2. Clean up the product surface with soft cloth or operating gauze that is dipped with 75% medical alcohol twice, which lasts for 3 min until the surface of the device is visually clean.
3. In case of visual inspection of residual dirt and if residual dirt is visible on the equipment, repeat the above steps with a new soft cloth soaked in above-mentioned disinfectant.

After-sales maintenance, troubleshooting, and disposal

► Routine maintenance

Perform routine maintenance every other week:
(1) Make sure the power cable is intact and not damaged.
(2) Clear memory timely or regularly as insufficient memory space may disable photo and video functions. The system only supports U disks with a capacity of no more than 128G. If the U disk cannot be read by the system, the capacity information will not be displayed on the main interface, please replace it with another available U disk.

CAUTION

- The maintenance of the product should be conducted following the Instructions for Use. Any additional maintenance can be conducted only by the technical personnel or the authorized personnel of Ezisurg Medical Co., Ltd.
- Contact Ezisurg Medical Co., Ltd. for after-sales repair and maintenance.
- Address: No. 399, Miaojiao Rd., Pudong New District, 201315 Shanghai, P.R. China
- Tel: +0086-4000906176

► Troubleshooting

Inspect the power cable and other connecting line first in case of any fault. Consult the checklist below before contacting the service provider. It is suggested to have a standby device during the surgery. If the fault cannot be solved, shut down the failed device, pull out its power cable, and connect the standby device to continue the surgery. In case of fault removal failure, contact the authorized representative of Ezisurg Medical Co., Ltd.. For common faults and troubleshooting methods, refer to Table.

Common faults and troubleshooting methods

Fault	Possible cause	Troubleshooting
The indicator does not light up after the power switch is turned on.	· The power cord is not connected, · The image processor is damaged.	· Check that the power cable is connected · Contact the manufacturer
Unclear output image	· The exposed surface of the electronic lens gets contaminated · The endoscope gets damaged	· Clean the exposed surface of the electronic lens with the alcohol cotton · Contact the manufacturer to replace the endoscope
Dim output image	· The brightness of the illuminator for the endoscope is not adjusted properly · The brightness of the monitor is not adjusted properly	· Adjust LED brightness of the Video Image Processor by operating the corresponding button · Adjust the brightness of the monitor to a proper scope
There is an image, but photo and video functions are unavailable.	· The memory capacity is used up. · Functional failure	· Delete some hard disk storage files; If you cannot delete it, please contact the manufacturer · If the function fails, please contact the manufacturer
No output of the external monitor	· The power supply of the external monitor is not connected, or the power switch is not turned on · The power cable of the processor is not connected to the socket, or the power switch is not turned on · The HDMI cable is not connected properly	· Check that the power supply of the external monitor is connected and that the power switch is turned on · Check that the power cable is connected to the socket and that the power switch is turned on · Re-connect the HDMI cable properly

Fault	Possible cause	Troubleshooting
No output of the external monitor	· HDMI cable fault	· Replace the HDMI cable
	· The signal plug of the endoscope is not connected in place	· Pull out the plug and re-insert it in place
	· The connecting channel of the monitor is wrong	· Connect the monitor via the correct channel (Some monitor models have separate output and input ports)
	Wrong signal input source is selected for the monitor	Select the correct signal input source matching the monitor
Abnormal image color	· Abnormal monitor color setting	· Reset the monitor color
	· Abnormal color setting of the processor	· Reset the processor color
	· HDMI cable fault	· Replace the HDMI cable
	· Absence of white balance (WB) operation	· Conduct white balance (WB) operation

If the troubleshooting methods do not work when the Video Image Processor develops a fault, contact the authorized representative of Ezisurg Medical Co., Ltd. and follow the rules below:
(1) Cleaning and disinfection the failed product before packaging.
(2) As the product is a precise instrument, use rigid packaging with anti-shock measures to prevent any damage during transportation.
(3) Detail the fault in writing. For any other requirements, contact the authorized representative of Ezisurg Medical Co., Ltd.

► Disposal

Dispose of the product, its accessories and packaging material according to the regulations of the hospital, the administration and/or the local government.
When the device is scrapped, users must delete the data on the TF card and format the data.

Essential performance

1. Be able to display real-time images.
2. Be able to display image with the correct orientation.
3. There is no noise in image which can affect doctor's diagnosis, noise may manifest as grains or specks, horizontal or vertical stripes, flickering, and snowy screen.

Electromagnetic compatibility

CAUTION

- For this equipment, special precautions regarding electromagnetic compatibility (EMC) must be taken, and it must be installed and used in accordance with the electromagnetic compatibility information specified in this manual.
- This product conforms to the relevant requirements of IEC 60601-1-2 standard electromagnetic compatibility.
- Portable and mobile radio frequency communication devices may hinder the normal operation of the Video Image Processor or the electronic endoscope. Please use the processor and the endoscope under the recommended electromagnetic environment.
- The use of accessories and cables other than those sold as spare parts for the product by the manufacturer may increase electromagnetic emission or lower electromagnetic immunity.
- The following cable types (included but not limited to) should be used to meet disturbing radiation and immunity standards.

Cable information

Cable name	Length
Power cable	2.0 m
HDMI cable	1.5 m

- The Video Image Processor and the electronic endoscope should not be used when they are placed close to, or stacked with, another device. If they have to be used, observe to check that they work normally in the given configuration.
- WARNING
- Devices or systems should not be close to or stacked on top of other devices. If they must be close to or stacked on top of other devices, they should be observed and verified to operate normally in the configuration in which they are used.
 - Class A equipment intended for use in industrial environments may have potential difficulties in ensuring EMC in other environments due to product conduction and radiation harassment.
 - Except for cables sold by the manufacturer as spare parts for internal components of the product, accessories and cables other than those specified for use may result in increased product emission or reduced disturbance immunity.
 - Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the product, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result in .

Manufacturer's Declaration – Electromagnetic Emissions		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions (CISPR 11)	Group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic environment.
RF emissions (IEC/CISPR 11)	Class A	The equipment is suitable for use in all establishments, other than non-domestic or domestic establishments directly connected to the public low-voltage power supply network that supplies buildings used.
Harmonic emissions (IEC 61000-3-2)	Not Applicable	
Voltage fluctuations/ flicker emissions (IEC 61000-3-3)	Not Applicable	

Manufacturer's Declaration – Electromagnetic Immunity		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Immunity test	Compliance Level	Electromagnetic Environment – Guidance
Electromagnetic discharge(ESD) (IEC61000-4-2)	± 8 kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV Air	Floors should be made of wood, concrete or ceramic tile that hardly produces static. If floors are covered with synthetic material that tends to produce static, the relative humidity should be at least 30%.
Electromagnetic fast transient/ burst (IEC61000-4-4)	± 2 kV for power supply lines ± 1 kV for signal lines	Mains power quality should be of a typical commercial (original condition feeding the facilities) or hospital environment.
Surge (IEC61000-4-5)	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines (IEC61000-4-11)	0% U _r , 0.5cycle (At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°) 0% U _r , 1 cycle 70% U _r , 25/30 cycles Single phase: at 0° 0% U _r , 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment with highly reliable service. If the user of the equipment requires continued operation during power mains interruptions, it is recommended that the equipment be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field (IEC 61000-4-8)	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Proximity fields From RF wireless communications equipment (IEC 61000-4-3)	Refer to next table	Portable and mobile RF communications equipment should be used no closer to any part of the product
Proximity magnetic fields (IEC 61000-4-39)	134.2 kHz, PM, 2.1 kHz, 65 A/m 13.56 MHz, PM, 50 kHz, 7.5 A/m	The sources of proximity magnetic fields should be used no closer than 0.15 m to any part of the product
Conducted RF (IEC 61000-4-6)	3 Vrms, 150 kHz ~80 MHz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2\sqrt{P}$
Radiated RF (IEC 61000-4-3)	3 V/m 80 MHz ~2.7 GHz	$d=2.3\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey(a), should be less than the compliance level in each frequency range(b). Interference may occur in the vicinity of equipment marked with the following symbol: 
Note : 1: U_i is the AC mains voltage before application of the test level. 2: At 80 MHz and 800 MHz, the higher frequency range applies. 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. 4: To established the proximity threshold 0.15 m for Proximity magnetic fields, the IEC Subcommittee (SC) 62A considered the types of proximity magnetic field disturbance sources expected: – induction cooking appliances and ovens operating at frequencies up to 30 kHz; – RFID readers operating at both 134.2 kHz and 13.56 MHz; – electronic article surveillance (EAS) systems; – sponge detection systems; – equipment used for position detection (e.g. in catheter labs); – wireless power transfer charging systems for electrical vehicles that operate in the frequency range of 80 kHz to 90 kHz. These frequencies and applications are representative examples based on sources of magnetic field disturbance in use at the time of publication of the collateral standard IEC 60601-1-2:2014+A1:2020 . a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the equipment. b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Recommended minimum separation distances:

Nowadays, many RF wireless equipment have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. This Systems has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014+A1:2020 . The customer and/or user should help keep a minimum distance between RF wireless communications equipment and this Systems as recommended below.

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	430-470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9
5500						
5785						

Note: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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Shanghai International Holding Corp. GmbH (Európa)