
easyUS™ Ultrasonic Surgical System Generator User Manual

REF

ES01

(for use with software version V1)

Ezisurg Medical Co., Ltd.

Please read all information carefully before use.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Important!

Please read all information carefully. Failure to properly follow the instructions may lead to serious surgical consequences.

1. This manual is designed to provide instructions for use of the generator ES01 and its accessory footswitch FSU01. It is not a reference to surgical techniques.
2. The ES01 generator should be used with easyUS™ handpiece eHP01 and eHP02 and easyUS™ shears SHC series and SHF series manufactured by Ezisurg Medical. Please refer easyUS™ Single Use Ultrasonic Shears Instructions for Use, easyUS™ Single Use Ultrasonic Scissor-type Shears Instructions for Use and easyUS™ Ultrasonic Surgical System Handpiece Instructions for Use for more information.
3. This manual contains important information for operation of ES01. It should be kept where it may be referenced during usage. Please print or copy pages as necessary to keep nearby.
4. Ezisurg Medical Co., Ltd reserves the right to change the electrical/electronic or mechanical configurations and components without notice. Schematic revisions may not match PCB revisions in the unit being maintained; consult Ezisurg Medical Co., Ltd before making any changes.
5. This user manual is only for software version 1 of the generator ES01. Software version other than version 1 may not match the diagram in the manual. Please examine the software version of the generator ES01 before starting maintenance; consult Ezisurg Medical Co., Ltd for the proper service manual for the using software version.
6. easyUS™ is the registered trademark of this product.

WARNING: No modification of this equipment is allowed.

System Description

The Ultrasonic Surgical System converts high frequency current into ultrasonic oscillation to mechanically fragment soft-tissue cells upon contact with a vibrating tip for cutting and coagulating tissue during surgery. It mainly consists of a generator (with software embedded), a footswitch, a handpiece and an easyUS™ Ultrasonic shears. The shears and the handpiece have a torque wrench and/or a test tip as spare parts. The torque wrench is to secure the test tip or shears to the handpiece. The test tip is to connect the handpiece and is used to help to calibrate and diagnose the handpiece performance. Among them, the generator outputs electrical energy (controlled by the software); the

footswitch is used to active MAX or MIN output; the handpiece converts the electrical energy to mechanical motion; the shears connects to the handpiece and pass the mechanical motion to the distal end of the blade in the shears.

The Ultrasonic Surgical System does not include an integrated irrigation/suction system.
The Ultrasonic Surgical System does not support the function of patient physiological parameter monitoring.

The shaft of the easyUS™ single use Ultrasonic shears, including blade, clamp, tissue pad, trigger tube and outer tube are applied parts.

See Figure 1-1 for the diagram of Ultrasonic Surgical System.


Figure 1-1 Diagram of Ultrasonic Surgical System

- 1、 generator (including software) 2、 footswitch 3、 handpiece
4、 single use ultrasonic shears

See Table 1-1 for the specifications of Ultrasonic Surgical System.

Table 1 -1 Specifications of Ultrasonic Surgical System

Category	Product Code	Description
Generator	ES01	
Handpiece	eHP01	
Shears	SHC14A	14cm shaft, curved blade
Shears	SHC23A	23cm shaft, curved blade
Shears	SHC36A	36cm shaft, curved blade
Shears	SHC45A	45cm shaft, curved blade
Footswitch	FSU01	

Category	Product Code	Description
Test Tip	TP01	Spare part of handpiece
Torque Wrench	TW01	Spare part of handpiece and shears

See Figure 1-2 for the diagram of Ultrasonic Surgical System.


Figure 1-2 Diagram of Ultrasonic Surgical System

1、 generator (including software) 2、 footswitch 3、 handpiece 4、 single use ultrasonic scissor-type shears

See Table 1-2 for the specifications of Ultrasonic Surgical System.

Table 1-2 Specifications of Ultrasonic Surgical System

Category	Product Code	Description
Generator	ES01	
Handpiece	eHP02	
Shears	SHF9A	9cm length, curved blade, handswitch color is black
Shears	SHF17A	17cm length, curved blade, handswitch color is black
Footswitch	FSU01	
Test Tip	TP02	Spare part of handpiece eHP02
Torque Wrench	TW03-01	Spare part of handpiece eHP02 and shears SHF9A
Torque Wrench	TW03-02	Spare part of handpiece eHP02 and shears SHF17A

Generator

Generator is to generate high frequency alternating current to drive the handpiece. In addition, generator provides Human Machine Interface to have user control, monitor, and diagnose the status of Ultrasonic Surgical System.

Handpiece

Handpiece is to connect a generator and a shears. There is an ultrasonic transducer inside the handpiece to convert electrical energy to mechanical energy. The handpiece has a cable connects to a generator. Then the handpiece converts high frequency alternating current to high frequency mechanical motion, and transmit the motion to shears.

A handpiece eHP01 is packaged with a torque wrench and a test tip. A handpiece eHP02 is packaged with a test tip. The handpiece, the torque wrench and the test tip are packed non-sterile when manufactured. The handpiece must be cleaned and sterilized prior to each use. Please refer *easyUS™ Ultrasonic Surgical System Handpiece Instructions for Use* for details.

Shears

easyUS™ Ultrasonic shears is handled by surgeons. The shears has multiple functions like clasp, dissection, cutting and coagulating. The proximal end of the blade in shears receives high frequency mechanical motion and transmits the motion to the distal end of the blade.

A shears is packed with a torque wrench in sterilized package and sterilized using ethylene oxide when manufactured. Please refer *easyUS™ Single Use Ultrasonic Shears Instructions for Use* for details.

Footswitch

Footswitch is the foot control part of Ultrasonic Surgical System. It includes two pedals, active MIN and MAX output accordingly.

Torque Wrench

Torque wrench is the spare part of shears and handpieces. It is used to secure test tips with handpiece or shears with handpiece in appropriate torque.

When torque wrench is packed with easyUS™ Ultrasonic shears it is sterilized using ethylene oxide. When torque wrench is packed with handpiece it is non-sterilized. Please refer *easyUS™ Single Use Ultrasonic Shears Instructions for Use* and *easyUS™ Ultrasonic Surgical System Handpiece Instructions for Use* for details.

Test Tip

Test tip is the spare part of handpiece. It is used to help to diagnose the handpiece performance.

Test tip is packed with handpiece and non-sterilized. Please refer *easyUS™ Ultrasonic Surgical System Handpiece Instructions for Use* for details.

Theory of Operation

easyUS™ Ultrasonic Surgical System is for hemostatic cutting and/or coagulation. When activated, the blade in shears vibrates longitudinally at 55.5 kHz. This ultrasonic vibration is transmitted to the distal end of the blade, causing the blade to vibrate and generate heat. The coupling causes collagen molecules within the tissue to vibrate and become denatured, forming an irreversible cross-linking, hence occludes the blood vessels for a permanent fusion. At the same time, the cavitation effect caused by high-frequency vibration accelerates the vaporization of heated water in the tissue cells in contact with blade, further breaks the protein hydrogen bond, and disintegrates the cell, so the tissue is divided.

Intended Purpose

The instrument is intended to seal and cut vessels and/or to dissect, coagulate and cut soft tissues via ultrasonic oscillation.

Device Performance

The Generator ES01 supplies energy to the handpiece eHP01 or eHP02. The Generator ES01 uses a touchscreen display interface and has a unique receptacle port that accepts the handpiece eHP01 or eHP02.

The Generator ES01 is not designed to operate unauthorized devices. Such use is not in accordance with the design & use parameters, instructions and guidelines for the product.

Indications for Use

easyUS™ Ultrasonic Surgical System is indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The Ultrasonic Surgical System instruments can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels. When used with the compatible instruments, the Ultrasonic Surgical System can be used to coagulate isolated vessels up to and including 5 mm in diameter.

Contraindications

1. The instruments are not indicated for incising bone.
2. The instruments are not intended for contraceptive tubal occlusion.
3. The Ultrasonic Surgical Systems instruments are contraindicated when, in the judgment of the physician, ultrasonic surgery would be contrary to the best interest of the patient.

Preparation at the Point of Use

Turn off the generator ES01 before processing. Remove gross soiling on the device with lint free towel wetted in cold water (<40°C) immediately after use. Don't use a fixating detergent or hot water (>40°C) as this can cause the fixation of residuals which may influence the result of the reprocessing process.

Transportation

If applicable, safely transport the device to the reprocessing area to avoid any damage and contamination to the environment.

Target Patient Group

The intended patient population consists of adult and pediatric patients.

Intended Users

Intended users are medical professionals qualified for the ultrasonic surgical technique and surgical procedure to be performed.

Clinical Benefits

Support for successful cutting and coagulation.

Undesirable Side Effects / Residual risk

Undesirable side effects and risks associated with ultrasonic devices include the potential for bleeding, tissue injury via mechanical or thermal damage, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, electrical shock, foreign body or magnetic resonance incompatibility, and property or environmental damage. Also, unintended harm, extended surgery, or altered surgical approach may result from issues related to device activation, damaged devices, electromagnetic interference, audible noise due to misassembly, misuse of the torque wrench, or an attempt to alter the device.

Instructions for Reprocessing

The generator ES01 can be cleaned and disinfected by wiping the device with disinfectant. An automated reprocessing or immersion cleaning and disinfection is not possible. For this Instructions, routine decontamination should be performed on a daily basis.

Reprocessing procedures have only limited implications to the ES01. The limitation of the numbers of reprocessing procedures is therefore determined by the function / wear of the ES01. From the processing side, there is no maximum number of allowable reprocessing. The ES01 should no longer be reused in case of dysfunction and/or signs of material

degradation, e.g. damaged insulation or other change in the material, etc.

Overview

Device Description

The generator ES01 supplies energy to the handpiece. The handpiece converts electrical energy to mechanical motion, then transmits the motion to distal end of the blade in the shears to cut/coagulate soft tissue. ES01 uses a touchscreen display interface and has a unique receptacle port that accepts the handpiece.

How Supplied

ES01 is supplied in a semi-ready-to-use state. The shipping box contains the generator ES01, power cord and user manual. The Ultrasonic Surgical System instruments like handpieces and shears are not included in this packaging and must be purchased separately.

Illustration and Nomenclature

Front Panel of the Generator


Figure 2 Front Panel of ES01

1. Standby touch key	Hold Standby touch key to switch between Standby Mode and Ready Mode. The key shows amber under Standby mode and green under Ready Mode.
2. Error light	Normally off. Show red when system has fatal error.
3. Display/Touch screen	Displays system information and serves as interface for adjusting controls and settings.
4. Handpiece receptacle	Receptacle used to attach the handpiece connector.

Back Panel of the Generator


Figure 3 Back Panel of ES01

5. Air vents	For the cooling fan.
6. 1 st Footswitch receptacle	Receptacle used to connect the 1 st footswitch to the generator.
7. 2 nd Footswitch receptacle	Receptacle used to connect the 2 nd footswitch to the generator.
8. On/Off switch	Power on/off the generator.
9. Potential equalization terminal	Provides means for connection to a potential equalization conductor.
10. Fuse holder	Use to hold fuse for the generator.
11. Power cord receptacle	Receptacle used to attach the power cord to the generator.

Footswitch

Figure 4 Footswitch

12. MIN (left pedal)	Activate generator for the minimum output.
13. MAX (right pedal)	Activate generator for the maximum output.

Useful Life

The useful life of generator ES01 depends on parts and accessories connects the generator if ES01 is under normal use and preventive maintenance. If parts or accessories failed or are beyond their useful life there needs to replace them to continue use.

Customer Service

Warranty

This warranty and the rights and obligations hereunder shall be construed under and governed by the laws of China.

Ezisurg Medical Co., Ltd warrants this product to be free from defects in material and workmanship under normal use and preventive maintenance for the respective warranty period shown below. Ezisurg Medical Co., Ltd's obligation under this warranty is limited to

the repair or replacement, at its option, of any product, or part thereof, which has been returned to Ezisurg Medical Co., Ltd or its distributor within the applicable time period shown below and which examination disclosed, to Ezisurg Medical Co., Ltd's satisfaction, to be defective.

This warranty does not apply to any product, or part thereof, that has been:

- 1) Adversely affected due to use with devices manufactured or distributed by parties not authorized by Ezisurg Medical Co., Ltd;
- 2) Repaired or altered outside Ezisurg Medical Co., Ltd's factory in a way so as to, in Ezisurg Medical Co., Ltd's judgement, affect its stability or reliability;
- 3) Subjected to improper use, negligence or accident;
- 4) Used other than in accordance with the design and use parameters, instructions and guidelines for the product or with functional, operational or environmental standards for similar products generally accepted in the industry.

Ezisurg Medical Co., Ltd's products are warranted for the following periods after delivery to the original purchaser:

Generator ES01 (including power cord)	One (1) year, parts and labor
Footswitch FSU01	One (1) year, parts and labor

Unless superceded by applicable local law, this warranty is in lieu of all other warranties, express or implied, including the warranties of merchantability and fitness for a particular purpose, and of all other obligations or liabilities on the part of ezisurg medical co., ltd and is a purchaser's exclusive remedy. In no event shall ezisurg medical co., ltd be liable for special, incidental or consequential damages including, without limitation, damages resulting from loss of use, profits, business or goodwill, other than as expressly provided by a specific law. Ezisurg medical co., ltd neither assumes nor authorizes any other person to assume for it any other liability in connection with the sale or use of any of Ezisurg Medical Co., Ltd's products. There are no warranties that extend beyond the terms hereof. Ezisurg Medical Co., Ltd reserves the right to make changes to products built and/or sold by them at any time without incurring any obligation to make the same or similar changes on products previously built and/or sold by them.

Theory of Operation of System

In Ultrasonic Surgical System, the generator ES01 generator electrical signals, which is passed to handpiece eHP01 or eHP02 via the handpiece cable, then the handpiece converts the electrical energy to mechanical motions. The shears receives the mechanical motions from handpiece, amplifies the motions and delivers the energy to tissues.

Generator

The generator converts the AC line voltage to a controlled DC level. The DC level is then modulated at the resonant frequency of the hand piece. The modulated signal is then filtered and delivered to the hand piece, where it resonates the acoustic drive train. Below is the basic function diagram of the generator:


Figure 5 The Basic Functional Diagram of Generator ES01

Power Entry Module

Consist of switch, power receptacle, filter, fuse, AC-DC power. It is to provide 48 VDC to the whole system.

Power Regulation Circuit

Converts 48 VDC to different types of powers used in the generator.

DSP

DSP contains software for running. The software drives the display and touchscreen component to control HMI as well as to drive FPGA to control peripheral devices like to adjust output frequency and power dynamically, to play voice and other components.

FPGA

FPGA is an auxiliary component in circuit.

Frequency Synthesis Circuit

Generate high frequency AC current. The frequency is dynamically adjustable along with the changes of the load to resonate the whole system.

Power Adjust Circuit

Adjust output current and voltage according to feedback impedance to have system work on its best status.

Power Amplify Circuit

Amplify the high frequency AC current via the transformer so that the current is to drive ultrasonic the transducer in handpiece.

Patient Isolation Circuit

Provide a safety isolation barrier to patient and/or user. The voltage generated on the secondary of the transformer is isolated from the primary and thereby not referenced to earth ground.

Display and Touchscreen Component

The display and touchscreen component communicates with DSP and draws the screen as DSP requires as well as passes user touch action to DSP.

Touch Key Component

The touch key component communicates with FPGA. It senses user actions on the touch key and illuminate different color LEDs according to current status.

Handpiece

The handpiece houses several major components that generate, amplify, and deliver ultrasonic energy to the instrument end-effector. When each component is attached to the other and tuned, an acoustic drive train is formed. The ultrasonic transducer in handpiece converts electrical energy to mechanical motions. When an AC waveform is applied to the transducer, the piezoelectric material expands and contracts to produce longitudinal motion.

Shears

The blade in the shears, as known as the wave guide, is to amplify the mechanical motions and deliver the energy to tissues. The distal end of the blade is in contact with human tissues for cutting and/or coagulation. There is handswitch on the shears. User may use handswitch. The action of handswitch is transmitted to DSP in generator. The DSP controls minimum power output or maximum power output according to the handswitch pressed.

Footswitch

User may use footswitch. The action of footswitch is transmitted to DSP in generator. The DSP controls minimum power output or maximum power output according to the pedal stepped down.

Instructions for Use

Setup

WARNING: To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

Caution: Do not block the air vents of the generator to avoid generator overheating.

1. Examine the generator ES01 with Ultrasonic Surgical System handpiece and shears for damage. Do not use damaged devices and instruments.
2. Secure the generator on a cart or any other suitable fixture in the appropriate position.
3. Connect the power cord to the power cord receptacle on the back panel of the generator and to a grounded electrical outlet. The power requirements for the generator are listed on the label on its back panel.

4. Connect the footswitch connector to the footswitch receptacle on the back panel of the generator if applicable.
5. Secure the Ultrasonic Surgical System handpiece and shears by the torque wrench. Connect the handpiece connector to the handpiece receptacle on the front panel of the generator.

Caution: Verify that the generator is secured to a cart or any other suitable fixture before powering up.

6. Turn on the generator using the On/Off switch on the back panel of the generator.
7. All lights on the front panel illuminates and the system will run its initiation sequence. A tone is played while an animation is shown on the screen during the initiation.
8. After the initiation sequence is done the system enters Standby mode. User may replace handpieces and/or shears under Standby mode.
9. Hold the Standby touch key for 3 seconds and the test prompt screen appears. Follow the instructions on the screen to run the Pre-Run test.
10. The ready screen appears when the Pre-Run test is finished and the generator is ready for use.

Ready Screen

The ready screen has the following features:

Figure 6 Ready Screen

1. MIN	Indicates the user-adjustable minimum power level setting. When this power level is activated (by footswitch or handswitch), the 'MIN' indicator will flash. The system defaults to 'MIN' power level 3.
2. MAX	Indicates the maximum power level setting. This setting is always 5. When this power level is activated (by footswitch or handswitch) the 'MAX' indicator will flash.
3. POWER LEVEL DECREASE	Touch this button to decrease the minimum (MIN) power setting to the desired level (from 1 to 5). The level chosen will be shown on the screen.
4. POWER LEVEL	Touch this button to increase the minimum (MIN)

INCREASE 5. SETTINGS 6. VOLUME	power setting to the desired level (from 1 to 5). The level chosen will be shown on the screen. Touch this button to go to the system settings. Touch the plus (+) or minus (-) button to adjust the volume of the activation tones (from level 1 to level 4).

Basic Functionality

Mode Switch

After the generator is setup and powered on, it can be switched between Standby mode and Ready mode by holding the Standby touch key for 3 seconds. In Ready mode, the Standby touch key shows green. In Standby mode, the Standby touch key shows amber and the brightness is lower than Ready mode. In Standby mode, the touch screen rejects any touch operation from user as well as the handswitch and footswitch are void.

Generally the generator shall be set to Standby mode when replace instruments.

Activation

In Ready mode, user can activate system for minimum power output or maximum power output by footswitch or handswitch. The power output is activated when user press the handswitch or step down the footswitch and stopped when user release the handswitch or footswitch.

The right pedal on footswitch is for MAX power output. The left pedal on footswitch is for MIN power output. The handswitch is under knob on shears and controlled by forefinger. The upper button is for MAX power output. The lower button is for MIN power output. The default power level for MIN power is 3. Touch the Power Level Increase or Power Level Decrease button to increase or decrease the MIN power setting.

NOTE: Only MIN power level is adjustable between level 1 and level 5. The MAX power level is always level 5.

When the system switched from the Standby mode to the Ready mode, press the button of handswitch or step down the footswitch, the generator will start the Pre-Run test. During the test, the jaws of the scissors should be kept open. Do not touch other instruments, surgical sheets, patients or other objects. After the test passes, the generator is ready for use. Otherwise, the screen will display an error message to help in troubleshooting.

When activate MAX power output, the generator display flash on the right area where shows the MAX power level. A higher and rapid tone is played continuously simultaneously. When activate MIN power output, the generator display flash on the left area where shows the MIN power level. A lower and slow tone is played continuously simultaneously.

The generator ES01 supports two footswitches. The two footswitches are same. Each of the two footswitches may activate system for power output.

Replace Instruments

When replace Ultrasonic Surgical System handpiece and/or shears, the generator ES01 shall be set to Standby mode. In Standby mode, the Standby touch key shows amber and neither handswitch nor footswitch can activate the system. The user may re-connect the handpiece to the generator, or re-secure the handpiece and the shears.

When instruments installation is done, switch the generator to Ready mode. In Ready mode, the Standby touch key shows green and the user can perform System Test (see System Test section on [Settings](#) section) to verify the installation.

Shutdown

Turn off the generator using the On/Off switch on the back panel of the generator.

Settings

Volume

Adjust the volume using the plus or minus button (from level 1 to level 4). Default volume setting is 2.

Minimum Power Level

Touch the power level increase or the power level decrease button to increase or decrease the level of minimum power (from level 1 to level 5). Default minimum power setting is 3.

Settings Navigation

Touch the 'Settings' button to make adjustments in the system settings and access system information. On the settings screen, touch any options to highlight this one then down to the specification setting screen on this option. Touch the close button 'X' to exit from anywhere within the Settings menu. Touch the back arrow to go to the previous menu from the submenu screens.


Figure 7 Setting Screen

Display Brightness

Select 'Display Brightness' to adjust the brightness.


Figure 8 Display Brightness Screen

Touch the plus or minus button to adjust the brightness. Touch the back arrow or close button after adjusting the settings to exit. Once a new brightness setting is selected that setting becomes the new system default.

Activation Set

In the activation set, there are two submenus: "Hand/Foot Activation" and "Intelligent tissue perception".


Figure 9 Activation Set Screen

Hand/Foot Activation

Select "Hand/Foot Activation" to select activation mode among "Both Hand & Foot Active", "Only Hand Active" and "Only Foot Active". A dot and the text highlighted indicate current selection.

Figure 10 Hand/Foot Activation Screen

System Test

There needs to run System Test to ensure products operates properly when new shears or new handpiece is installed to use. During the system test, please keep the jaw opened, ensure the blade is not in contact with any other instruments, drapes, patients and other objects, Safety measures (in accordance with hospital protocol) taken in the presence of aerosols should be in effect during the system test. The Testing screen is like below:


Figure 11 System Test Screen (Running)

If all tests passed, the "Test passed" is showed. Press the OK button generator goes back to Settings screen. If any test failed during test, the generator enters Error screen. See Error Messages and Troubleshooting section.


Figure 12 System Test Screen (Done)

Language

Select 'Language' to adjust the language displayed on the generator. A dot and the text highlighted indicate current selection.

Figure 13 Language Screen

When user changes the language, for example, from Simplified Chinese to English, the system shows the Language Change Confirm Screen:

Figure 14 Language Change Confirm Screen

Press '✓' to confirm to switch to the new language. Press 'X' to cancel switch and the language is not changed.

System Information

Select "System Information" to check information of the generator, for example, software version.


Figure 15 System Information Screen

Error Messages and Troubleshooting

The generator ES01 enters Error mode with audible indicator when system failed. In Error mode, the screen shows an error notification on left top and the error information in the middle. The error information usually contains brief description, resolving suggestion and digital code of the error. Some error screen may have an OK button. Press the OK button the generator will go back to normal state. The digits error code helps to identify the root cause by referring this manual. Here is an example:


Figure 16 Error Screen

See table below for the information and troubleshooting suggestion of each error code:

Table 2 Error Code Messages and Troubleshooting

Error Code (3 digits)	Root Cause	Troubleshooting Suggestions
007	The activation button of the handswitch is pressed when installing the handpiece.	Release the activation button of the handswitch and reinstall the handpiece.
008	The activation button of the handswitch is detected pressed while the connection signal of the shears is not detected.	Check the connection between the handpiece and generator, Please reconnect it if the connection is loose,
009	The activation button of the 2 nd footswitch is pressed when installing the 2 nd footswitch.	Release the activation button of the 2 nd footswitch and reinstall the footswitch.
010	The activation button of the 1 st footswitch is pressed when installing the 1 st footswitch.	Release the activation button of the 1 st footswitch and reinstall the footswitch.
011	The activation button of the 2 nd footswitch is detected pressed while the connection signal of 2 nd the footswitch is not detected.	Check the connection between the 2 nd footswitch and generator. Please reconnect it if the connection is loose.
012	The activation button of the 1 st footswitch is detected pressed while the connection signal of the 1 st footswitch is not detected.	Check the connection between the 1 st footswitch and generator. Please reconnect it if the connection is loose.
015	Both handswitch and footswitch activation signals are detected.	Release the activation button or pedal.
016	Two footswitch activation signals are detected simultaneously.	Release the activation pedals.

017	Two handswitch activation signals are detected simultaneously.	Release the activation buttons.
018	The activation button of handswitch is stuck	Please retry activation, if the problem persists, replace the shears.
022	Power amplifier circuits fault.	Please retry activation, if the problem persists, please contact your representatives.
029	Impedance is bigger than a threshold during prerun-test.	Tighten assembly and reactivate.
030	Impedance is over limit.	Try to re-activate to continue. If the error persists, use the test tip to locate the cause of the error. Use the test tip instead of the shears for the system test. If the test is passed, replace the shears. Otherwise, replace the handpiece.
031	Multiple occurrences of impedance exceeding limit or failure to lock resonant frequency.	Replace the shears.
056	Handpiece disconnected during testing process or activation process.	Reconnect the handpiece to the generator.
081	Handpiece inner circuit is detected shorted during the test.	Retry the test. If the error persists, replace the handpiece.
082	The handpiece cable circuit is detected open.	Replace the handpiece.
083	Handpiece inner circuit is detected shorted during activation.	Try to re-activate to continue. If the error persists, replace the handpiece.
111	Power MOSFET is over temperature.	Power off the generator, ensure the air vent is not blocked by anything, then power on the generator again.
112	Main board is over temperature.	Power off the generator, ensure the air vent is not blocked by anything, then power on the generator again.
152	Power MOSFET bias voltage initialization error.	Restart the generator.
166	Power supply failure.	Restart the generator.
171	Output over current.	Restart the generator.
180	The fan is blocked.	Restart the generator, if the problem persists, please contact your representatives.
181	Failed to lock resonant frequency.	Tighten assembly and reactivate.

184	Too much pressure on the blade.	Release the pressure on the blade. Try to activate to continue.

If error still occurs after instructions listed above, contact your local Ezisurg Medical Co., Ltd sales representative.

Cleaning and Disinfection Instructions

Before cleaning, thoroughly inspect the device(s) for any signs of damage, cracks, or improper mechanical function. Do not use the device(s) if there are signs of damage. Discard and replace device or send device to Ezisurg Medical Co., Ltd or an authorized Ezisurg Medical Co., Ltd service facility for repair where appropriate, if damage or degradation is present.

1.Wipe the surfaces of the device with one CaviWipes™1 towelette or a cloth wetted in other effective disinfectant which is approved for its efficacy by FDA and/or Health Canada approval, VAH/DGHEM-listing, CE marking, etc., until the device surface is free of visible soils. Prior to disinfection, confirm that no visible soil remains on the device.

2.Use a second CaviWipes™1 towelette to thoroughly wet the surfaces. Repeated use of the towelette may be required to ensure that the surfaces remain visibly wet for 1 minute at room temperature.

3.Leave the surfaces to dry on its own.

The above-mentioned reprocessing process (cleansing disinfection) has been successfully validated. Refer to test report: Ezisurg Manual Chem Dis Report, Report No.: MDS-RECD-251105-246

Notes: Operators should refer to Standards & Recommended Practices for additional guidance on cleaning; e.g., ISO 15883. Please follow the manufacturer's instructions of disinfectants.

Other Instruments and Accessories

For Ultrasonic Surgical System handpiece, shears and other accessories, please see their user manual for instructions on cleaning, disinfection and sterilization.

Maintenance and Repair

Visual inspection for cleanliness of the products.

After the cleansing disinfection, a thorough inspection and maintenance ensures that the device is fit for use.

1.Check that the device has no dents, cracks, deformations, scratches, etc.;

2.Check all markings on the product for clear visibility.

Discard and replace any components as necessary.

Do not use the device with following defects: material deformation, damaged insulation or other change in the material, etc.

Test and maintain the device acc. to manufacturer's instruction.

Electrical Safety Checks

The hospital shall contact Ezisurg Medical Co., Ltd or an authorized Ezisurg Medical Co.,

Ltd service facility for periodically (at least once a year recommended) electrical safety check performed by qualified service personnel.

The generator ES01 contains a Potential Equalization Terminal on the back panel. This is provided for compatibility with other medical systems requiring such connections. This conductor is not intended for protective earthing.

Do not open the generator performed by personnel out of Ezisurg Medical Co., Ltd or an authorized Ezisurg Medical Co., Ltd service facility. Open the generator voids the generator warranty.

Replace Main Fuses

If the generator ES01 failed to power on, please check the main fuse status. Locate the fuse holder on the back panel of the generator (see Figure). Remove the fuse holder and examine the fuse. If the fuse is broken, replace the fuse with the appropriate type and rating. Then insert the fuse holder. See System Specifications for fuse type and rating.

Software Update

The generator ES01 software version may be displayed on the System Information screen. If the software needs update, contact Ezisurg Medical Co., Ltd or an authorized Ezisurg Medical Co., Ltd service facility. Do not update the software by personnel out of Ezisurg Medical Co., Ltd or an authorized Ezisurg Medical Co., Ltd service facility.

Disposing of generator ES01 (Environmental Protection)

The generator ES01 and accessories must not be disposed of at the end of life with other waste. To recycle the waste equipment, contact Ezisurg Medical Co., Ltd or your local sales representative to discuss local waste disposal solutions and processes. The generator ES01 poses disposal risks similar to consumer electronics such as computers. There are no radioactive substances, batteries, or hazardous liquids that may leak in the generator ES01.

System Information

Main Power Wiring Diagram


Figure 17 Main Power Wiring Diagram of Generator ES01

1. Universal AC Mains Input
2. Power Entry Module (power receptacle, filter)
3. Protective Earth Connection to Chassis
4. Mains Fuses, 3.15A, 250V
5. Mains Switch
6. 48VDC Medical Grade Power Supply
7. Generator Main PCB
8. Isolated Transducer Drive / Patient Connection
9. Isolation Barriers
10. Functional Earth Connections to Chassis

Assembly Drawing


Figure 18 Assembly Drawing of Generator ES01

1. 48VDC Medical Grade Power Supply
2. Generator Main PCB
3. Filter (with Fuse Holder)
4. Mains Switch
5. Potential Equalization Terminal
6. Footswitch Receptacle
7. Cooling Fans
8. Speaker
9. Handpiece Receptacle
10. LED and touch key PCB
11. Display and Touchscreen PCB with fixed metal
12. Display and Touchscreen fixed metal
13. Display and Touchscreen
14. Front Plastic Panel
15. Ornament Panel for Display and Touchscreen
16. Cover metal of Chassis
17. Back metal of Chassis
18. Bottom metal of Chassis
19. Footswitch

Precautions

1. Minimally invasive procedures should be performed only by persons have adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
2. Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure. Please check carefully on the invasive parts, especially on the surface roughness, sharp edges or outcrop, which might cause safety hazard of the surgery.
3. Handpieces are not sterilized when manufactured. Medical institutions should sterilize the handpiece prior to use. Please refer *Ultrasonic Surgical System Handpiece User Manual* for details.
4. The torque wrench and test tip are not sterilized ,and are only used in non-sterile environments. They must not be used in sterile environments.
5. Please use the torque wrench provided in the shears package to secure the shears and handpieces. Using torque wrench provided other than Ezisurg Medical Co., Ltd might not be able to secure shears and handpieces in appropriate torque, which might cause products broken or decrease the system performance.
6. Operators should perform system test when shears are first in use. During the system test, please keep the jaw opened, ensure the blade is not in contact with any other instruments, drapes, patients and other objects, Safety measures (in accordance with hospital protocol) taken in the presence of aerosols should be in effect during the system test.
7. When replace shears or re-secure shears and handpiece, please switch generator to Standby mode.

8. Ultrasonic Surgical System is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. It is possible to create sparks by hitting other metal instruments. Sparks may ignite flammable gases such as bowel gas.
9. After acting on tissues for a long time, the distal end of the blade, the clamp and the distal end of outer tube might have higher temperature than 41 degree Celsius. Please be careful to operate and avoid to in contact with tissue, drapes, and other parts in patients.
10. When tissue stuck on the distal end of the blade hence decrease the effect of cutting and/or coagulating, please remove the shears from patients, open the jaw, put the distal end of the blade inside normal saline. Then activate system with MAX output to clear the stuck tissue.
11. Operator should activate the system after the jaw is fully closed (clasp the trigger until sense the force feedback and/or hear the notification tone) when using shears to cut tissue. Stop activating immediately after cutting and/or coagulating operation is done. Otherwise the tissue pad might be worn out soon.
12. When using shears to grip tissue to cut and/or coagulate, please ensure the tissue is not beyond the best working range. Deeply grip tissue might decrease the effect of cutting and/or coagulating.
13. Generally, it is better not to keep activating system longer than 15 seconds.
14. The easyUS™ Ultrasonic shears is disposable instrument. It is designed to be used in one surgery. The instrument should be disposed after surgery. Resterilization is not allowed.
15. Do not use if sterilized package is damaged.
16. When generator prompts error messages and/or error tones, it is an indicator that the shears or handpiece is not operating properly. It may be due to the handpiece is beyond its useful life or the shears has not been attached properly, which may result in abnormally high temperatures and user or patient injury.
17. If the handpiece temperature becomes uncomfortable to hold, stop use to prevent burn injury,
18. As with all energy sources (Electrosurgery, Lasers, or Ultrasound) there are concerns about the carcinogenic and infectious potential of the by-products such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
19. Do not activate an electrosurgical device in close proximity to the Ultrasonic Surgical System, it may cause the touch key or the touchscreen misoperation by electromagnetic interference. This situation may cause user or patient injury, should be avoid.
20. The aerosols creates by the activations of Ultrasonic Surgical System shears in fatty tissue are potentially flammable.

21. To avoid user or patient injury in the event that accidental activation occurs, the shears' blade should not be in contact with the patient, drapes or flammable materials while not in use.
22. During and following activation in tissue, the shears' blade may become hot. Avoid unintended blade contact with tissue, drapes, surgical gowns, or other unintended sites after activation.
23. The handpiece meets the international safety standard EN 60601-1 for user contact and is not intended for patient contact.
24. Handle the handpiece carefully. If the binder bolt is damaged the handpiece resonant frequency may be shifted and decrease the performance.
25. Do not bang or drop the handpiece.
26. Do not clean the handpiece electrical connector with alcohol.
27. After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
28. Stop the operation immediately and inspect the shear or handpiece to confirm whether they have been assembled well, if the shear or handpiece does not work properly or a buzzing is heard during surgery.
29. Do not operate the generator ES01 in a moist environment as a shock hazard may exist. If liquids have entered the generator ES01 the unit must be returned to the manufacturer for testing prior to use.
30. Do not operate the unit in close proximity to volatile solvents such as methanol or alcohol as explosion may occur.
31. Avoid use of generator ES01 adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, monitor the generator ES01 and the other equipment to assure normal operation.
32. Use of accessories and cables other than those specified may result in unpredictable performance, increased electromagnetic emissions, or decreased electromagnetic immunity.
33. No internal user serviceable parts. For service, return the generator to Ezisurg Medical Co., Ltd or an authorized Ezisurg Medical Co., Ltd service facility.
34. To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
35. This device is not MR safe and is not MR compatible.
36. In case of system failure, ensure availability of the appropriate backup equipment relevant to the specific procedure.
37. The touch screen display of the generator is very sensitive. Do not use sharp metal objects on the touch screen.
38. Opening of the generator ES01 invalidates the warranty and could create hazardous condition.
39. Read instructions prior to use and follow the hospital's clinical practice guidelines for

ultrasonic surgery, electrosurgery, gynecology and laparoscopic procedures.

40. Replace fuses only with the appropriate type and rating. See System Specifications.
41. Do not sterilize the generator ES01. Sterilization will damage the unit.
42. Do not restrict the openings on the bottom and the back panel of the generator ES01, as they provide the required airflow for cooling.
43. If electromagnetic interference with other equipment is suspected, reorient the device or remove possible sources of interference (for example, cellular phones, radios, etc.) from the room.
44. When Ultrasonic Surgical System is used together with endoscopic equipment, the patient leak current increase hazard may exist.
45. There is no specific requirement for interconnection when Ultrasonic Surgical System is used together with endoscopic equipment.
46. The distal end of the blade may have higher temperature than 41 degree Celsius. Be careful to operate and avoid damaging unintentional objects.
47. The Ultrasonic Surgical System shall not be serviced or maintained while in use with a patient.
48. Do not position the generator ES01 so that it is difficult to operate the disconnection device while the power cord separated.
49. No modification of this equipment is allowed.
50. Do not reboot the generator immediately (in less than 0.5s) after powering off the generator, it may cause the screen can not start up normally.
51. Successful hemostasis may require adjunct measures when the device is used on solid organs. Due to the difficulty of visualizing internal structures, proceed slowly and do not attempt to transect large masses of tissue in one activation. Avoid the division of large vascular/biliary bundles when using the instrument under these conditions.
52. Device has not been evaluated in main vessels of the central circulatory system and is not intended for use in the following named vessels: arteriae pulmonales, aorta ascendens, arcus aortae, aorta descendens to the bifurcatio aortae, arteriae coronariae, arteria carotis communis, arteria carotis externa, arteria carotis interna, arteriae cerebrales, truncus brachiocephalicus, venae cordis, venae pulmonales, vena cava superior, and vena cava inferior.

System Specifications

Output Characteristics

See table below for the ultrasonic output characteristics of the shears SHC14A、SHC23A、SHC36A、SHC45A、SHF9A and SHF17A :

SHC14A SHC23A SHC36A SHC45A	Reference primary tip vibration excursion: 30µm~120µm Primary acoustic output area: 2.4 mm ² (±30%) Drive frequency: 55.5kHz±1kHz Frequency modulation type: Non-modulating
SHF9A SHF17A	Reference primary tip vibration excursion: 30µm~120µm Primary acoustic output area: 1.6mm ² (±30%) Drive frequency: 55.5kHz±1kHz Frequency modulation type: Non-modulating
NOTE: Shears SHC14A、SHC23A、SHC36A and SHC45A shall be used together with generator ES01 and handpiece eHP01; Shears SHF9A、SHF17A shall be used together with generator ES01 and handpiece eHP02.	

Conformance to Standards

The Ultrasonic Surgical System conforms to the following international standards:

EN (IEC) 60601-1	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.
EN (IEC) 60601-1-2	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests.

System Specifications

Product Codes	Required Components for System Operation: Generator: ES01 Handpiece: eHP01 Instruments: SHC14A、SHC23A、SHC36A、SHC45A Optional Components: Footswitch: FSU01
	Required Components for System Operation: Generator: ES01 Handpiece: eHP02 Instruments: SHF9A、SHF17A Optional Components: Footswitch: FSU01
Degree of Protection Against Electric Shock	Type CF Applied Part
Class of Protection Against Electric Shock	Class I
Degree of Protection Against Harmful Ingress of Water	ES01: IPX1 FSU01: IPX8
Main Fuse	F3.15AL250V (Fast-acting, 3.15Amp, glass body, 5x20 mm package size, quantity: 2)
Main Input	Voltage: 110-240 VAC Frequency: 50/60 Hz Current Consumption: 3 Amp
Output	Voltage: 125 VAC RMS maximum Power: 35 watts maximum Frequency: 55.5 kHz

Ambient Operating Conditions	Temperature: 10°C to 28°C Humidity: 30%-75% non-condensing Atmospheric Pressure Range: 700hPa-1060hPa
Power Cord	International removable power cord set with the following characteristics: Plug Style: as needed by particular country requirements Receptacle: IEC 60320 C13 with straight non-angled cord entry Cord Length: 2.5m maximum Current Rating: 10A Minimum conductor size cross-sectional area: 1.0 mm ² copper Voltage Rating: 250 VAC minimum Wiring: International
Weight	ES01: 8.2 Kg FSU01: 1.9 Kg
Overall Dimensions	ES01: 370 mm x 380 mm x 150 mm FSU01: 350 mm x 190 mm x 65 mm

Electromagnetic Compatibility (EMC)

Electromagnetic Emissions

Guidance and manufacturer's declaration – electromagnetic emissions		
Ultrasonic Surgical System is intended for use in the electromagnetic environment specified below. The customer or the user of Ultrasonic Surgical System should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Ultrasonic Surgical System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The Ultrasonic Surgical System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Electromagnetic Immunity

Guidance and manufacturer's declaration – electromagnetic immunity			
Ultrasonic Surgical System is intended for use in the electromagnetic environment specified below. The customer or the user of Ultrasonic Surgical System should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
ESD IEC 61000-4-2	±8kV contact ±15kV air	±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line to line ±2kV line to earth	±1kV line to line ±2kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0%U _T , 0.5cycle (At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315°) 0%U _T , 1 cycle 70%U _T , 25/30 cycles Single phase: at 0° 0%U _T , 250/300 cycles	0%U _T , 0.5cycle (At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315°) 0%U _T , 1 cycle 70%U _T , 25/30 cycles Single phase: at 0° 0%U _T , 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment.

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of Ultrasonic Surgical System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80MHz to 800MHz $d = 2.3\sqrt{P}$ 800MHz to 2.7GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m 80 MHz to 2,7 GHz	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Ultrasonic Surgical System is used exceeds the applicable RF compliance level above, Ultrasonic Surgical System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating Ultrasonic Surgical System.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.


Recommended separation distances between portable and mobile RF communications equipment and Ultrasonic Surgical System

Recommended separation distances between portable and mobile RF communications equipment and Ultrasonic Surgical System			
Ultrasonic Surgical System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Ultrasonic Surgical System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Ultrasonic Surgical System as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Storage and transportation Conditions:

Temperature:

Humidity:

Atmospheric Pressure Range:

Storage of processed products in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.

Contact information

Ezisurg Medical Co., Ltd.

399, Miaoqiao Rd., Pudong New District, 201315 Shanghai, P.R. China

Website: www.ezisurg.com

Tel.1: 4000906176

Tel.2: +86-21-50456176

Fax: +86-21-50676156

Shanghai International Holding Corp. GmbH (Europe)

Add: Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175

Fax: +49-40-255726

For a patient/ user/ third party in the European Union and in countries with identical

regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/ or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

For additional information, contact your Ezisurg Medical Co., Ltd. Sales Representative or Ezisurg Medical Co., Ltd. Customer Service.

No electronic version of the IFU provided.

Symbols

Caution

Affix on the device label and the front panel

Keep dry

Affix on the device package

Fragile
handle with care

Affix on the device package

This way up

Affix on the device package

Keep away from
sunlight

Affix on the device package

Temperature limit

Affix on the device package

Humidity limitation

Affix on the device package

Atmospheric
pressure limitation

Affix on the device package

Use no hooks

Affix on the device package

Stacking limit by
number

Affix on the device package

Model number

Affix on the device label and bar code label

Serial Number

Affix on the bar code label

Date of manufacture
Country of
manufacture: CHINA

Affix on the device label and
bar code label

Manufacturer

Affix on the device label

Authorized
representative in the
European Community

Affix on the device label

Return waste to a
collection system or
treatment and
recycling facilities.
Applicable in EU.
Follow
decontamination
instructions before
returning waste.

Affix on the device label

Refer to instruction
manual/booklet

Affix on the device label


Device meets the
requirements of
Medical Device
Directive 93/42/EEC,
as amended by
2007/47/EC. 0197 is
the code of notified
body TUV Rheinland.

Affix on the device label

Equipotentiality

Affix on the device label

Fuse

Affix on the fuse label

Foot switch

Affix on the device label

on (power)

Affix on the power switch

off (power)

Affix on the power switch

Stand-by

Affix on the front panel

Type CF Applied Part

Affix on the front panel

Unique device
identifier

Affix on the bar code label

www.ezisurg.com

Consult instructions
for use or consult
electronic instructions
for use

Affix on the transport
package label

Medical device

Affix on the bar code label
and transport package label

Package quantity

Affix on the transport
package label

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