**EU REP** Shanghai International Holding Corp. GmbH (Europe)  
Add: Eiffeistrasse 80, 20537 Hamburg, Germany  
Tel: +49-40-2513175  
Fax: +49-40-255726

 Ezisurg Medical Co., Ltd.  
No.399, Miaoqiao Rd., Pudong New District, Shanghai, P.R. China  
Website: www.ezisurg.com  
Tel: +86-21-50456176  
Fax: +86-21-50676156

CE 0197

Rev.E  
2026-07-01

# easyEndo™ UniPlus

|       |                                                                  |
|-------|------------------------------------------------------------------|
| en    | Single Use Powered Endoscopic Staplers                           |
| pl    | Jednorazowy Zasilany Stapler Endoskopowy                         |
| cs    | Jednorázové napájené endoskopické svorkovače                     |
| pt-pt | Agrafadores Endoscópicos Eléctricos Descartáveis                 |
| es    | Grapadoras Endoscópicas Motorizadas de Un Solo Uso               |
| it    | Cucitrici Endoscopiche Alimentate Monouso                        |
| de    | Einweg-betriebene endoskopische Klammernahtgeräte                |
| fr    | Agrafeuse endoscopique électrique à usage unique                 |
| nl    | Aangedreven Endoscopische Nietmachine voor eenmalig gebruik      |
| sv    | Engångsdriven Endoskopisk Häftapparat                            |
| ro    | Capsatoare Electrice Endoscopice de Unică Folosință              |
| el    | Ηλεκτροδοτούμενα ενδοσκοπικά συρραπτικά μίας χρήσης              |
| hu    | Egyszer használatos, elektromos működtetésű endoszkópos varrógép |


## IMPORTANT!


Please read the following instructions carefully before using this product. This document is designed to provide instructions for using this product. It is not a reference to surgical techniques. Failure to follow the instructions may lead to serious surgical consequences, such as anastomotic potential for bleeding, leakage, or disruption.

## Device Overview


1.

2.


(a)

(b)


(c)

3.


1


4.


5.


6.

7.


8.


2


9.


10.


11.


12.


13.


14.

15.

16.

## English

### Device Description

easyEndo™ UniPlus Single Use Powered Endoscopic Staplers are sterile, single patient use instruments that simultaneously cut and staple tissue. The instruments can be used in open or minimally invasive surgery. A primary lithium Battery Pack is included in the stapler package, which must be installed prior to use. Refer to the G. Battery Pack Disposal section for specific requirements for disposing of the Battery Pack. The Powered Stapler is packaged without a Reload, which must be loaded prior to use. Refer to section B for both the Battery Pack and Reload installation. The Curved Tip is packed with the stapler with specific models (refer to Table 2).

Table 1. List of device features

|                                     |                           |                              |
|-------------------------------------|---------------------------|------------------------------|
| 1. Pin                              | 11. Power Indicator       | 21. Cutting Knife            |
| 2. Stapler Loading Indicating Arrow | 12. Forward Indicator     | 22. Curved Tip               |
| 3. Shaft                            | 13. Backward Indicator    | 23. Cartridge                |
| 4. Articulation Knob                | 14. Safety Indicator      | 24. Anvil Jaw                |
| 5. Rotating Knob                    | 15. Reset Wrench          | 25. Tissue Stop Wedge        |
| 6. Close/Fire Button                | 16. Load/Lock Button      | 26. Loading Indicating Line  |
| 7. Open/Back Button                 | 17. Reset Cover           | 27. Loading Indicating Arrow |
| 8. Safety Button                    | 18. End of Staple Line    | 28. Shipping Wedge           |
| 9. Handle                           | 19. End of Cut Line       | 29. Staple                   |
| 10. Battery Pack                    | 20. Knife Blade Indicator |                              |

Table 2. Stapler Product Codes Chart

| REF    | Shaft Length (L1)<br>(mm) | Outer Diameter of Shaft (D1)<br>(mm) | Description                                               |
|--------|---------------------------|--------------------------------------|-----------------------------------------------------------|
| PU12L  | 250                       | 12                                   | Single Use Powered Endoscopic Stapler, Long               |
| PU12M  | 160                       | 12                                   | Single Use Powered Endoscopic Stapler, Medium             |
| PU12S  | 70                        | 12                                   | Single Use Powered Endoscopic Stapler, Short              |
| PU12LG | 250                       | 12                                   | Single Use Powered Endoscopic Stapler, Long, Curved tip   |
| PU12MG | 160                       | 12                                   | Single Use Powered Endoscopic Stapler, Medium, Curved tip |
| PU12SG | 70                        | 12                                   | Single Use Powered Endoscopic Stapler, Short, Curved tip  |

Table 3. Reload Product Codes Chart

| REF   | Cartridge Color | Staple Line Length (L2)<br>(mm) | Staples Rows<br>(pattern) | QTY /<br>row | Staple height H (mm)                    |                     | Thickness range*(mm) |
|-------|-----------------|---------------------------------|---------------------------|--------------|-----------------------------------------|---------------------|----------------------|
|       |                 |                                 |                           |              | Open                                    | Closed              |                      |
| P30M  | Grey            | 30                              | 6(staggered)              | 8            | 2.0                                     | 0.75                | 0.75 - 1.0           |
| P30W  | White           |                                 |                           |              | 2.5                                     | 1.0                 | 1.0 - 1.5            |
| P30B  | Blue            |                                 |                           |              | 3.5                                     | 1.5                 | 1.5 - 2.0            |
| P30VM | Brown           |                                 |                           |              | Inner: 2.0<br>Middle: 2.5<br>Outer: 3.0 | 0.75<br>1.0<br>1.25 | 0.88-1.5             |

|       |        |    |              |    |                                         |                     |            |
|-------|--------|----|--------------|----|-----------------------------------------|---------------------|------------|
| P45M  | Grey   | 45 | 6(staggered) | 11 | 2.0                                     | 0.75                | 0.75 - 1.0 |
| P45W  | White  |    |              |    | 2.5                                     | 1.0                 | 1.0 - 1.5  |
| P45B  | Blue   |    |              |    | 3.5                                     | 1.5                 | 1.5 - 2.0  |
| P45C  | Gold   |    |              |    | 3.8                                     | 1.7                 | 1.7 - 2.2  |
| P45G  | Green  |    |              |    | 4.1                                     | 2.0                 | 2.0 - 2.5  |
| P45T  | Black  |    |              |    | 4.4                                     | 2.3                 | 2.3-3.0    |
| P45VM | Brown  | 45 | 6(staggered) | 11 | Inner: 2.0<br>Middle: 2.5<br>Outer: 3.0 | 0.75<br>1.0<br>1.25 | 0.88-1.5   |
| P45MT | Purple |    |              |    | Inner: 3.5<br>Middle: 3.8<br>Outer: 4.1 | 1.25<br>1.5<br>1.75 | 1.5-2.25   |
| P45XT | Black  |    |              |    | Inner: 3.8<br>Middle: 4.1<br>Outer: 4.4 | 1.75<br>2.0<br>2.25 | 2.25-3.0   |
| P60M  | Grey   |    |              |    | 2.0                                     | 0.75                | 0.75 - 1.0 |
| P60W  | White  |    |              |    | 2.5                                     | 1.0                 | 1.0 - 1.5  |
| P60B  | Blue   |    |              |    | 3.5                                     | 1.5                 | 1.5 - 2.0  |
| P60C  | Gold   | 60 | 6(staggered) | 15 | 3.8                                     | 1.7                 | 1.7 - 2.2  |
| P60G  | Green  |    |              |    | 4.1                                     | 2.0                 | 2.0 - 2.5  |
| P60T  | Black  |    |              |    | 4.4                                     | 2.3                 | 2.3-3.0    |
| P60VM | Brown  |    |              |    | Inner: 2.0<br>Middle: 2.5<br>Outer: 3.0 | 0.75<br>1.0<br>1.25 | 0.88-1.5   |
| P60MT | Purple |    |              |    | Inner: 3.5<br>Middle: 3.8<br>Outer: 4.1 | 1.25<br>1.5<br>1.75 | 1.5-2.25   |
| P60XT | Black  |    |              |    | Inner: 3.8<br>Middle: 4.1<br>Outer: 4.4 | 1.75<br>2.0<br>2.25 | 2.25-3.0   |

\*Tissue thickness range: Porcine tissue measured at 8g/mm2 prior to firing.

Compatibility

All easyEndo™ UniPlus Single Use Powered Endoscopic Staplers are compatible with all easyEndo™ UniPlus series Reloads.

Type CF Applied part

Single Use Reloads and Curved Tip (refer to figure2)

Intended Use/Purpose

The instruments are intended for transection, resection, and/or creation of anastomoses in open and minimally invasive general abdominal, and thoracic surgical procedures.

Indications

The stapler, when used with easyEndo™ UniPlus series Reloads, has applications in abdominal, and thoracic surgery for resection, transection, and creation of anastomosis. It may be used for transection and resection of liver substance, hepatic vasculature, and biliary structures and for transection and resection of the pancreas.

The stapler, when used with Curved Tip, can be used to blunt dissect or separate target tissue from other certain tissue.

Contraindications

- 1. Do not use the instruments on the aorta.
- 2. The device should not be used to staple tissues that are necrotic, friable, or have altered integrity, e.g., ischemic or edematous tissues.
- 3. Do not use any linear cutter on major vessels without making provision for proximal and distal control.
- 4. Tissue thickness should be carefully evaluated before firing any stapler. Refer to Table 3. Reload Product Codes Chart for tissue compression requirement (closed staple height) for each staple size. If tissue cannot comfortably compress to the closed staple height, or easily compresses to less than the closed staple height, the tissue is contraindicated as it may be too thick or too thin for the selected staple size.
- 5. These instruments are not intended for use when surgical stapling is contraindicated.

Target Patient Group

The intended patient population is adults.

Intended Users

The Staplers and Reloads for Single Use should only be used by professional medical staff who have been fully trained in the use of The Staplers and Reloads for Single Use techniques and surgeries.

Clinical Benefits

- Reduced anastomotic bleeding and anastomotic leakage rate;
- Reduced hospital stay and operative time.

Undesirable Side Effects/Residual Risks

Undesirable side effects and risks associated with powered surgical stapling include the potential for bleeding, tissue injury, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, hypersensitivity, leakage, Biological toxicity, property or environmental damage, electrical shock, and foreign body magnetic resonance incompatibility. Also, unintended harm, extended surgery, or altered surgical approach may result from staple line failures, inability to cut, or damaged devices.

Safety Information

The Use of Warning, Caution, and Note Statements

Information relative to the completion of a task in a safe and thorough manner will be supplied in the form of a Warning, a Caution, or a Note statement. These statements are found throughout the documentation.

These statements should be read before continuing to the next step in a procedure.

- **WARNING:** A Warning statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in personal injury or loss of life.
- **CAUTION:** A Caution statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in damage to or destruction of the equipment.
- **NOTE:** A Note statement indicates an operating or maintenance problem, practice, or condition that is necessary to accomplish a task efficiently.

Warnings and Precautions

1. Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
2. easyEndo™ UniPlus Staplers can only be used with the easyEndo™ UniPlus series Reloads manufactured by Ezsurg. Verify compatibility of all staples and accessories prior to using the Stapler (refer to Table 2, Table 3).
3. This product has never been verified for cross-assembly with reloads of similar devices of other brands. Do not cross-assemble the reloads of other brands with this product during the operation. Otherwise, it may cause stapler damage, failure to fire, staple malformation, anastomosis failure, etc.
4. Any type of battery other than that provided with the Stapler may cause unpredictable performance problems, increased electromagnetic radiation, or reduced electromagnetic immunity of the Stapler.
5. Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers (such as trocar) are employed together in a procedure, verify compatibility prior to initiation of the procedure.
6. Radiation therapy before surgery may cause tissue changes, resulting in tissue thickness exceeding the applicable range of the selected staple. Therefore, when choosing the staple specification, consider the patient's various preoperative treatments.
7. This device is packaged and sterilized for single use only. Do not re-use, reprocess, or re-sterilize. Re-use, reprocessing, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure, that in turn may result in patient injury, illness, or death. Also, reprocessing or re-sterilizing single-use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to the patient's injury, illness, or death.
8. Prepare a full set of products as a backup during the operation. This will enable the operation to continue on if one Stapler malfunctions during the operation.
9. Do not use if sterile package is damaged.
10. Do not use after expiry date.
11. Make sure that the Curved Tip is correctly and securely connected to the easyEndo™ UniPlus Single Use Reloads Anvil Jaw.
12. Always inspect the tissue thickness and select an appropriate staple size prior to specifying the use of the easyEndo™ UniPlus Single Use Powered Endoscopic Staplers and Reloads. Refer to Table 3. Reload Product Codes Chart for proper cartridge selection. Overly thick or thin tissue may result in unacceptable staple formation.
13. The use of Curved Tip Stapler is contraindicated on tissue or structures that cannot fit completely within the jaws proximal to the transitional angle of the Curved Tip.
14. When the tissue is too thick, it may trigger the safety mechanism causing the motor to stall and the Cutting Knife to stop. Press the "Open/Back" Button to open the jaws, and replace the Reload as described in section **E. Reloading the Instrument**. Then, position the instrument around a smaller tissue section or use a more appropriate Reload.
15. The Reload is packaged in the open position. Do not attempt to close the Reload.
16. Do not attempt to remove the Shipping Wedge until the Reload is loaded into the instrument.
17. Do not clamp instrument prior to removing Shipping Wedge.
18. Ensure that the stapler is loaded firmly with the Reload. The application of the instrument loaded not firmly during the operation may lead to separation of the Stapler and the Reload, in such case, an open surgery may be necessary to remove the separated Reload.
19. It is strictly prohibited to use the Reload if Staple is exposed to the surface of Cartridge or the colored drivers can be seen.
20. The Reload cannot be removed from the Stapler only by rotation or pull, otherwise, it is necessary to reload the Reload.
21. If the Reload jaws remain closed (or incompletely open) after pressing the Open/Back Button, unload the Reload and replace a new one.
22. If the jaws of the Reload does not articulate after manipulating the Articulation Knob, unload the Reload and replace with a new one.
23. Do not press the Safety Button while the jaws are closed and the Stapler is not ready to fire. Avoid inadvertently pressing the Close/Fire Button, which may cause the Stapler to be partially or completely fired.
24. Before inserting or removing the stapler from trocar cannula, ensure the jaws are closed and aligned linearly with the Shaft (i.e., straight). Failure to have the jaws in the straight position will result in difficult insertion or withdrawal of the Stapler and may result in damage to the instrument.
25. Do not move the Stapler in a forcible manner at any time (i.e., when it is inserted, removed from the trocar cannula, and during the entire operation), otherwise, it may cause Stapler failure or damage.
26. The instrument Jaw must be completely visible past the trocar sleeve before opening within the body cavity.
27. Do not manipulate the Articulation Knob after grasping the tissue.
28. Clamping and unclamping delicate structures such as venous and bile ducts may damage tissue regardless of whether the stapler is fired.
29. Prior to positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc., are incorporated into the instrument jaws. Firing over obstacles may result in incomplete cutting action and/or improperly formed staples.
30. Placement of tissue proximal to the Tissue Stop Wedge (on the jaws) may result in instrument malfunction. Any tissue extending beyond the cut mark on the jaws, may be transected without Staples and may lead to serious surgical consequences, such as bleeding, leakage, or disruption.
31. When using the Reload with a Curved Tip to manipulate tissue, avoid exerting excessive pressure on fragile structures with the Curved Tip of the device.
32. When using the easyEndo™ UniPlus Single Use Reloads with Curved Tip to manipulate tissue, do not apply excessive force to the Curved Tip, excessive force may separate it from the anvil Jaw.
33. When using the Curved Tip, it is necessary to ensure that the Curved Tip can be seen in the field of vision.
34. In case of a Curved Tip Reload, ensure that the tissue/vessel to be transected does not extend beyond the End of Cut Line on the Reload. The device will only cut to the cut line. Tissue contained within the jaws distal to this line will not be transected.
35. If the clamping mechanism becomes inoperative and the jaws do not clamp on the tissue, do not fire the Stapler. Remove and do not continue to use Stapler.
36. If a Stapler malfunction occurs while applying Staples across a blood vessel, then the user should clamp or ligate the vessel before releasing the Stapler, while the Stapler is still closed on the tissue.
37. Ensure that the tissue lies flat and is positioned properly between the jaws. Any "bunching" of tissue along the Cartridge, particularly in the crotch of the jaws, may result in an incomplete staple line.
38. Before firing, check that no unexpected tissues and organs are clamped by the jaws.
39. Incomplete firing may result in malformed Staples, incomplete cut line, bleeding, and/or difficulty removing the device.
40. Examine the staple lines for pneumostasis/hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.
41. Before removing the instrument, be sure tissue is cleared from the jaws, and then close the jaws. Failure to gently clear tissue from the jaws may result in bleeding or damage to adjacent tissues or structures. Additionally, it may interfere with removal of the instrument through the trocar.
42. When the Stapler needs to be used multiple times in a single operation, used Reloads must be unloaded and replaced with new ones. The device's lock-out feature is designed to prevent an empty single use reload from being fired a second time.
43. The Stapler and Reload which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.

- 44. There is a risk of increased leak rates when staple lines are crossed, even if there may be clinical circumstances when a surgeon may deem it necessary or appropriate to do so.
- 45. In a single operation, the Stapler can be fired at most 25 times by replacing Reload. Each Reload can be fired only once.
- 46. After the manual override system is used, the Stapler is disabled and cannot be used for any subsequent firings.
- 47. Handle and dispose of the Battery Pack in accordance with the Battery Pack Disposal section.
- 48. Do not use autoclaves, ethylene oxide, or radiation to re-sterilize or disinfect the Battery Pack.
- 49. If handled incorrectly, the Battery Pack may cause a fire. Do not disassemble, heat up, autoclave, squeeze, puncture, short-circuit external contact, or charge Battery Pack.
- 50. Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- 51. Portable radio and mobile RF communications equipment can affect medical devices. If it is suspected that there is electromagnetic interference between the Stapler and other equipment, reposition the Stapler or remove possible sources of interference from the room (e.g., Portable radio, mobile RF communications equipment, etc.).
- 52. Use of accessories, transducers and cables other than those specified or provided by the MANUFACTURER of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- 53. PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of this equipment, including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could result.
- 54. easyEndo™ UniPlus Single Use Powered Endoscopic Stapler is considered medical electrical equipment. Medical electrical equipment requires special precaution regarding electromagnetic compatibility (EMC), and needs to be installed and put into service according to the EMC information provided in this document.
- 55. Do not use the Stapler in an aircraft environment.
- 56. Cobalt classified as a CMR substance is present in the device in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
- 57. Hazards can result from unauthorized modification of the ME equipment.
- 58. If the Safety Indicator remains solid red with an electronic beep for approximately 2 seconds, press the Open/Back Button to open the instrument and replace the Reload with a new one. If the Safety Indicator flashes red, the device must be replaced.

## Use instructions

### A. Prior to Using Instrument

- 1. Assess the tissue prior to using the Stapler. The Stapler should not be used to staple tissues that are necrotic, friable, or have altered integrity, e.g., ischemic or edematous tissues, as the staples may be less likely to securely approximate tissue.

### B. Assembling the Stapler

- 1. Using sterile technique, remove the Stapler, Battery Pack, and Reload from their packages. To avoid damage, carefully place the Stapler, Battery Pack, and Reload in the sterile field.
- 2. Install the Battery Pack prior to use. Insert the Battery Pack, with the buckle facing upward, into the Stapler along the slots on the rear of the instrument (refer to Figure 4). An audible click will be heard when the Battery Pack is fully inserted.

#### Note:

- After installing the Battery Pack, the power indicator lights up green for a long time and the electronic prompt tone "beep" is heard. If the Power Indicator light flashes red and the electronic prompt sound "beep" is heard, it indicates that the Battery Pack is low, and it is recommended to replace the device. If the Power Indicator lights up red for a long time and the continuous prompt tone "beep" is heard or the light does not light up, it means that the Battery Pack has been discharged, you must replace the device.

- 3. Select a new Reload with the appropriate staple size for the tissue thickness. Refer to Table 3. Reload Product Codes Chart.

#### WARNING:

Carefully evaluate the tissue thickness prior to use. Overly thick or thin tissue may result in unacceptable staple formation.

#### CAUTION:

- The Reload is packaged in the open position. Do not attempt to close the Reload.
- Do not attempt to remove the Shipping Wedge until the Reload is loaded into the instrument.
- Do not clamp instrument prior to removing shipping wedge.

**NOTE:** The color of the Cartridge identifies the staple size, as described in Table 3. Reload Product Codes Chart.

**NOTE:** The intended tissue in the label are recommended approximate intended tissues.

- 4. Insert the Curved Tip into the Anvil Jaw, as needed.
  - a. Firstly, insert the distal end of Anvil Jaw into the Curved Tip opening cavity above the spring plate (refer to Figure 2a);
  - b. Push the Curved Tip horizontally, and the spring plate on the Curved Tip follows the slot on the back of Anvil Jaw (refer to Figure 2b,2c);
  - c. An audible click will be heard when the Curved Tip is fully inserted.

**Caution:** Make sure that the Curved Tip is correctly and securely connected to the easyEndo™ UniPlus Single Use Reloads Anvil Jaw.

#### NOTE:

- When using the easyEndo™ UniPlus Single Use Reloads, the Curved Tip can be installed and removed as needed.
- When using the Curved Tip Reload, the Curved Tip of the Anvil Jaw can be used to manipulate the tissue and will be clearly visible when the Reload is clamped on the tissue/vessel that has to be transected.
- The Powered Staplers with Curved Tips can be used for blunt dissection.
- 5. Prior to inserting the Reload, ensure that the Stapler's Articulation Knob is neutral position (0° relative) to the instrument.
- 6. To load the staplers with the appropriate Reload (Refer to Figure 5).
  - a. Insert the Pin located at the distal end of the instrument shaft into the Reload.
  - b. Ensure that the Loading Indicating Arrow on the reload aligns with the Stapler Loading Indicating Arrow on the Stapler.
  - c. Push the Reload in and twist clockwise relative to the Stapler, so that the Reload will lock into place and an audible click is heard.
  - d. Ensure that the Loading Indicating Line aligns with the Stapler Loading Indicating Arrow.
  - e. When the Reload is correctly loaded on the stapler, the light blue Load/Lock Button fuse will reset.
- 7. Remove the Shipping Wedge from Reload as depicted in Figure 5. Cycle the instrument to confirm proper loading of the Reload. To do so, close the jaws by pressing Close/Fire Button. Press the Open/Back Button to confirm that the jaws open fully.

#### WARNING:

- Ensure that the Stapler is loaded firmly with the Reload. The application of the instrument loaded not firmly during the operation may lead to separation of the Stapler and the Reload, in such case, an open surgery may be necessary to remove the separated Reload.
- It is strictly prohibited to use the Reload if Staples are exposed to the surface of Cartridge or the colored drivers can be seen (refer to operation of E. Reloading the Instrument below).

#### CAUTION:

- If the Reload jaws remain closed (or incompletely open) after pressing the Open/Back Button, unload the Reload and replace a new one (refer to operation of E. Reloading the Instrument below).
- If the jaws of Reload does not articulate after manipulating the Articulation Knob, unload the Reload and replace with a new one (refer to operation of E. Reloading the Instrument below).

**NOTE:** Check the function of the Stapler. The Forward Indicator will light up when the Close/Fire Button is pressed, and the Backward Indicator will light up when the Open/Back Button is pressed. If this does not occur, replace the Stapler.

- 8. The Stapler is now loaded and ready for use.

### C. Articulating the Stapler

- 1. Introduce Stapler to body cavity, as follows:

- a. Close the jaws by pressing Close/Fire Button.

**CAUTION:** Do not press the Safety Button while the jaws are closed and the Stapler is not ready to fire. Avoid inadvertently pressing the Close/Fire Button, which may cause the Stapler be partially or completely fired.

- b. Ensure the Articulation Knob is in the neutral position, so that the jaws are straight.

- c. Introduce the Stapler into the body cavity through a trocar of the appropriate size, or through an incision (refer to Figure 6).

**CAUTION:** Before inserting, ensure the jaws are closed and linearly aligned with the Shaft (i.e., straight). Failure to have the jaws in the straight position will result in difficult insertion or withdrawal of the Stapler and may result in damage to the instrument.

- 2. Once the jaws have entered body cavity completely, press Open/Back Button to open the jaws. When using a trocar, the jaws must be visible past the trocar sleeve before opening the jaws.

- 3. Position the jaws at an appropriate angle before grasping tissue by using the Rotating Knob, and Articulation Knob, as needed.

- The Rotating Knob can rotate the jaws 360 degrees (refer to Figure 7).
- The Articulation Knob can articulate the jaws left and right, with three gears on the left and right, and the maximum bending angle is 60 degrees (refer to Figure 8).

#### WARNING:

Do not manipulate the Articulation Knob and Rotating Knob after grasping the tissue.

### D. Firing the stapler

- 1. Position the jaws around the tissue to be transected.

#### WARNING:

Prior to positioning the Stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed Staples, and/or inability to open the instrument jaws.

**Note:** The stapler will not cut tissue beyond the cut mark indicated on the jaws. More than one application of the stapler may be necessary for tissue exceeding the length of the Reloads (30mm, 45mm, 60mm).

- 2. Close the jaws across the tissue to be transected by pressing the Close/Fire Button (refer to Figure 9). Holding the jaws in place for 15 seconds after closing, and prior to firing, may result in better compression and staple formation.

#### WARNING:

Placement of tissue proximal to the Tissue Stop Wedge (on the jaws) may result in instrument malfunction. Any tissue extending beyond the cut mark may be transected without Staples and may lead to serious surgical consequences, such as bleeding, leakage, or disruption.

- When using the Reload with a Curved Tip to manipulate tissue, avoid exerting excessive pressure on fragile structures with the Curved Tip of the device.
- When using the easyEndo™ UniPlus Single Use Reloads with Curved Tip to manipulate tissue, do not apply excessive force to the Curved Tip, excessive force may separate it from the Anvil Jaw.
- When using the Curved Tip, it is necessary to ensure that the Curved Tip can be seen in the field of vision.
- In case of a Curved Tip Reload, ensure that the tissue/vessel to be transected does not extend beyond the End of Cut Line on the Reload. The device will only cut to the cut line. Tissue contained within the jaws distal to this line will not be transected.

**CAUTION:** If the clamping mechanism becomes inoperative and the jaws do not clamp on tissue, do not fire the Stapler. Remove and do not continue to use the Stapler.

#### NOTE:

- The jaws may be repositioned on the tissue prior to firing by pressing the Open/Back Button, which will reopen the jaws.
- During the 15 seconds pre-pressing of tissue, the Safety Indicator light flashes. After completion, the Safety Indicator light goes out and an electronic prompt sound "beep" can be heard.
- 3. Firing and releasing stapler:
  - a. Press the green Safety Button to enable the firing function. It is ready to be fired when the Safety Indicator light will remain on and two electronic beeps are heard. (refer to Figure 10).
  - b. Press and hold the Close/Fire Button to fire and slide the Knife Blade Indicator to the End of Cut Line (refer to Figures 11 and 12).

#### WARNING:

Incomplete firing may result in malformed Staples, incomplete cut line, bleeding, and/or difficulty removing the Stapler.

#### NOTES:

- The instrument will not fire the Staples and cut tissue unless the green Safety Button is pressed.
- Make sure the Stapler is completely fired. The time to complete the firing depends on length of the Reloads (30mm, 45mm or 60mm) of Cartridges and tissue status.
- If firing stalls and fails, the system will stop. Try to fire again, if still unable to complete the firing, stop firing, press and hold the Open/Back Button to move backward. If it is still unable to move backward, refer to section "F. Manually Overriding System".
- When Knife Blade Indicator slides to the End of Cut Line, Forward Indicator and Safety Indicator turn off, only Power Indicator is on, and an electronic prompt sound "beep" can be heard.
- Users can continuously or pulse press the Close/Fire Button during the firing process to complete the firing.
- c) Once the Stapler has been completely fired, release the Close/Fire Button to return the Cutting Knife to the starting position (i.e., before firing).
- d) Press the Open/Back Button to open the jaws (refer to Figure 13). Gently pull the Stapler away from the transected tissue and ensure the tissue is released from the jaws.

**CAUTION:** Examine the staple lines for pneumostasis/hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.

- 4. To remove the Stapler from the body cavity:

- a. Ensure the Articulation Knob is in the neutral position, so that the jaws are straight (refer to Figure 14).

**CAUTION:** Failure to have the jaws in the straight position will result in difficult insertion or withdrawal of the Stapler and may result in damage to the Stapler.

- b. Press and hold the Close/Fire Button to close the jaws, and completely remove the Stapler while in the closed position.

**CAUTION:** When the Stapler needs to be used multiple times in a single operation, used Reloads must be unloaded and replaced with new ones. The device's lock-out feature is designed to prevent an empty single use reload from being fired a second time.

### E. Reloading the Instrument

- 1. Confirm Stapler position:

- Ensure that the jaws are open. Press the Open/Back Button, if needed.
  - Ensure that the Articulation Knob is in the neutral position.
- Pull the light blue Load/Lock Button back towards the instrument, twist the Reload counterclockwise 45° and remove the Reload from the Shaft of the instrument (refer to Figure 15).
  - Using sterile technique, remove the new Reload from the package. To avoid damage, carefully place the Reload in the sterile field.
  - To load the new Reload, repeat section B. Assembling the Stapler, steps 3 to 7.

#### F. Manual override system

- Before using the manual override, if the Cutting Knife is unable to return to home position and the jaws will not open, try the following:
  - Ensure the Battery Pack is securely installed (refer to Figure 4) and the Stapler has power, then try pressing the Open/Back Button again (refer to Figure 13).
- If the Cutting Knife still does not return, use the manual override:
  - Remove the Reset Cover labeled "MANUAL OVERRIDE" on the top of the Stapler Handle (refer to Figure 16).
  - Insert the Reset Wrench into the square hole ① (see Figure 16). Turn the Reset Wrench counterclockwise until it stops.
  - Insert the Reset Wrench into the square hole ② (see Figure 16), and turn the Reset Wrench back and forth counterclockwise until the Cutting Knife completely returns to its home position. The jaws should now be open, releasing the tissue.
  - After releasing the tissue, move the jaws away, then remove the Stapler from the body cavity.
  - Discard the Stapler.

**CAUTION:** After the manual override system is used, the Stapler is disabled and cannot be used for any subsequent firings.

#### G. Battery Pack Disposal

**NOTE:** The Battery Pack contains a built-in battery drain that will turn on once installed in the Stapler. Once it is first turned on, the Battery Pack will continue to drain even after it has been removed from the Stapler. Once removed, depending on your local regulations, it can be put directly into the battery recycling bin or normal waste stream.

- Before installing into the Stapler:  
If Battery Pack needs to be disposed of prior to installation into instrument (e.g. product is beyond expiration date indicated on the package, Battery Pack is dropped), first install the Battery Pack into the Stapler, then remove, to trigger the built-in battery drain.
- After using the Stapler:  
The Battery Pack must be removed from the instrument prior to disposal.  
If the Battery Pack requires decontamination prior to disposal, follow hospital protocol.  
To remove the Battery Pack, press the buckle spring plate and pull the Battery Pack backwards (refer to Figure 4). There is no need to disassemble the Battery Pack.

#### WARNING:

Do not use autoclaves, ethylene oxide, or radiation to sterilize or disinfect Battery Pack.

### Technical Specifications

#### Electrical safety performance:

- The general electrical safety performance of the easyEndo™ UniPlus Single Use Powered Endoscopic Staplers complies with the requirements of IEC 60601-1:2005+AMD1:2012+AMD2:2020. See Table 5 below for details.

Table 5 Electrical safety characteristics

This product has following safety characteristics:

| Serial number | Safety characteristics                                                                                                                                                     | Parameter                           |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------|
| 1             | Classified according to type of protection against electric shock                                                                                                          | Internally powered equipment        |
| 2             | Classified according to degree of protection against electric shock                                                                                                        | Type CF applied part                |
| 3             | Classified according to degree of protection against ingress of water                                                                                                      | IPX0                                |
| 4             | Classified according to safety degree when used in the environment with flammable anesthetic mixture with air or flammable anesthetic mixture with oxygen or nitrous oxide | Non-AP/APG equipment                |
| 5             | Classified according to mode of operation                                                                                                                                  | Continuous operation                |
| 6             | Input voltage and frequency                                                                                                                                                | 12VDC                               |
| 7             | Applied part with protection against the effect of the defibrillation discharge:                                                                                           | No                                  |
| 8             | Signal output or input part                                                                                                                                                | No                                  |
| 9             | Permanently installed equipment or non-permanently installed equipment                                                                                                     | non-permanently installed equipment |

#### 2. Electromagnetic compatibility requirements:

For this instrument, special precautions regarding electromagnetic compatibility (EMC) must be taken, and it must be installed and used in accordance with the electromagnetic compatibility information specified in this instruction.

Medical electrical equipment requires special EMC precautions and must be installed and used in the EMC environment specified in this document.

Portable and mobile radio frequency communication equipment may affect this instrument.

**NOTE:** The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

#### WARNINGS:

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the MANUFACTURER of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of this equipment, including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could result.

#### Essential performance:

- No Cutting Knife movement without activation of the Close/Fire Button and Open/Back Button.
- When there is a Stapler failure, manual forced reset can return the Cutting Knife and open the jaws of the Stapler.
- No tissue will be cut without activation of the Safety Button.

Electromagnetic Emissions Guidance:

|                                                                    |
|--------------------------------------------------------------------|
| Guidance and manufacturer's declaration – electromagnetic emission |
|--------------------------------------------------------------------|

The easyEndo™ UniPlus Single Use Powered Endoscopic Staplers is intended to be used in the electromagnetic environment specified below. The customer or the user should guarantee that it is used in such an environment:

| Emission test                                      | Compliance     | Electromagnetic environment - guidance                                                                                                                                                                                                                   |
|----------------------------------------------------|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emission CISPR 11                               | Group 1        | The easyEndo™ UniPlus Single Use Powered Endoscopic Staplers uses radio frequency energy only for its internal function. Therefore, its radio frequency emission is very low and is not likely to cause any interference to nearby electronic equipment. |
| RF emission CISPR 11                               | Class A        |                                                                                                                                                                                                                                                          |
| Harmonic emission IEC 61000-3-2                    | Not applicable | The easyEndo™ UniPlus Single Use Powered Endoscopic Staplers can't be connected to the public low-voltage power supply network                                                                                                                           |
| Voltage fluctuation/flicker emission IEC 61000-3-3 | Not applicable |                                                                                                                                                                                                                                                          |

Electromagnetic Immunity Guidance:

| Guidance and manufacturer's declaration - electromagnetic immunity                                                                                                                                                        |                                                                      |                                                                                                                                                |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| The easyEndo™ UniPlus Single Use Powered Endoscopic Staplers is intended to be used in the electromagnetic environment specified below. The customer or the user should guarantee that it is used in such an environment: |                                                                      |                                                                                                                                                |
| Immunity test                                                                                                                                                                                                             | Compliance level                                                     | Electromagnetic environment - guidance                                                                                                         |
| Electrostatic discharge IEC 61000-4-2                                                                                                                                                                                     | ± 8 kV contact discharge<br>±2 kV/±4 kV, ±8 kV, ±15 kV air discharge | Floor should be wood, concrete or ceramic tile. If the floor is covered with synthetic material, the relative humidity should be at least 30%. |
| Radiated RF EM fields IEC 61000-4-3                                                                                                                                                                                       | 3V/m, 80MHz – 2.7 GHz<br>80%AM at 1kHz                               | Portable and mobile RF communications equipment should be used no closer to any part of Staplers                                               |
| Proximity fields From RF wireless communications equipment                                                                                                                                                                | Refer to next page                                                   |                                                                                                                                                |
| Power frequency magnetic field (50 Hz/60 Hz) IEC 61000-4-8                                                                                                                                                                | 30 A/m                                                               | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.      |
| Proximity magnetic fields IEC 61000-4-39:2017                                                                                                                                                                             | 134.2 kHz, PM, 2.1 kHz, 65 A/m<br>13.56 MHz, PM, 50 kHz, 7.5 A/m     |                                                                                                                                                |

Recommended separation distance between portable and mobile RF communication equipment and the powered endoscopic linear cutting staplers and components for single use

Nowadays, many types of RF wireless equipment are being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. The easyEndo™ UniPlus Single Use Powered Endoscopic Staplers have been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2020. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and The easyEndo™ UniPlus Single Use Powered Endoscopic Staplers recommended below.

| Test frequency(MHz) | Band(MHz) | Service                                                         | Modulation                         | Immunity test level(V/m) |
|---------------------|-----------|-----------------------------------------------------------------|------------------------------------|--------------------------|
| 385                 | 380-390   | TETRA 400                                                       | Pulse modulation<br>18Hz           | 27                       |
| 450                 | 430-470   | GMRS 460<br>FRS 460                                             | FM<br>± 5 kHz deviation 1 kHz sine | 28                       |
| 710                 | 704-787   | LTE Band 13, 17                                                 | Pulse modulation<br>217Hz          | 9                        |
| 745                 |           |                                                                 |                                    |                          |
| 780                 |           |                                                                 |                                    |                          |
| 810                 |           |                                                                 |                                    |                          |
| 870                 | 800-960   | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5          | Pulse modulation<br>18Hz           | 28                       |
| 930                 |           |                                                                 |                                    |                          |
| 1720                |           |                                                                 |                                    |                          |
| 1845                |           |                                                                 |                                    |                          |
| 1970                | 1700-1990 | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS | Pulse modulation<br>217Hz          | 28                       |
| 2450                | 2400-2570 | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7            | Pulse modulation<br>217Hz          | 28                       |
| 5240                | 5100-5800 | WLAN 802.11 a/n                                                 | Pulse modulation<br>217Hz          | 9                        |
| 5500                |           |                                                                 |                                    |                          |
| 5785                |           |                                                                 |                                    |                          |

#### Operating Environment Conditions

The Stapler and Reload for Single Use are generally used in professional healthcare environments such as hospitals and clinics.

- Temperature: 10°C – 40°C;
- Relative humidity: 30% – 75%;
- Range of atmospheric pressure: 700hPa – 1,060hPa;

#### Storage conditions

Store at room temperature, and storage environment parameters are as follows.

- Ambient temperature: 5°C – 30°C;
- Relative humidity: 10% – 80%;
- Range of atmospheric pressure: 700hPa – 1,060hPa;

#### Transportation conditions

This product should be prevented from pressure, direct sunlight, rain and snow during transportation.

- Ambient temperature: -20°C – 60°C;
- Relative humidity: 5% – 80%;
- Range of atmospheric pressure: 700hPa – 1060hPa;

## Dispose after use

This device is for single use only. After use, please dispose of it according to hospital regulations. Do not discard it arbitrarily.

## Contact information

Ezisurg Medical Co., Ltd.  
No.399, Miaoqiao Rd., Pudong New District, Shanghai, P.R. China  
Website: www.ezisurg.com  
Tel: +86-21-50456176  
Fax: +86-21-50676156

**EU REP** Shanghai International Holding Corp. GmbH (Europe)  
Add: Eiffestrasse 80, 20537 Hamburg, Germany  
Tel: +49-40-2513175  
Fax: +49-40-255726

After launch of European Database on Medical Devices/EUDAMED, the summary of safety and clinical performance (SSCP) will be made available to the public at <https://ec.europa.eu/tools/eudamed>.

For a patient/ user/ third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/ or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: [https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en)

For additional information, contact your Ezisurg Medical Co., Ltd. Sales Representative or Ezisurg Medical Co., Ltd. Customer Service Tel: +86-400906176.

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