Shanghai International Holding Corp. GmbH (Europe)

Add: Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175

Fax: +49-40-255726

Ezisurg Medical Co., Ltd.

No.399, Miaojiao Rd., Pudong New District, Shanghai, P.R. China

Website: www.ezisurg.com

Tel: +86-21-50456176


Fax: +86-21-50676156

CE 0197

Rev.D
2026-04-02

easyEndo™ UniPlus


en	Single Use Reloads
pl	Wkłady Jednorazowego Użytku
cs	Jednorázové náboje
pt-pt	Recargas Descartáveis
es	Recargas de un solo uso
it	Ricariche Monouso
de	Einweg-Nachlademagazine
fr	Chargeur à usage unique
nl	Wegwerpnavullingen
sv	Engångsreloads
ro	Reîncărcări de Unică Folosință
el	Κασέτες συνδετήρων μίας χρήσης
hu	Egyszer használatos újratöltők

IMPORTANT!

Please read the following instructions carefully before using this product. This document is designed to provide instructions for using this product. It is not a reference to surgical techniques. Failure to follow the instructions may lead to serious surgical consequences, such as anastomotic potential for bleeding, leakage, or disruption.

Device Overview

Device Description

easyEndo™ UniPlus Single Use Reloads are sterile, single patient use instruments that simultaneously cut and staple tissue. The instruments can be used in open or minimally invasive surgery. The Reload is not preloaded with a Curved Tip, which can be installed as needed.

Table 1. List of device features

1. End of Staple Line	2. End of Cut Line	3. Knife Blade Indicator
4. Cutting Knife	5. Curved Tip*	6. Cartridge
7. Anvil Jaw	8. Tissue Stop Wedge	9. Loading Indicating Line
10. Loading Indicating Arrow	11. Shipping Wedge	12. Staple

*Curved Tip: The easyEndo™ UniPlus Reloads is not preloaded with a Curved Tip, which is placed inside the easyEndo™ UniPlus Staplers packaging and can be installed as needed (refer to section B and Step 3).

Table 2. Reload Product Codes Chart

REF	Cartridge Color	Staple Line Length (L2) (mm)	Staples Rows (pattern)	QTY / row	Staple height H (mm)		Thickness range* (mm)
					Open	Closed	
P30M	Grey	30	6(staggered)	8	2.0	0.75	0.75 - 1.0
P30W	White				2.5	1.0	1.0 - 1.5
P30B	Blue				3.5	1.5	1.5 - 2.0
P30VM	Brown	45	6(staggered)	11	Inner: 2.0 Middle: 2.5 Outer: 3.0	0.75 1.0 1.25	0.88-1.5
P45M	Grey				2.0	0.75	0.75 - 1.0
P45W	White				2.5	1.0	1.0 - 1.5
P45B	Blue				3.5	1.5	1.5 - 2.0
P45C	Gold				3.8	1.7	1.7 - 2.2
P45G	Green				4.1	2.0	2.0 - 2.5
P45T	Black				4.4	2.3	2.3-3.0
P45VM	Brown				Inner: 2.0 Middle: 2.5 Outer: 3.0	0.75 1.0 1.25	0.88-1.5
P45MT	Purple				Inner: 3.5 Middle: 3.8 Outer: 4.1	1.25 1.5 1.75	1.5-2.25
P45XT	Black				Inner: 3.8 Middle: 4.1 Outer: 4.4	1.75 2.0 2.25	2.25-3.0

P60M	Grey	60	6(staggered)	15	2.0	0.75	0.75 - 1.0
P60W	White				2.5	1.0	1.0 - 1.5
P60B	Blue				3.5	1.5	1.5 ~ 2.0
P60C	Gold				3.8	1.7	1.7 ~ 2.2
P60G	Green				4.1	2.0	2.0 ~ 2.5
P60T	Black				4.4	2.3	2.3-3.0
P60VM	Brown				Inner: 2.0 Middle: 2.5 Outer: 3.0	0.75 1.0 1.25	0.88-1.5
P60MT	Purple				Inner: 3.5 Middle: 3.8 Outer: 4.1	1.25 1.5 1.75	1.5-2.25
P60XT	Black				Inner: 3.8 Middle: 4.1 Outer: 4.4	1.75 2.0 2.25	2.25-3.0
*Tissue thickness range: Porcine tissue measured at 8g/mm² prior to firing.							

Compatibility

The easyEndo™ UniPlus Single Use Reloads can be used with easyEndo™ UniPlus Single Use Powered Endoscopic Staplers. Refer to the Instructions for Use for the easyEndo™ UniPlus Single Use Powered Endoscopic Staplers for specific indications, contraindications, Warnings and Precautions, and operating instructions.

Refer to the Instructions for Use for the easyEndo™ UniPlus Single Use Powered Endoscopic Staplers for articulation angles. Compatibility of other manufacturer's staples and staplers has not been determined.

Intended Use/Purpose

The instruments are intended for transection, resection, and/or creation of anastomoses in open and minimally invasive general abdominal, and thoracic surgical procedures.

Indications

The easyEndo™ UniPlus Single Use Reloads have applications in abdominal, and thoracic surgery for resection, transection and creation of anastomosis. It may be used for transection and resection of liver substance, hepatic vasculature, and biliary structures and for transection and resection of the pancreas.

Contraindications

- Do not use the instruments on the aorta.
- The device should not be used to staple tissues that are necrotic, friable, or have altered integrity, e.g., ischemic or edematous tissues.
- Do not use any linear cutter on major vessels without making provision for proximal and distal control.
- Tissue thickness should be carefully evaluated before firing any stapler. Refer to Table 2. Reload Product Codes Chart for tissue compression requirement (closed staple height) for each staple size. If tissue cannot comfortably compress to the closed staple height, or easily compresses to less than the closed staple height, the tissue is contraindicated as it may be too thick or too thin for the selected staple size.
- These instruments are not intended for use when surgical stapling is contraindicated.

Target Patient Group

The intended patient population is adults.

Intended clinical benefits

Support for Successful Anastomosis.

Clinical Benefits

- Reduced anastomotic bleeding and anastomotic leakage rate;
- Reduced hospital stay and operative time.

Undesirable Side Effects/Residual Risks

Undesirable side effects and risks associated with powered surgical stapling include the potential for bleeding, tissue injury, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, hypersensitivity, leakage, Biological toxicity, property or environmental damage, electrical shock, and foreign body magnetic resonance incompatibility. Also, unintended harm, extended surgery, or altered surgical approach may result from staple line failures, inability to cut, or damaged devices.

Safety Information

The Use of Warning, Caution, and Note Statements

Information relative to the completion of a task in a safe and thorough manner will be supplied in the form of a Warning, a Caution, or a Note statement. These statements are found throughout the documentation.

These statements should be read before continuing to the next step in a procedure.

- WARNING:** A Warning statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in personal injury or loss of life.
- CAUTION:** A Caution statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in damage to or destruction of the equipment.
- NOTE:** A Note statement indicates an operating or maintenance problem, practice, or condition that is necessary to accomplish a task efficiently.

Warnings and Precautions

- Ensure that the Reloads are compatible with the Staplers.
- Depending upon the stapler handle used, also refer to the applicable Instructions for Use for the compatible Stapler for specific indications, contraindications, warnings, precautions, and operating instructions.
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers (such as trocars) are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- Radiation therapy before surgery may cause tissue changes, resulting in tissue thickness exceeding the applicable range of the selected staple. Therefore, when choosing the staple specification, consider the patient's various preoperative treatments.
- This device is packaged and sterilized for single use only. Do not re-use, reprocess, or re-sterilize. Re-use, reprocessing, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure, that in turn may result in patient injury, illness, or death. Also, reprocessing or re-sterilizing single-use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to the patient's injury, illness, or death.
- Do not use if sterile package is damaged.
- Do not use after expiry date.
- Always select a reload with the appropriate Staple size for the tissue thickness. Refer to Table 2. Reload Product Codes Chart for proper Cartridge selection. Overly thick or thin tissue may result in unacceptable Staple formation.
- The use of Curved Tip Stapler is contraindicated on tissue or structures that cannot fit completely within the jaws proximal to the transitional angle of the Curved Tip.
- The Reload is packaged in the open position. Do not attempt to close the Reload.
- Do not attempt to remove the Shipping Wedge until the Reload is loaded into the instrument.
- Do not clamp instrument prior to removing Shipping Wedge.
- When using the Reload with a Curved Tip, ensure that the tissue/vessel to be transected does not extend beyond the End of Cut Line on the Reload. The device will only cut to the cut line. Tissue contained within the jaws distal to this line will not be

transected.

15. Ensure that the Stapler is loaded firmly with the Reload. The application of the instrument loaded not firmly during the operation may lead to separation of the Stapler and the Reload, in such case, an open surgery may be necessary to remove the separated Reload.
16. It is strictly prohibited to use the Reload if Staple is exposed to the surface of Cartridge or the colored drivers can be seen.
17. The instrument Jaw must be completely visible past the trocar sleeve before opening within the body cavity.
18. Clamping and unclamping delicate structures such as venous and bile ducts may damage tissue regardless of whether the stapler is fired.
19. Prior to positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc., are incorporated into the instrument jaws. Firing over obstacles may result in incomplete cutting action and/or improperly formed Staples.
20. Placement of tissue proximal to the Tissue Stop Wedge (on the Jaws) may result in instrument malfunction. Any tissue extending beyond the cut mark on the Jaws, may be transected without Staples and may lead to serious surgical consequences, such as bleeding, leakage, or disruption.
21. When using the Reload with a Curved Tip to manipulate tissue, avoid exerting excessive pressure on fragile structures with the Curved Tip of the device.
22. If a Stapler malfunction occurs while applying staples across a blood vessel, then the user should clamp or ligate the vessel before releasing the stapler, while the Stapler is still closed on the tissue.
23. Ensure that the tissue lies flat and is positioned properly between the jaws. Any "bunching" of tissue along the Cartridge, particularly in the crotch of the jaws, may result in an incomplete staple line.
24. Before firing, check that no unexpected tissues and organs are clamped by the jaws.
25. Incomplete firing may result in malformed staples, incomplete cut line, bleeding, and/or difficulty removing the device.
26. Examine the staple lines for pneumostasis/hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.
27. Make sure that the Curved Tip is correctly and securely connected to the Anvil Jaw.
28. When using the Curved Tip to manipulate tissue, do not apply excessive force to the Curved Tip, excessive force may separate it from the anvil Jaw.
29. Before removing the instrument, be sure tissue is cleared from the jaws, and then close the jaws. Failure to gently clear tissue from the jaws may result in bleeding or damage to adjacent tissues or structures. Additionally, it may interfere with removal of the instrument through the trocar.
30. The Reload which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
31. There is a risk of increased leak rates when staple lines are crossed, even if there may be clinical circumstances when a surgeon may deem it necessary or appropriate to do so.
32. Cobalt classified as a CMR substance is present in the device in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

Use instructions

A. Prior to Using Instrument

1. Assess the tissue prior to using the instrument. The instrument should not be used to staple tissues that are necrotic, friable, or have altered integrity, e.g., ischemic or edematous tissues, as the staples may be less likely to securely approximate tissue.

B. Assembling the Instrument

1. Using sterile technique, remove the Reload from their packages. To avoid damage, carefully place the Reload in the sterile field.
2. Select a new Reload with the appropriate staple size for the tissue thickness. Refer to Table 2. Reload Product Codes Chart.

WARNING:

Carefully evaluate the tissue thickness prior to use. Overly thick or thin tissue may result in unacceptable staple formation.

CAUTION:

- The Reload is packaged in the open position. Do not attempt to close the Reload.

- Do not attempt to remove the Shipping Wedge until the Reload is loaded into the instrument.
- Do not clamp instrument prior to removing shipping wedge.

NOTE: The color of the Cartridge identifies the staple size, as described in Table 2. Reload Product Codes Chart.

NOTE: The intended tissue in the label are recommended approximate intended tissues.

3. Insert the Curved Tip into the Anvil Jaw, as needed.
 - a. Firstly, insert the distal end of Anvil Jaw into the Curved Tip opening cavity above the spring plate (refer to Figure 2a);
 - b. Push the Curved Tip horizontally, and the spring plate on the Curved Tip follows the slot on the back of Anvil Jaw (refer to Figure 2b,2c);
 - c. An audible click will be heard when the Curved Tip is fully inserted.

WARNING:

When using the Reload with a Curved Tip, ensure that the tissue/vessel to be transected does not extend beyond the End of Cut Line on the Reload. The device will only cut to the cut line. Tissue contained within the jaws distal to this line will not be transected.

Caution: Make sure that the Curved Tip is correctly and securely connected to the Anvil Jaw.

NOTE:

- The Curved Tip is placed inside the easyEndo™ UniPlus Staplers packaging and can be installed as needed.
- When using the Reload with a Curved Tip, the Curved Tip of the Anvil Jaw can be used to manipulate the tissue and will be clearly visible when the reload is clamped on the tissue/vessel that has to be transected.
- When using the Reload with a Curved Tip, the Curved Tip of the Anvil Jaw can be used for blunt dissection.
- 4. Dependent upon the stapler used, refer to the applicable Instructions for Use for the compatible Stapler for specific indications, contraindications, warnings and precautions, and operating instructions.

C. Reloading the Instrument

1. Dependent upon the stapler used, refer to the applicable Instructions for Use for the compatible Stapler for specific indications, contraindications, warnings and precautions, and operating instructions.

Technical Specifications

Table 3 MRI Safety Information

 MR Conditional	
Non-clinical testing of the worst-case implants has demonstrated that the implant matrix "Staples" is MR Conditional. A patient with staples can be safely scanned in an MR system meeting the following conditions. Failure to follow these conditions may result in injury to the patient.	
Item Name/Identification	Staples
Static Magnetic Field Strength [T]	1.5 T or 3.0 T
Maximum Spatial Field Gradient [T/m]	55 T/m (5,500 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Whole body transmit coil (short body coil), Head RF receive coil or Extremity RF receive coil
Maximum Whole Body SAR [W/kg]	whole body SAR < 2.0 W/kg at 1.5 T and 3.0 T
Maximum B1 ⁺ _{rms} [μT]	B1 ⁺ _{rms} < 27.1 μT at 1.5 T Not reported for 3.0 T
Scan Duration	for 15 minutes of continuous RF (a sequence or back to back series/scan without breaks)
MR Image Artifact	The presence of this implant may produce an image artifact of 1.92 mm.
If information about a specific parameter is not included, there are no conditions associated with that parameter.	

NOTE: All parameters are connected by AND logical conjunction and have to be within the limits always at the same time. For MR systems that use SAR only, SAR values shall be used. For MR systems that use B1⁺_{rms} and SAR, B1⁺_{rms} can be used.

Note: Staples are safe as implant components under specific conditions in 1.5T and 3.0T MRI environment.

Operating Environment Conditions

The Reload for Single Use are generally used in professional healthcare environments such as hospitals and clinics.

- a) Temperature: 10°C – 40°C ;
- b) Relative humidity: 30%–75%;
- c) Range of atmospheric pressure: 700hPa–1,060hPa;

Storage conditions

Store at room temperature, and storage environment parameters are as follows.

- a) Ambient temperature: 5°C – 30°C ;
- b) Relative humidity: 10%–80%;
- c) Range of atmospheric pressure: 700hPa–1,060hPa;

Transportation conditions

This product should be prevented from pressure, direct sunlight, rain and snow during transportation.

- a) Ambient temperature: -20°C – 60°C ;
- b) Relative humidity: 5%–80%;
- c) Range of atmospheric pressure: 700hPa–1060hPa;

Dispose after use

This device is for single use only. After use, please dispose of it according to hospital regulations. Do not discard it arbitrarily.

Contact information

Ezisurg Medical Co., Ltd.

No.399, Miaoqiao Rd., Pudong New District, Shanghai, P.R. China

Website: www.ezisurg.com

Tel: +86-21-50456176

Fax: +86-21-50676156

Shanghai International Holding Corp. GmbH (Europe)

Add: Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175

Fax: +49-40-255726

After launch of European Database on Medical Devices/EUDAMED, the summary of safety and clinical performance (SSCP) will be made available to the public at <https://ec.europa.eu/tools/eudamed>.

For a patient/ user/ third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/ or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

For additional information, contact your Ezisurg Medical Co., Ltd. Sales Representative or Ezisurg Medical Co., Ltd. Customer Service Tel: +86-4000906176.

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