easyUS™ Single Use Ultrasonic Shears Instructions for Use

REF

SHC14A、SHC23A、SHC36A、SHC45A

Ezisurg Medical Co., Ltd.

Please read all information carefully before use.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician

Important!

Please read all information carefully.

1. This manual is designed to provide instructions for use of the easyUS™ Single Use Ultrasonic Shears SHC14A、SHC23A、SHC36A and SHC45A of Ultrasonic surgical system. It is not a reference to surgical techniques.
2. The shears SHC14A、SHC23A、SHC36A and SHC45A is designed for single use only. To reuse or reassemble the shears might introduce functional failure of the instruments and create a risk on patient of the damage as well as the infectious potential. Please don't reuse, reassemble and re-sterilize the instruments.
3. The easyUS™ Single Use Ultrasonic Shears should be used with the generator and handpiece of Ultrasonic surgical system. Please refer *easyUS™ Ultrasonic surgical system Generator User Manual* and *easyUS™ Ultrasonic surgical system Handpiece Instructions for Use for eHP01* for more information.

Device Description

The Ultrasonic surgical system converts high frequency current into ultrasonic oscillation to mechanically fragment soft-tissue cells upon contact with a vibrating tip for cutting and coagulating tissue during surgery. It mainly consists of a generator (with software embedded), a footswitch, a handpiece and a single use shears. The shears and the handpiece have a torque wrench and/or a test tip as spare parts. The torque wrench is to secure the test tip or shears to the handpiece. The test tip is to connect the handpiece and is used to help to calibrate and diagnose the handpiece performance. Among them, the generator outputs electrical energy (controlled by the software); the footswitch is used to active MAX or MIN output; the handpiece converts the electrical energy to mechanical motion; the shears connects to the handpiece and pass the mechanical motion to the distal end of the blade in the shears.

The Ultrasonic surgical system does not include an integrated irrigation/suction system.

The Ultrasonic surgical system does not support the function of patient physiological parameter monitoring.

The easyUS™ Single Use Ultrasonic Shears had different models and the surgeon can choose the appropriate model by clinical preference.

The diameter of the shears shaft (see D1 in Figure 1) is 5.5mm. The shears can use with

5mm trocar or a larger one with an adaptor. It is also able to use the shears directly through the incision. The trocar is surgical accessory and not included in the system.

The shaft of the easyUS™ Single Use Ultrasonic Shears, including blade, clamp, tissue pad, trigger tube and outer tube are applied parts.

See Figure 1 for the diagram of the easyUS™ Single Use Ultrasonic Shears.



Figure 1 Diagram of easyUS™ Single Use Ultrasonic Shears

See Table 1 for the product codes and dimensions of the shears.

Table 1 Specifications of easyUS™ Single Use Ultrasonic Shears

Product Code	Shaft Length L1	Shaft Diameter D1
SHC14A	140 mm	5.5 mm
SHC23A	230 mm	5.5 mm
SHC36A	360 mm	5.5 mm
SHC45A	450 mm	5.5 mm

The product code of shears consists of the shears type, shaft length of the shears and the character code.

Note 1: Character Code indicates different types for the shears with same dimensions. For example, A is for easyUS™ Single Use Ultrasonic Shears.

Note 2: Shears type indicates the shape of the distal end of the blade in the shears. For example, SHC is for Shears Curved.

Theory of Operation

The Ultrasonic surgical system is for hemostatic cutting and/or coagulation. The blade in shears vibrates longitudinally at 55.5 kilohertz. This ultrasonic vibration at the blade enhances its cutting ability. The same vibration seals small vessels with coagulated blood and tissue proteins. Hemostasis occurs when tissue couples with the instrument. This coupling causes collagen molecules within the tissue to vibrate and become denatured, forming a coagulum.

Device performance

The easyUS™ Single Use Ultrasonic Shears can output ultrasonic energy stably to incise soft tissue and seal small vessels when activate the hand switch or footswitch after connecting to the handpiece and generator.

Indications for Use

easyUS™ Single Use Ultrasonic Shears, when used with the Ultrasonic Surgical System devices, are indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The instruments can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels. easyUS™ Single Use Ultrasonic Shears can be used to coagulate isolated vessels up to and including 5 mm in diameter.

Contraindications

1. The instruments are not indicated for incising bone.
2. The instruments are not intended for contraceptive tubal occlusion.
3. The Ultrasonic Surgical Systems instruments are contraindicated when, in the judgment of the physician, ultrasonic surgery would be contrary to the best interest of the patient.

Target Patient Group

The intended patient population consists of adult and pediatric patients.

Intended Users

Intended users are medical professionals qualified for the ultrasonic surgical technique and surgical procedure to be performed.

Clinical Benefits

Support for successful cutting and coagulation.

Undesirable Side Effects/ Residual Risks

Potential bleeding, tissue injury, introduction of non-sterile surfaces, property or environmental damage, unintended harm, extended surgery or altered surgical approach may result from incorrect use or damaged devices.

Operation

① See Figure 2 for each part of easyUS™ Single Use Ultrasonic Shears


Figure 2 Construction of easyUS™ Single Use Ultrasonic Shears

The shears consists of a handle, a handpiece interface, a hand switch, a trigger, a knob and a shears shaft. The shears shaft consists of a blade, an inner tube, an outer tube and a jaw which is consisted of a clamp and a tissue pad.

② Installation

Clasp the trigger to close the jaw. Insert the torque wrench from the distal end of the shears. Insert handpiece into the handpiece interface on shears. Rotate the torque wrench clockwise to secure the shears and the handpiece until to hear at least two notification tone “Ka-”. Clasp the trigger to close the jaw again to remove torque wrench from the distal end of the shears.

To install the handpiece and the generator, please make sure the handpiece connector is inserted into the receptacle along the location line. Fully insert the handpiece connector.

③ Pre-Run Test

Power up the generator, hold the Standby touch key for 3 seconds to switch the generator to Ready mode. Ensure the shears jaw is open and not in contact with any objects. Activate system by hand switch or footswitch. The system will perform Pre-Run Test and play a test tone until the test is passed and system enters Running mode. If test failed the system enters Error mode and prompts failed error information. Please refer *Ultrasonic surgical system Generator User Manual* for troubleshooting.

④ Use

Clasp the trigger to close the jaw. Insert the shears into a 5mm Trocar or a larger one with appropriate adapter.

In endoscopic surgery, the shears needs to be used in conjunction with the endoscope. The operator locates the target requiring cutting or coagulation through the endoscopic imaging.

Generally, operator should activate the system after the jaw is fully closed (clasp the trigger until sense the force feedback and/or hear the notification tone) when using shears to cut tissue. Another way is to keep the jaw open and use another side of the blade to cut tissue

(known as Back Cut).


Figure 3 Diagram of Shears Jaw

When using shears to grip tissue to cut and/or coagulate, please ensure the tissue is not beyond the best working range (see Figure 3). Deeply grip tissue might decrease the effect of cutting and/or coagulating.

Stop activating immediately after cutting and/or coagulating operation is done. Otherwise the tissue pad might be worn out soon.

When tissue stuck on the distal end of the blade hence decrease the effect of cutting and/or coagulating, please remove the shears from patients, open the jaw, put the distal end of the blade inside normal saline. Then activate system with MAX output to clear the stuck tissue.

Note: During activating, don't have the blade in contact with any metal object.

⑤ Finishing

After surgery is done, close the jaws by pressing trigger and remove the shears from trocar. Detach handpiece from shears with torque wrench. The easyUS™ Single Use Ultrasonic Shears are disposable. Disposal should be performed according to hospital requirements and regulations. The handpiece is able to reuse. Please refer *easyUS™ Ultrasonic surgical system Handpiece Instructions for Use* to get information about how to clean, disinfect and sterilize the handpiece.

Precautions

1. Minimally invasive procedures should be performed only by persons have adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
2. Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure. Please check carefully on the invasive parts, especially on the surface roughness, sharp edges or outcrop, which might cause safety hazard of the surgery.

3. Please use the torque wrench provided in the shears package to secure the shears and handpieces. Using torque wrench provided other than Ezisurg Medical Co., Ltd might not be able to secure shears and handpieces in appropriate torque, which might cause products broken or decrease the system performance.
4. Operators should perform system test when shears are first in use. During the system test, please keep the jaw opened, ensure the blade is not in contact with any other instruments, drapes, patients and other objects, Safety measures (in accordance with hospital protocol) taken in the presence of aerosols should be in effect during the system test.
5. When replace shears or re-secure shears and handpiece, please switch generator to Standby mode.
6. Ultrasonic surgical system is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. It is possible to create sparks by hitting other metal instruments. Sparks may ignite flammable gases such as bowel gas.
7. After acting on tissues for a long time, the distal end of the blade, the clamp and the distal end of outer tube might have higher temperature than 41 degree Celsius. Please be careful to operate and avoid to in contact with tissue, drapes, and other parts in patients.
8. When tissue stuck on the distal end of the blade hence decrease the effect of cutting and/or coagulating, please remove the shears from patients, open the jaw, put the distal end of the blade inside normal saline. Then activate system with MAX output to clear the stuck tissue.
9. Operator should activate the system after the jaw is fully closed (clasp the trigger until sense the force feedback and/or hear the notification tone) when using shears to cut tissue. Stop activating immediately after cutting and/or coagulating operation is done. Otherwise the tissue pad might be worn out soon.
10. When using shears to grip tissue to cut and/or coagulate, please ensure the tissue is not beyond the best working range. Deeply grip tissue might decrease the effect of cutting and/or coagulating.
11. Generally, it is better not to keep activating system longer than 15 seconds.
12. After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
13. As with all energy sources (Electrosurgery, Lasers, or Ultrasound) there are concerns about the carcinogenic and infectious potential of the by-products such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
14. It is not allowed to bend, sharpen or change the shape of the blade in anyway. Doing that might cause the shears broken and bring damage to either operators or patients.
15. Reuse, reprocessing or resterilization of easyUS™ Single Use Ultrasonic Shears is strictly PROHIBITED! The easyUS™ Single Use Ultrasonic Shears is disposable instrument. It is designed to be used in one surgery. The instrument should be

disposed after surgery. Reuse, reprocessing or resterilization will damage the structural integrity of the instrument and make it malfunction, which may result in injury, illness or death. In addition, it will cause contamination, infection or cross-infection, including, but not limited to, transmission of infectious diseases from one patient to another, which has high risk of illness or death. Ezisurg Medical Co., Ltd does not bear the liability for instrument misuse.

16. Do not use if sterilized package is damaged. Only assemble the shears in the sterile area before the procedure.
17. The Ultrasonic surgical system shall not be serviced or maintained while in use with a patient.
18. One or more components of this device contains the Cobalt (CAS No.7440-48-4; EC No.231-158-0) defined as CMR 1B in a concentration above 0.1% weight by weight.
19. Ensure the availability of backup shears in case of the device failure.
20. The shears belongs to the type CF applied part and meets the requirements for interconnection with endoscopes.
21. Before each use, it is necessary to confirm whether the specifications of the Trocar are suitable for the shears and the endoscope, and the puncture positions of the shears and the endoscope should be kept at an appropriate distance from each other.
22. Before each use, the outer surface of the portions of the shears should be checked to ensure there are no unintended rough surfaces, sharp edges or protrusions
23. When using the shears in conjunction with the endoscope, it is necessary to avoid any contact between the shears and the endoscope.
24. Successful hemostasis may require adjunct measures when the device is used on solid organs. Due to the difficulty of visualizing internal structures, proceed slowly and do not attempt to transect large masses of tissue in one activation. Avoid the division of large vascular/biliary bundles when using the instrument under these conditions.
25. Device has not been evaluated in main vessels of the central circulatory system and is not intended for use in the following named vessels: arteriae pulmonales, aorta ascendens, arcus aortae, aorta descendens to the bifurcatio aortae, arteriae coronariae, arteria carotis communis, arteria carotis externa, arteria carotis interna, arteriae cerebrales, truncus brachiocephalicus, venae cordis, venae pulmonales, vena cava superior, and vena cava inferior.

Sterilization

easyUS™ Single Use Ultrasonic Shears with Torque Wrench product is sterilized using ethylene oxide, The use-by date is marked on the packaging. The shelf life is 3 years. The product should not be used after expiration.

Disposal

easyUS™ Single Use Ultrasonic Shears with Torque Wrench product is one use only. Disposal should be performed according to hospital requirements and regulations.

Storage and transportation Conditions:

Avoid heavily stacking, exposing under sunlight, and raining/snowing. Alternatively, follow the order contract term.

During transportation, the temperature and humidity should be maintained within -20°C ~ 55°C and 10% ~ 80%, respectively. Avoid exposing under high-temperature for long time.

Storage Condition:

The product should be stored at the temperature of 5°C ~ 30°C and humidity of 10% ~ 80%.

No condensation; Atmospheric pressure: 700hPa to 1060hPa.

Ambient Operation Conditions:

Used in the operating rooms of large hospitals, specialized hospitals and other medical institutions that can perform open surgery and laparoscopic surgery, under normal operating room conditions.

Ezisurg Medical Co., Ltd.
Site: No. 399, Miaoqiao Rd., Pudong New District, 201315 Shanghai, P.R. China
Website: www.ezisurg.com


Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffestrasse 80, 20537, Hamburg, Germany
Tel: +49-40-2513175
Fax: +49-40-255726

For a patient/ user/ third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/ or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

For additional information, contact your Ezisurg Medical Co., Ltd. Sales Representative or Ezisurg Medical Co., Ltd. Customer Service.

Symbols

	Do not re-use	Affix on the sterile dialyzing paper label and package label
	Caution	Affix on the sterile dialyzing paper label and package label
	Contains hazardous substances	Affix on the bar code label and package label
	Keep dry	Affix on the transport package label
	Keep away from sunlight	Affix on the transport package label
	Use no hooks	Affix on the transport package label
	Stacking limit by number	Affix on the transport package label
	Use-by date	Affix on the package label, transport package label and bar code label
	Batch code	Affix on the package label, transport package label and bar code label
	Model number	Affix on the sterile dialyzing paper label, bar code label, package label and transport package label
	Do not re-sterilize	Affix on the sterile dialyzing paper label and package label

	Do not use if package is damaged and consult instructions for use	Affix on the sterile dialyzing paper label and package label
	Date of manufacture Country of manufacture: CHINA	Affix on the bar code label, package label and transport package label
	Manufacturer	Affix on the package label and transport package label
	Authorized representative in the European Community	Affix on the sterile dialyzing paper label and package label
	Device meets the requirements of Medical Device Directive 93/42/EEC, as amended by 2007/47/EC. 0197 is the code of notified body TUV Rheinland.	Affix on the sterile dialyzing paper label and package label
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Affix on the sterile dialyzing paper label and package label
	Unique device identifier	Affix on the bar code label, package label and transport package label
	Sterilized using ethylene oxide Single sterile barrier system with protective packaging outside	Affix on the package label
	Sterilized using ethylene oxide Single sterile barrier system	Affix on the sterile dialyzing paper label
 www.ezisurg.com	Consult instructions for use or consult electronic instructions for use	Affix on the sterile dialyzing paper label and package label
	Medical device	Affix on the package label

	Temperature limit	Affix on the transport package label
	Humidity limitation	Affix on the transport package label
	Atmospheric pressure limitation	Affix on the transport package label
	Type CF Applied Part	Affix on the sterile dialyzing paper label and package label
	Refer to instruction manual/booklet	Affix on the sterile dialyzing paper label and package label
	Package quantity	Affix on the package label and transport package label

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